Avior

User Guide

AVIOR

EN	User Guide	2
FR	Guide de l'utilisateur	40
DE	Bedienungsanleitung	78
IT	Guida per il paziente	116
ES	Guía del usuario	154
NL	Gebruikershandleiding	192

Contents2

1 Description and Intended Purpose4

2 Safety Information7

3 Avior User Controls and Charging Connections9

 3.1 AC-DC Battery Charger Adapter 10

4 Key Points 11

5 Battery Charging 13

6 How to Don and Doff the Device 14

7 Transportation and Storage 14

8 Device Operation 15

 8.1 Controls and Indicators 15

 8.2 Device Functions and Responses 16

9 Program Operation 25

 9.1 Operating Modes: User Modes and Functions 26

 9.2 Safety Mode 28

 9.3 How to Engage Stand-to Sit Mode in Normal Operation 28

10 SMARTSTEP—Wearer App 29

 10.1 App Overview 29

Appendix 1 Installing the SMARTSTEP Programming App 31

Appendix 2 Wearer Registration on SMARTSTEP 32

Appendix 3 SMARTSTEP—General Information 34

Manufacturer’s Declarations and Legal Information 37

About the Manual

Terminology

The following terms used throughout these instructions are defined as follows:

The terms *device* and *device type* are used to refer to Avior.

The term *tablet/phone* is used to refer to a tablet or mobile phone running an Android or iOS operating system.

The term *app* is used to refer to *SMARTSTEP*, the Blatchford programming application which can be installed on a *tablet/phone*.

The terms *user* and *wearer* are used to refer to the amputee who wears the device.

The term *clinician* is used in these instructions to refer primarily to the prosthetist, i.e., the *practitioner* qualified to prescribe and fit the device.

For each user access level to *SMARTSTEP*, a separate term is used, e.g., the term *clinician app* refers to *SMARTSTEP* as it appears and functions when accessed as an authorised practitioner/clinician. Similarly, the term *wearer app* refers to the app as it appears and functions when accessed respectively by an authorised wearer.

In this document, the *wearer app* is the primary app interface described

Text emphasis in bold italics is used to highlight a device/app button, selection option or switch name, e.g., ***SILENT*** mode; it is also used to highlight the ***Bluetooth®*** tradename.

This warning symbol highlights important safety information which must be followed carefully.

Scope

This User Guide provides instructions for use intended for the user only, and includes information relevant to the user about the device and about the *SMARTSTEP* wearer app.

Information intended specifically for the practitioner can be found in the Clinician's Manual 938474.

Note... As the app is designed for multiple screen aspect ratios and different operating system, and is being continuously improved, the layout and content of images of screens shown in this document may differ from how they are seen on your tablet/phone.

1 Description and Intended Purpose

These instructions for use are intended for the user. **Please keep these instructions.**

Please read and ensure that you understand these instructions, in particular all safety and maintenance information.

Application

The device is a Class 1 medical device according to the EU Medical Device Regulation 2017/745 to be used exclusively as part of a lower limb prosthesis.

Intended for a single user only (Single Patient – multiple use).

It is for use solely as a prosthetic knee joint for lower extremity use in accordance with the manufacturer's recommendations.

The device is a microprocessor-controlled knee which automatically modulates the resistance to flexion and extension at the knee joint to provide a locked or fixed swing resistance/walking cadence prosthesis as required by the wearer.

The device is suitable for use by transfemoral (or knee disarticulation) amputees weighing between 44 kg (97 lb) and 125 kg (275 lb).

Device Features

- Microprocessor-controlled stance and swing
- **Standing Support** mode that controls knee flexion to support standing
- Stumble recovery in all phases of gait
- **Bluetooth®** communications between programming app and the device
- **SMARTSTEP** wearer app interface
- User modes:
 - **Wheelchair** mode
 - **Manual Lock** (Flexion Lock) mode
 - **Silent** mode
- Up to two days of normal use from the internal rechargeable Li-ion battery pack
- Clinician-adjustable volume and frequency of audio indicator tones
- Clinician-adjustable intensity of haptic response vibrations
- Device* dust proof and splash proof to IP65

*Note... The battery charger is not water resistant and so charging must be carried out in a dry environment.

Device Set-up

The device will be fitted and initially configured for you by your prosthetist using the clinician's app. Your prosthetist may make adjustments to the performance of the device during your rehabilitation.

User Modes

In addition to configuring the normal automatic operating modes and functions and fine tuning the performance of the device to suit your needs and capabilities, your practitioner may enable the following additional automatic **User** modes:

- **Donning the Limb** mode—a mode which allows the device to flex to aid donning while seated.
- **Standing Support**—to provide support while standing still.
- **Wheelchair** mode—only operates while seated and when the thigh is horizontal. It provides assistance to support the knee at a chosen angle when the knee is extended. The knee flexion angle where at which the mode will engage is between 45° and 10° (relative to fully extended = 0°).

The device has two manually-selected **User** modes, which you can select using specific key presses (see Section 8) on the **Mode** button on the device:

- **Manual Lock** mode—allows the device to be manually locked against flexion.

Note... this mode is only available if it has been enabled by your clinician.

- **Silent** mode—allows non-critical audio indications to be muted.

Note... the practitioner can also control silent mode from their app if connected to your device, and can adjust the sound and vibration levels and frequencies.

The device provides operating status and warning information directly via:

- Single LED indicator
- The Battery Indicator (LED strip)
- Audio indicators (beeps)
- Haptic indicators (vibrations)

SMARTSTEP—Wearer's App Interface

The Blatchford app **SMARTSTEP** has a wearer-level interface, or "Wearer app", which can provide you with statistical information on usage and other information.

If you would like to, or have been asked to, subscribe to the wearer app, first download the **SMARTSTEP** app to your tablet/phone (see Appendices) and register yourself on the app to apply for access to the wearer app functions.

When you have been granted authorization for wearer-level access, this will allow you to login and open the wearer app and access certain screens and information appropriate to your level of access.

To see live data from your device you will need to connect the app to your Avior via **Bluetooth®**. To do this, switch on **Bluetooth®** on both your tablet/phone and on your Avior, then your app will connect your tablet/phone to your Avior automatically.

When connected, you will be able to use **SMARTSTEP** to see your **SMARTSTEP** Home screen with options to navigate to different information screens and functions.

Status information shown directly on the device such as battery level indication is reproduced on the app when connected via **Bluetooth®** to the device, however device warnings (that would appear as LED light changes on Avior) are not pushed to the app (i.e., warnings, alerts or messages are not sent to overlay the current screen view while connected via **Bluetooth®**). Audible and haptic warnings and indications given by the device are not replicated on the app.

The app is available in two platform variants, one for iOS tablets/phones, one for Android tablets/phones—both variants have near identical screens, controls and levels of authorised access. Layout of information depends on the programming screen aspect ratio.

Contraindications

This device is not suitable for Activity Level 3 or Activity Level 4 users.*

Not suitable for wearers whose weight and/or activity levels fall outside of the recommended levels, i.e., weigh over 125 kg and/or categorized as Activity Level 3 or Activity Level 4 users.

Not to be immersed in water or exposed to salt or chlorinated water.

Not suitable for users who do not have suitable cognitive and physical ability to use and operate the device.

Not suitable for sports use.

***Activity Level**

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Contact your practitioner if your condition changes.

For use only as part of a lower extremity prosthetic limb build.

Device must only be fitted by suitably trained practitioners or technicians.

Be aware of finger trap hazard at all times.

Only power down the device when you are seated.

Immediately stop using the device and report to your clinician, service provider or Blatchford any breakage, damage, crack, leakage of hydraulic fluid or any changes in the performance or function of the device, e.g. restricted movement, non-smooth motion or unusual noises. Only continue use of your device once it has been confirmed safe to do so by your clinician, service provider or Blatchford.

Device contains Lithium-ion battery, do not tamper/puncture.

To prevent irreparable damage to the battery, do not allow exposure to temperatures below -20°C (-4°F), or above $+50^{\circ}\text{C}$ (122°F). Only charge the battery between temperatures of 0°C (32°F) and $+45^{\circ}\text{C}$ (113°F).

Assembly, programming, maintenance and repair of the device (including handling and replacement of the battery unit) must only be carried out by a suitably qualified practitioner or technician that has attended an approved training course.

Only for use in accordance with manufacturers recommendations; not for sports use.

Magnetic charger is solely for use with this device, do not use charger to charge any other device.

Do not place near any heat source. Do not leave in direct sunshine or inside a vehicle in hot weather.

Do not immerse device in water. It may be worn in a freshwater shower (avoid soaping the device). Rinse and wipe dry immediately after showering, or if the device gets wet through any other means.

Avoid exposing the device to abrasive or corrosive environments, such as those containing sand, salt or chemically treated water (e.g., swimming pools). After any exposure to such environments, rinse the device in fresh water, taking care not to exceed the water resistance limitations (see Section 4) and dry thoroughly.

Ensure no water enters the charging socket on the device and the cover is on at all times, except when charging.

Never connect the charger to the device or charge the device while the device is being worn.

The battery charger is not water resistant. Charging the device must be carried out only in a dry environment.

Avoid strong magnetic fields, sources of electrical interference, atmospheres containing liquids and/or powders.

Do not tamper with the mechanical set-up of the prosthesis. If any component screws or nuts become loose, stop using the device straight away and report this to your practitioner who can adjust and recalibrate the device as necessary.

-  Use of accessories, chargers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
-  If the device indicates a fault or critical warning, or goes into *Safety* mode for any reason, stop using the device as soon as is practicable and safe to do so and investigate the nature of the fault/warning. Refer to Section 8.2.8. If you cannot clear the fault, contact your practitioner or service provider for further advice or to arrange a service.
-  Portable RF (radio frequency) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer; this may otherwise result in the degradation of performance of the device.
-  Magnets can impact the function of pacemakers and implantable cardioverter defibrillators (e.g., actuation of reed switch). Keep a minimum distance of 20 cm (8 inches) between the magnetic connectors and the implanted devices to prevent malfunction and danger to health. Contact your practitioner for further advice on whether the device can be used by a user fitted with any other specific electronic medical devices.
-  Inform your practitioner/service provider of any changes in body weight and/or activity level.
-  You must not use the device to operate heavy machinery while wearing the device.
-  Always use a handrail when ascending and descending stairs and at any other time if available.
-  Leg-over-leg stairs descent should only be undertaken if you have received suitable training and have been approved to do so. Otherwise a 'step-to' gait should be used when descending stairs.
-  Ensure any vehicle used has suitable adaptation where necessary for the user's condition. When operating a motor vehicle, all persons are required to observe all applicable driving laws where they are driving. The prosthetic limb must not be used to control the vehicle or its components (e.g., clutch pedal, brake pedal, accelerator).
-  When starting to walk or ascend stairs, you should always lead with your unaffected leg. When descending stairs you should always lead with the prosthetic leg.
-  Avoid, where possible, using the device on slippery, uneven or loose surfaces (sand, snow, mud, cobbles etc.). If unavoidable, proceed with extreme caution.
-  Before wearing footwear with a different heel height, contact your practitioner to arrange the necessary alignment and programming changes.

3 Avior User Controls and Charging Connections


Fig. 1 Avior Component Identification

3.1 AC-DC Battery Charger Adapter

Battery Charger with various plug in blade adapters.

Input 100-240VAC, 50-60Hz.

Output 5-5.2 V DC.

3.1.1 Fitting the Battery Charger Plug Blade

1. Slide down battery charger clip and hold clip open.
2. Insert latch on required plug blade into notch in battery charger unit recess.
3. Push plug blade into battery charger unit recess and slide up clip to close and retain plug blade.

Note....If using the IEC C18 ("kettle" type) plug blade adapter with this battery charger, an additional cable is required (not supplied) that has an IEC C17 socket at one end and an appropriate mains plug at the other. The cable and connectors must comply with the local electrical regulations where used.


Fig. 2 Battery Charger

4 Key Points

This device is a water resistant (IP65) microprocessor-controlled multiaxial knee with a hydraulic control cylinder that automatically controls knee flexion and extension resistances to aid walking standing and sitting.

The device is configured by your practitioner for your individual gait needs and adapts in real time to your current motion situation. It has various modes of operation and features that include adjustable stance flexion and extension, adjustable swing extension, intuitive standing support mode, supported stand-to-sit, ratchet sit-to-stand. Additionally, yielding down slopes and stairs with enhanced support, and **Wheelchair** mode, may be enabled and/or configured for you if deemed suitable by your practitioner.

A **Power** ON/OFF button and a **Mode** selector button on the device provide manual controls (to power the device ON or OFF, to enable/disable certain user modes, to perform a battery check, etc.). Specific button-press durations and combinations (either short, long or very long presses of one or other button, or both buttons together) are used to select/deselect individually each of the available functions/modes. A single LED indicator, a battery indicator, audible sounder and haptic vibrator provide device feedback for operating modes, battery charge status and multi-level alarm conditions (see Section 8).

If a high priority alarm is raised, the device operation switches automatically to **Safety** mode, which limits the motion of the device so that you can still walk, but with a restricted gait, to a safe place to address the fault/remove the device.

Do not tighten or adjust any of the bolts and screws on the device yourself or make any adjustments to the mechanical set-up of the prosthesis after setup and calibration by your clinician as this could adversely affect operation, and may require the device to be recalibrated.

You cannot make any adjustments to the programmed performance of the device, these can only be carried out by your practitioner. However, you can, using the **Mode Selector Switch**, switch ON or OFF the following user modes:

- **Manual Lock** mode (if this function has been enabled for you by your practitioner), and
- **Silent** mode.

The device is powered by a Lithium-Ion battery and should be charged daily using the provided charger. When the battery is fully charged the device will provide up to two days of normal use.

The device functionality addresses the following user activities:

- | | |
|---|--|
| • Walking forwards (all forward directions) | • Sitting-to-standing (ratchet function) |
| • Walking backwards (all backward directions) | • Ability to reposition limb after sitting still for 2 seconds |
| • Slope ascent | • Kneeling |
| • Slope descent | • Donning and doffing* |
| • Stairs descent | • Wheelchair use* |
| • Stairs ascent | • Manual Lock* |
| • Standing | • Getting into a car |
| • Standing-to-sitting | • Getting out of a car |
| • Sitting (and free to reposition knee) | • Stumble recovery |

* If associated user mode or function has been enabled for you.

If you are granted permission to use the wearer app (wearer's interface to **SMARTSTEP**) with your device, you can keep track of your activity and usage. If you are subscribed to the app, your practitioner may also ask you complete on-line tasks, questionnaires, etc., via the app if this functionality is made available to you.

Maintenance

There is no planned service requirement during the 3 year warranty period, unless stated otherwise as a requirement of any extended warranty purchased by your provider.

The device must be inspected by the clinician every year.

Any maintenance must be carried out by competent personnel. Do not carry out maintenance on this device yourself. Instead, return it to your practitioner/service provider for maintenance.

Do not lubricate any components of the device

If a warning/critical warning occurs and you cannot fix the fault directly (e.g. by recharging the battery if the battery empty warning is given), then contact your practitioner who may be able to diagnose the fault and explain what remedial action you can take to clear the fault. If you cannot clear the fault, then the device will need to be returned to the practitioner/service provider for a service or repair.

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase or decrease in knee stiffness
- Instability
- Reduced knee support (free movement)
- Any unrecognized beeps or warning lights
- Any unusual noise*

* Due to the nature of the hydraulics there may be some slight air noise from the cylinder during the first few steps. This is not detrimental to the function of the unit and should dissipate quickly. If symptoms persist, please contact your practitioner.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces; IPA (Isopropyl Alcohol (Isopropanol)) wipes, Clinell wipes and Lysol may also be used to clean outside surfaces. DO NOT use aggressive cleansers.

When rinsing the device, only use low pressure showers or water jets that do not exceed a flow rate of 12.5 litres per minute and pressure of 30 kPa (4.4 psi). Thoroughly drain and dry before use.

To clean the magnetic mating surfaces of the charging connector and charging terminal use a clean and dry cotton swab or a soft bristle brush. Make sure that no cotton fibres or bristles are left on the surfaces.

Never clean or expose the device to moisture while charging.

Never apply any cleaning liquid, detergents or products to the magnetic mating surfaces.

Limitations on Use

Environment

The device is dustproof and water resistant to IP65¹ and can withstand temporary exposure to a low-pressure jet of fresh water with no harmful effects. The device can be worn while showering in fresh water (avoid soaping the device while showering), and it must be rinsed and dried immediately after showering. The device is not to be used for bathing, paddling/walking in water, swimming or prolonged submersion.

Avoid exposing the device to corrosive elements such as salt water, chlorinated water, acids, bleach, aggressive cleansers and detergents (except for approved cleansers for cleaning the device as described in *Cleaning* in Section 4). If exposed, it should be rinsed with fresh water and dried immediately.

Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

The charging surfaces of the connectors should be kept clean and kept away from other magnetic or conductive materials.

Exclusively for use between -10 °C and 50 °C (14 °F and 122 °F).

Suitable for shower use

Water resistant but not suitable for submersion

¹ Water resistance tested from all directions using a water jet from a 6.3 mm (1/4") dia. nozzle at a distance of 3 m (9.8 ft) with fresh water @12.5 litres/minute and 30 kPa (4.4 psi) for a maximum of 3 minutes. Dust tight.

5 Battery Charging

The device has rechargeable Lithium-ion batteries and should be charged daily using the charger supplied with the device (see Section 3.1 for charger description).

Never connect the charger to the device or charge the device while the device is being worn.

Only charge between battery temperatures of 0 °C (32 °F) and +45 °C (113 °F). If, in use or during charging, the device battery temperature exceeds these maximum and minimum temperature limits then a charging error is flagged the device LED will flash yellow (medium priority). If the battery temperature goes above +60 °C (140 °F) then **Safety** mode will be triggered and the LED will flash red (high priority).

Note... The charger itself has no indicator LED.

To fully charge from empty will take a maximum of 6 hours. When fully charged, the device will provide up to two days of normal use.

Ensure the interfacing connections are clean, dry and free from debris.

Note... Connectors are magnetic so may attract iron or steel particles, filings etc.

If charging from a mains wall power outlet only the provided wall charger should be used.

The device can be charged using the provided cable in combination with a CE Marked, 5V USB 1.0 or above compliant portable battery pack or in-car charger.

Blatchford takes no liability for products damaged by faulty or unapproved USB chargers. If your device begins to behave erratically during or after charging, please contact your service provider or Blatchford.

The device's battery status/warning indications in use are described in Section 8.2.6. The device's charging sequence indications (LED lights, sounds and haptic signals), and battery level indications while charging are described in in Section 8.2.5.

6 How to Don and Doff the Device

Donning and doffing should be carried out seated, however it is possible to don or doff from a standing position (with support).

If the limb is powered ON and held horizontally for two seconds while you are seated it will activate **Donning the Limb** mode (if this mode has been enabled by the clinician) and allow the knee to flex easily, to facilitate donning from a seated position.

7 Transportation and Storage

Storage

Ensure that the device is in a straight (fully extended) position when not in use.

If the limb is to be stored upright then it must be secured to prevent it falling over or flexing.

Transportation

If the device is shipped either on its own or built into a limb it should be packaged in a suitable box, ideally with the device fully flexed. For storage below -10°C , store in original packaging with knee fully flexed.

For transport purposes the device must be switched off (using hardware shutdown, see Section 8.2.1) and packed securely so that it is

- a) thoroughly protected from damage in transit, and**
- b) cannot be inadvertently switched on and become operative while in transit.**

8 Device Operation

8.1 Controls and Indicators

Control Panel


Fig. 3 Control Panel

*Your clinician can select the audio volume level for you for non-safety critical audio indications and adjust the pitch of the audio tones and intensity of the haptic responses on their app. However you can mute these non-safety critica indications using **Silent** mode (See Section 9.1.5).

8.2 Device Functions and Responses

8.2.1 Powering the Device ON or OFF

To Power ON Avior:

1. Short press (press for 1 second and release) the **Power** ON/OFF button.


Avior responds with an audible chime (3 short tones rising in pitch) along with 3 vibrations in sync with the audio to acknowledge Avior is powering ON.

At the same time, the LED switches on and illuminates initially steady red, then either:

- Flashes green twice and emits two long beeps and vibrations in sync, or
- After start-up sequence is complete, and power ON is successful, changes to steady green; Avior is ready to use.
- If there is an error:
 - Flashes red.

The device takes up to 10 seconds to power ON and be ready for use, i.e., from the **Power** button press to when the system enters normal operation— either **Walking** mode (indicated by the LED illuminating a steady green), or **Donning the Limb** mode (indicated by the LED illuminating steady purple).

Note*... If device has not yet been programmed the LED will stay steady red on powering ON

To Power OFF Avior (software shutdown):

1. Long press (press for 4 seconds and release) the **Power** ON/OFF button.

Avior responds with an audible chime (3 short tones falling in pitch) to acknowledge power is OFF along with 3 vibrations in sync with the audio.

The LED illuminates initially steady red, then switches off when Avior is shut down.

If pressing for 4 seconds and releasing does not shut down Avior (e.g., if there is a fault with the program preventing shutdown), then carry out a hardware shutdown.

To Power OFF Avior (hardware shutdown):

2. Very long press (more than 8 seconds) the **Power** ON/OFF button. This will force Avior to shut down (a hardware shutdown). The device does not chime in this hardware shutdown.

The LED switches from flashing red then extinguishes when Avior is shut down.

Always use hardware shutdown before storing or transporting the device.

8.2.2 Switching **Bluetooth®** ON/OFF and Pairing with your Phone or Tablet (Option)


If you intend to use the **SMARTSTEP** app to monitor your Avior you will need this information, if not, then ignore this section.

You must also install **SMARTSTEP** on your phone or tablet and register to use **SMARTSTEP** with your Avior, then log in to the app and navigate to the **Bluetooth®** connection screen before attempting to set up the wireless connection between the app and your Avior .

Bluetooth® needs to be enabled on both your tablet/phone and on Avior to connect your tablet/phone to Avior via a **Bluetooth®** Low Energy wireless connection.

To switch ON **Bluetooth®** Low Energy on Avior

1. With Avior powered ON, short press BOTH **Power** and **Mode** buttons simultaneously for 1 second.


2. LED flashes blue while **Bluetooth®** is enabled on the device and is available to pair with a tablet/phone.
3. LED changes to steady blue when paired with the tablet/phone and **Bluetooth®** communication is successfully established.

To switch OFF **Bluetooth®** Low Energy on Avior

1. Short press both **Power** and **Mode** buttons simultaneously.

2. LED changes to previous state i.e., Green if powered ON and no error, (or red flashing if any error).

Note... There are no beeps or vibrations associated with switching **Bluetooth®** ON of OFF on Avior.

8.2.3 Battery Status Check

To check Avior battery status:

- 1. With Avior powered ON, press the **Power** button for 1 second.
- 2. The battery level indicator will illuminate to indicate the battery's current level of charge.
- 3. LED display remains illuminated for up to 10 seconds then extinguishes.

See Section 8.2.4 for battery status indications and their meanings.


8.2.4 Battery Status Indications with Device in Use and Battery Status Check Selected

On selection of a battery status check (see Section 8.2.3) the battery charge status is presented as shown below for up to 10 seconds and then extinguishes.

Battery Condition	Battery Charge Status (% full)	Audio Signal	Haptic Signal	Battery Level Indicator Display	
Full.	87.5–100 %	4 short beeps.	None.	4 LED bars illuminate steady green.	
Very good.	62.5–87.5 %	3 short beeps.	None.	3 LED bars illuminate steady green.	
Good.	37.5–62.5 %	2 short beeps.	None.	2 LED bars illuminate steady green.	
Low.	12.5–37.5 %	1 short beep.	1 short vibration.	1 LED bar illuminates steady green.	
Critically low.	5–12.5 %	1 short beep.	1 short vibration.	1 LED bar illuminates flashing red.	

8.2.5 Battery Charging Notifications while Charging

Note... If the device is powered off, the battery can be charged when connected to the charger, but no charging indications will be displayed. If device is powered ON, when the charger is disconnected, the device will give one long beep to notify the user of disconnection—useful if the device is accidentally disconnected from the charger. This audio function can be disabled if required.

State	Pre-Condition	Audio Signal	Haptic Signal	Battery Level Indicator Visual Signal	
Charging operating correctly.	Device powered on, charger connected.	None.	None.	LEDs illuminate, cycling to the right, then dim after 1 minute.	
Charger disconnected.	Charger connected and charging.	1 long beep (unless configured to mute).	In sync with audio.	Charge status displayed for 10 seconds.	
Error charging.	Charger connected but error charging.	2 short beeps repeated 6 times, pattern repeated every 20 seconds.	In sync with audio.	1 LED bar illuminated steady red.	

8.2.6 Battery Discharging Notifications

While using the device these indications appear without needing to press the power button.

State	Pre-Condition	Audio Signal	Haptic Signal	Battery Level Indicator Visual Signal	
Battery Low Warning. Battery status indications on Power ON or in use, when low/no charge.	Device powered ON. Battery level discharged to 25% charge.	1 set of 3 short beeps.	In sync with audio.	Battery low in use or on power up.	
Critical Warning. Stop using the device and recharge it.	Device powered ON. Battery discharged to 12.5% charge.	2 sets of 3 short beeps.	In sync with audio.	Battery critically low in use or on power up: 12.5% to 5% charge remaining.	
Final Critical Warning. Last warning before device enters <i>Safety</i> mode automatically.	Device powered ON. Battery discharged to <5% charge.	3 sets of 3 short beeps. Nine beeps in total, 12 seconds total duration*.	In sync with audio.	Battery critically low on in use or on power up: 5% to 0% charge remaining.	  LED flashes yellow See <i>Medium Priority Alarm-Battery Warning</i> in Section 8.2.8.

*On hearing the 12 seconds **Final Critical Warning** alarm start, the wearer should stop moving as soon as possible and stay still for at least 4 seconds—this action will initiate *Safety* mode (See Section 8.2.8) on the device to allow the wearer to move with a restricted gait to a safe position and remove the device.

8.2.7 User Modes—Selection and Indications

Note... The LED color purple is used to denote a user mode is active.

Silent Mode

With *Silent* mode active, all audio indications are turned off except for *High Priority* and *Medium Priority* alarms.

Task	Pre-Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
Enter <i>Silent</i> mode.	Device powered ON.	Short press <i>Mode</i> button twice. 	Short beep after each button press, mode selection confirmed with 2 short beeps.	Short vibrations in sync with audio beeps.	 Flashes purple twice, then extinguishes when <i>Silent</i> mode is switched on.
Exit <i>Silent</i> mode.	Device powered ON, in <i>Silent</i> mode.	Long press <i>Mode</i> button twice. 	Short beep after each button press, mode selection confirmed with 2 short beeps.	Long vibrations in sync with audio beeps.	 Flashes purple twice, then either extinguishes, or returns to previous mode indication.

Note... *Silent* Mode can also be switched on or off by your practitioner via their app.

Manual Lock Mode

In **Manual Lock** mode the knee can extend up to full extension but is locked against flexion and cannot flex. This mode cannot be entered if device knee flexion is greater than 45 degrees.

Task	Pre-Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
Enter Manual Lock mode.	Device powered ON.	Short press Mode button 3 times. 	Beep after each button press, mode selection confirmed with 3 short beeps.	Short vibrations in sync with audio beeps.	  Flashes purple 3 times, then illuminates steady purple when locked mode is switched on.
Exit Manual Lock mode.	Device powered ON, in Manual Lock mode.	Long press Mode button once. Release after device response is complete. 	2 long beeps.	2 long vibrations in sync with audio beeps.	 Flashes green twice, then extinguishes.

Wheelchair Mode

Task	Pre-Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
Enter Wheelchair mode.	Device power ON, User seated, thigh near horizontal.	Extend knee to a valid flexion angle and hold briefly to activate Wheelchair mode.	5 × short beeps.	Short vibrations in sync with audio beeps. Knee locks in flexion at chosen position.	  Flashes purple x 5, then illuminates steady purple.
Exit Wheelchair mode.	Device powered ON, in Wheelchair mode.	Shin raised and released, or knee knocked.	2 × long beeps.	Short vibrations in sync with audio beeps.	 Flashes green twice, then returns to previous mode indication.

8.2.8 Alarms and Safety Modes Indications

High Priority Alarm, Safety Mode

The device enters **Safety** mode when a serious event or fault has occurred and a high priority alarm, such as critically low battery warning or an internal device fault warning is raised.

In **Safety** mode the device is set to high flexion resistance but is free in extension to allow you to move with restricted motion carefully to a safe place and try to fix the fault or stop using the device. Carry out a power reset in the event of such a fault arising as this, in some cases, may clear the fault; but if it does not, or the fault reoccurs, then stop using the device as soon as is safe to do so and report the fault to your service provider or practitioner.

Alarm Type	Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
High Priority. Safety Mode.	Critical fault detected or Safety mode activated.	High priority alarm to instruct the user to stop using the device as soon as is practicable. User may need to change gait and/or activity.	Burst of dual-tone (low and high pitch) long beeps, 14 beeps in total, 0.5 second between pulses. At maximum volume.	Vibrations in sync with audio beeps.	 Flashes red until fault cleared or Safety Mode deactivated.

Medium Priority Alarm—System Too Hot or Too Cold

Alarm Type	Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
Medium Priority System too hot or too cold.	System too hot or too cold. Safety mode activated.	Medium priority alarm to instruct the user to limit activity if too hot to prevent overheating, or if too cold to move to warmer environment. User may need to change gait and/or activity.	3 × burst, where burst = 3 × long beeps, with 0.25 second between beeps, and 2.5 seconds between bursts. 9 beeps in total. Approx. 25 seconds duration.	Vibrations in sync with audio beeps.	 Flashes yellow until fault cleared or Safety Mode deactivated.

Medium Priority Alarm—Battery Empty

Alarm Type	Condition	Description	Audio Signal	Haptic Signal	LED Visual Signal
Medium Priority Battery Empty	Battery charge level <5% (empty).	Medium priority alarm to instruct the user that the battery is about to run out. Alarm is active for 12 seconds, then device enters Safety Mode , or earlier if wearer stands still for 4 seconds.	3 × 3 short beeps. Nine beeps in total. Approx. 12 seconds duration.	Vibrations in sync with audio beeps.	 Flashes yellow until device enters Safety mode, or charge is increased above 5%.  Battery Indicator 1 LED bar illuminates flashing red until charge is increased above 5%.
Action	 On hearing the alarm start, the wearer should stop moving and stay still for at least 4 seconds— this action will initiate Safety mode on the device to allow the wearer to move with a restricted gait to a safe position to remove the device.				

Low Priority Alarm

To pause a **Low Priority** alarm, short press the **Mode** button

Note... Only Low Priority alarms can be paused.

Alarm Type	Condition	Action	Audio Signal	Haptic Signal	LED Visual Signal
Low Priority	Service needed.	Alarm to instruct the user to send device for service. User may need to change plans on device use.	1 short beep repeated every 20 seconds up to 12 beeps total. Repeats when powered on or charger removed until service reset.	Vibrations in sync with audio beeps.	 Illuminates steady yellow until service reset, or paused.

9 Program Operation

The device's program automatically adapts the device operation to suit different activities (walking, sitting, standing, walk down stairs, etc.). The program follows operating mode and function settings configured by your practitioner and also any user modes enabled for you, or by you.

Depending on the base configuration chosen by your practitioner, your device will be set either to lock or to yield during the stance phase of walking, and set either to lock or to be free to swing during the swing phase of walking.

The configuration and fine tuning settings programmed by your practitioner will be chosen to best suit you, your condition and stage in your rehabilitation progress, as well as taking into account your ability and confidence in using the device.

The device is capable of adapting automatically to the following user activities and situations:

- Walking forwards
- Walking backwards
- Walking on uneven terrain, camber inside or outside, urban or country terrain
- Standing on a slope (up to 5 degrees)
- Walking up a slope (up to 5 degrees)
- Walking down a slope (up to 5 degrees)
- Ascending stairs ascent using a step-to gait only
- Descending stairs descent using a step-to gait or step-over-step
- Standing stationary
- Standing-to-sitting
- Sitting and adjusting limb position
- Sitting-to-standing
- Kneeling
- Donning and doffing
- Wheelchair use
- Manual Lock required
- Getting into a car
- Getting out of a car
- Stumbling

9.1 Operating Modes: User Modes and Functions

The following **User** modes and functions may be made available to you by your practitioner:

- **Standing Support**
- **Donning the Limb** mode
- **Manual Lock** mode
- **Wheelchair** mode

Note... **Silent** mode is a device function that can be enabled or disabled by both user and practitioner (see Section 9.1.5). **Silent** mode is disabled by default.

These **User** modes and functions, and **Safety** mode, are described in the next sections.

9.1.1 Standing Support

Standing Support mode is a clinician-enabled function. The device recognizes when you are stationary and automatically switches to **Standing Support** mode, if enabled for you. In this mode the device goes into **Lock** to provide support while standing. The sensitivity of the detection system that determines how long you need to be sufficiently stationary for **Standing Support** mode to engage can be adjusted by your clinician.

9.1.2 Donning the Limb Mode

Donning the Limb mode is a clinician-enabled user mode. If enabled for you, **Donning the Limb** mode will initiate after powering ON while you are seated.

After the powering ON start-up sequence, and while seated, hold the device horizontal and fully extended for two seconds. This will engage **Donning the Limb** mode allowing the knee to flex or extend. When the device detects you starting to load the device and beginning to move from sit to stand, it automatically exits **Donning the Limb** mode and enters **Sit-to-Stand** mode.

If **Donning the Limb** mode is not enabled for you, when the device is powered on and detects you are standing or sitting, it behaves normally, i.e., going into **Walk** mode or **Sit-to-Stand** mode respectively.

9.1.3 Manual Lock


Manual Lock mode is a clinician-enabled user mode. In **Manual Lock** mode, the device is locked against flexion but can still be extended until fully extended. Mode selection/de-selection is described in Section 8.2.7.

9.1.4 Wheelchair Mode


Wheelchair mode is a clinician-enabled user mode (for use if the clinician is confident you been trained in **Wheelchair** mode and can use that mode independently).

You must be in a seated position with the thigh almost horizontal before the mode can be activated (if enabled). To then activate the mode, proceed as shown in Fig. 4.


1. Extend the knee to the required position (to a knee flexion angle between 45° and 10° , where 0° degrees is knee fully extended) in one movement using either the sound limb beneath the device or by hand, and then hold it still at this position for one second; the knee will then maintain lock in flexion at this position (so the user then does not need to use the footrest on a wheelchair).

2. To release from **Wheelchair** mode flexion lock, extend the knee slightly, and release.

3. As a safety feature, the device can also be released from **Wheelchair** mode flexion lock by knocking the device against something (i.e., by applying an abrupt force under the foot, or to the front, back or sides of the foot). This will release the flexion lock.

4. On release from **Wheelchair** mode flexion lock, the knee becomes free to flex and the shin/foot will then descend gently until the foot meets the floor.


Fig. 4 **Wheelchair** Mode Function

9.1.5 Silent Mode

Silent mode, when enabled, disables non-critical audio feedback. It can be enabled or disabled by both the user and practitioner using the correct device mode button press sequence (see Section 8.2.7). The clinician may also enable or disable **Silent** mode using their app.

Note... **Silent** mode is disabled by default.

9.2 Safety Mode

In normal operation the device operates automatically, adapting to your movement and situation, whether you are walking, sitting, standing, ascending/descending slopes or stairs etc. However, in the event of any malfunction of the electronic or the hydraulic control system, or when the battery is fully discharged, a **High Priority** alarm (see Section 8.2.8) is raised and the device automatically switches to **Safety** mode.

In **Safety** mode the device is set to high flexion resistance but is free in extension to allow you to move with restricted motion carefully to a safe place and try to fix the fault or stop using the device.

After the malfunction has been resolved, or after a low/empty battery has been recharged, normal operation will resume when the device is next powered on.

9.3 How to Engage Stand-to Sit Mode in Normal Operation

To engage **Stand-to-Sit** mode, the user must stand with both limbs evenly loaded, reach behind for the arms of the chair or seat, move their body posteriorly, with nose over toes (which initially moves/inclines the shin posteriorly) and, if necessary, push forwards with their residual limb (see Fig. 5.)


Fig. 5 Engaging **Stand-to-sit** Mode

10 SMARTSTEP—Wearer App

10.1 App Overview

Note...*SMARTSTEP* is only available to users who have been granted access by Blatchford

SMARTSTEP allows you (when wearer access has been granted by Blatchford), access to your user-relevant data for the device which is assigned to you. It can provide a range of device usage statistics and give access to your account details.

Some options shown here may or may not be available at first software release. Check your app for updates and always update to the most recent software release to ensure you have the latest features and *SMARTSTEP* information.

Your assigned device's technical data is also accessible via *SMARTSTEP* by your authorized clinician.

For regulatory compliance, your personal data is only available to authorized personnel with your permission.

10.1.1 Installation and Registration

For information on installing, registering and logging in to *SMARTSTEP* with *Wearer-level* access, and how to connect your tablet/phone to Avior via *Bluetooth®*, refer to the Appendix sections and summary flow below.

1. When installed, the *SMARTSTEP* app icon will appear on your tablet/phone (see Fig. 6).
2. After successful registration and login to *SMARTSTEP*, follow the appendix instructions to connect your tablet/phone to your Avior via *Bluetooth®*.
3. On establishing successful *Bluetooth®* connection to your device, the *SMARTSTEP Avior Home* screen will open as shown in Fig. 7 .

From this screen you can see a summary of your activity, device connection status and battery level indication.

4. Select *DEVICE VIEW* to open the *Live Monitor* Screen to see summary data (such as total step count). For more detailed data select *VIEW ALL DEVICE DATA* .

Note...You may need to scroll down to see all the information.

5. To see your account information, select the *Account* Icon (see Fig. 9).

Fig. 6 *SMARTSTEP* App Icon


Fig. 7 *SMARTSTEP Home* Screen (when connected to your Avior)


Fig. 8 *SMARTSTEP Home Screen Functions*

Fig. 9 *SMARTSTEP Account Icon*

Connects to wearer account.

Shows connection status.

Indicates Avior battery level status.

Shows progress of data synchronization between device and *SMARTSTEP* cloud data as a percentage of syncing process complete.

DEVICE VIEW opens a screen view of your device data (see Live Monitor Screen).


Fig. 10 *Live Monitor Screen*


Fig. 11 *All Device Data Screen*

Appendix 1 Installing the SMARTSTEP Programming App

A1-1 Android Devices

1. Open the **Google Play Store** app on your Android programming device.
2. Search for the **SMARTSTEP** app in the search bar at the top of the screen.
3. Tap on the app name and verify the information of the distributor is 'Blatchford Products Ltd'.
4. Tap the **Install** button on the app.
This will begin the process of downloading and installing the **SMARTSTEP** app on your programming device.
5. When successfully installed, return to your **Home Screen** -> **App menu** and the **SMARTSTEP** app icon will appear among your other apps (see Fig. A1-1).

Fig. A1-1 **SMARTSTEP** App Icon

System Requirements

Android 11 Red Velvet Cake (API level 30) and above.


A1-2 Apple iOS Devices

1. Open the **App Store** app on your Apple programming device.
2. Tap the search icon at the bottom right of the screen.
3. Search for the **SMARTSTEP** app in the search bar at the top of the screen.
4. Tap on the app name and verify the information of the distributor is 'Blatchford Products Ltd'.
5. Tap the **GET** button on the **SMARTSTEP** app.
This will begin the process of downloading and installing the **SMARTSTEP** app on your programming device.
6. When successfully installed, return to your **Home Screen** -> **App menu** and the **SMARTSTEP** app icon will appear among your other apps (see Fig. A1-1).


System Requirements

iOS 14.0 and above.


Appendix 2 Wearer Registration on SMARTSTEP


1. On the **SMARTSTEP Welcome** screen, select **CREATE AN ACCOUNT**.


2. To register as a wearer, select **Create a wearer account** under **For product wearers**.


3. Enter the e-mail address you regularly use to communicate with Blatchford. Click on **SEND VERIFICATION CODE**.


4. Check the inbox of your e-mail account. Enter the verification code you received and click **VERIFY CODE**.


5. After successful verification of your e-mail address, click **NEXT**.


6. Complete the registration form.

7. Click **CREATE** to finalize your registration.


8. Select the preferred language for the application.

9. Select the region where you live.

10. After reading all app-related warnings and permission requests, click **ACCEPT AND CONTINUE**.

11. The app opens at your **SMARTSTEP** home page and is ready for use.

Appendix 3 SMARTSTEP—General Information

Section A3-1 describes the steps to carry out prior to configuring the device for the wearer, including how to install and connect **SMARTSTEP** to a specific device, how to establish **Bluetooth®** connection and select a specific Avior device to Program.

Section A3-2 describes how to assign a wearer ID to the device.

Note... Assigning a user to a device can be carried out before or after set-up and programming.

Section A3-2 describes the function of the **Menu** icons.

A3-1 Procedures Required Prior to Connect SMARTSTEP app to your Device

If the app has not already been installed on your device and login for wearer-level access provided during your training, then:

1. Download **SMARTSTEP** app to your tablet/phone (see Appendix 1).
2. Tap the **SMARTSTEP** icon in the menu screen of your tablet/phone.

The **SMARTSTEP Welcome** screen (Fig. A3-1) displays.

3. Select **CREATE AN ACCOUNT** and follow screen instructions to request wearer-level access to the app.

When wearer-level access is granted you will be sent login details by e-mail permitting log-in to the app.

4. Switch on **Bluetooth®** on your device.
5. Log in to the app again and enter the required information in the **SMARTSTEP Sign-in** screen (Fig. A3-3) .
6. If app does not connect successfully, check **Bluetooth®** is switched on (active) on both your mobile device (tablet/ phone) and your Avior device, and ensure you are within broadcast range (within 3 m). Then select **Retry** (see Fig. A3-4).
7. The app will open at the **SMARTSTEP (Wearer) Home** screen (Fig. A3-5) connected automatically to your device.

Fig. A3-1 **SMARTSTEP** App Icon


Fig. A3-2 **SMARTSTEP** App Welcome Screen


Fig. A3-3 SMARTSTEP Sign-in Screen


Fig. A3-4 SMARTSTEP Retry Screen


Fig. A3-5 SMARTSTEP (wearer) Home Screen

A3-2 SMARTSTEP Menu Icons

Note... Some icons may not be available for wearer-level access to *SMARTSTEP*.

 Menu	• Menu: Opens the list below of menu icons.
 Home	• Home: Opens the <i>Home</i> screen.
 Documents	• Documents: Opens the <i>Documents</i> screen with the list of accessible documents.
 Device	• Device: Opens the <i>Devices</i> screen with a list of devices that can be connected.
 Access	• Access: Shows a list of devices to which you have access, and provides facility to request access to other devices.
 Contact	• Contact: Opens <i>Contact Information</i> screen.
 About	• About: Presents information '<i>About the app</i>' and other information.
 Settings	• Settings: Opens the <i>Device Programmed Settings</i> screen
 Gait Visualiser	• Gait Visualizer: Opens the <i>Gait Visualizer</i> screen.
 Account	• Account: Displays <i>Account information</i> screen.
 Log out	• Log out: Log out from the app
 Blatchford Web	• Blatchford Web: Link to the Blatchford website

Fig. 12 SMARTSTEP Menu Icons

Manufacturer's Declarations and Legal Information

FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth®** module: FCC ID: QOQ-GM220 IC: BGM220PC22HNA2.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342

800-548-3534 (toll free) | 937-291-3636

customerservice@blatchfordus.com

RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchfordmobility.com

Medical Device

Single Patient – multiple use

Warranty

This device is warranted for 36 months.

Extended warranty of an additional 36 months is available.

A service between 30 and 36 months is required for all extended warranties.

For the maximum 36-month extended warranty, a second service must be carried out between 54 and 60 months.

We reserve the right to void extended warranties if the device is not serviced at the required intervals.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Trademark Acknowledgments

Avior, SMARTSTEP and Blatchford are registered trademarks of Blatchford Products Limited.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under license. Other trademarks and trade names are those of their respective owners.

Clinell is a trademark of Gama Healthcare Ltd.

Lysol is a registered trademark of Reckitt Benckiser.

Apple is a trademark of Apple Inc., registered in the U.S. and other countries.

Android is a trademark of Google LLC.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

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© Blatchford Products Limited 2025. All rights reserved.
Blatchford Europe GmbH