
# Avalon<sup>K2</sup>

## User Guide

AV22L—AV30R  
AV22LD—AV30RD

|    |                        |    |
|----|------------------------|----|
| EN | User Guide             | 2  |
| FR | Guide de l'utilisateur | 8  |
| DE | Bedienungsanleitung    | 14 |
| IT | Guida per il paziente  | 20 |
| ES | Guía del usuario       | 26 |
| NL | Gebruikershandleiding  | 32 |
| PL | Podręcznik użytkownika | 38 |
| PT | Guia do utilizador     | 44 |
| CS | Návod k obsluze        | 50 |

# Contents

|                                          |   |
|------------------------------------------|---|
| Contents .....                           | 2 |
| 1 Description and Intended Purpose ..... | 3 |
| 2 Safety Information .....               | 3 |
| 3 Construction .....                     | 4 |
| 4 Maintenance .....                      | 5 |
| 5 Limitations on Use .....               | 5 |
| 6 Technical Data .....                   | 6 |

---

# 1 Description and Intended Purpose

## Please keep these instructions.

The term *device* is used throughout this document to refer to AvalonK2.

### Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for single user.

The device dorsiflexes after mid stance and remains dorsiflexed at toe off through the swing phase to give greater toe clearance for improved confidence and security. The controlled plantar flexion provided by the device at heel contact improves time to foot flat which enhances stability on uneven surfaces and slopes.

---

## 2 Safety Information

**This warning symbol highlights important safety information which must be followed carefully.**

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user should contact their clinician if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

After continuous use the ankle casing may become hot to the touch.

Avoid exposure to extreme heat and/or cold.

The user must not adjust or tamper with the setup of the device.


Be aware of finger trap hazard at all times.

# 3 Construction

## Principal Parts

- |                                             |                                                    |
|---------------------------------------------|----------------------------------------------------|
| • Hydraulic Body Assembly Including Pyramid | Aluminum, stainless steel, titanium                |
| • Carrier Assembly                          | Aluminum, stainless steel                          |
| • Keel                                      | Thermoplastic composite                            |
| • Keel Attachment Screws                    | Stainless steel                                    |
| • Glide Sock                                | Ultra high molecular weight polyethylene (UHMW PE) |
| • Foot Shell                                | Polyurethane (PU)                                  |

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

We recommend a regular visual check at least annually. Signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise

The practitioner must also be informed of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. Do not use aggressive cleansers.

---

## 5 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Suitable for submersion


Foot products must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, it should be inverted and dried before further use.

Exclusively for use between -15 °C and 50 °C  
(5 °F to 122 °F).

# 6 Technical Data

|                                            |                                                        |
|--------------------------------------------|--------------------------------------------------------|
| Component Weight ( <i>size 26</i> )        | 780 g (1 lb 11 oz)                                     |
| Maximum User Weight                        | 150 kg (330 lb)                                        |
| Activity Level                             | 2                                                      |
| Size Range                                 | 22 to 30 cm                                            |
| Build Height<br>(See diagram below)        | 115 mm                                                 |
| Heel Height                                | 10 mm                                                  |
| Range of ankle movement                    | 6 degrees plantar flexion<br>to 6 degrees dorsiflexion |
| Operating and<br>Storage Temperature Range | -15 °C to 50 °C (5 °F to 122 °F)                       |
| Proximal connection                        | Male Pyramid (Blatchford)                              |

## Fitting Length

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Warranty

This device is warranted for 24 months—glide sock 3 months—foot shell 12 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

## Trademark Acknowledgements

Avalon and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

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© Blatchford Products Limited 2024. All rights reserved.
**Endolite India Ltd.**
Blatchford Europe GmbH