
Avalon^{K2}

Instructions for Use

AV22L—AV30R
AV22LD—AV30RD

EN	Instructions for Use	2
FR	Instructions d'utilisation	15
DE	Gebrauchsanweisung	28
IT	Istruzioni per l'uso	41
ES	Instrucciones de uso	54
NL	Gebruiksaanwijzing	67
PL	Instrukcje użytkowania	80
PT	Instruções de utilização	93
CS	Návod k použití	106

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	5
5 Maintenance	6
6 Limitations on Use	6
7 Bench Alignment.....	7
8 Biomimetic Alignment	8
9 Biomimetic Adjustment.....	9
10 Dynamic Adjustment	10
11 Technical Data	11
12 Assembly Instructions	12
13 Ordering Information	13

1 Description and Intended Purpose

These instructions are for the practitioner.

The term *device* is used throughout these instructions for use to refer to AvalonK2.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

This device dorsiflexes after mid-stance and remains dorsiflexed at toe-off through the swing phase to give greater toe clearance for improved confidence and security. The controlled plantar flexion provided by this device at heel contact improves time to foot flat which enhances stability on uneven surfaces and slopes.

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 2 who may benefit from enhanced stability and an increase in confidence on uneven surfaces and slopes.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in Activity Level 1 who would benefit from the enhanced stability offered by this device, but this decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

Due to the cushioned heel and limited energy return from the foot keel this device is not suitable for Activity Level 3 or 4 individuals. This type of user will be better served by a specially designed prosthesis optimized for their needs.

This device is not suitable for users with poor balance.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

Clinical Benefits

- Improved ground clearance reducing risk of trips and falls
- Increased security for transfemoral wearers
- Additional dorsiflexion range allows more posterior positioning of feet, aiding in sit-to-stand
- Reduced stress and loading rate on the residual limb
- Increased walking speed

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

After continuous use the ankle casing may become hot to the touch.

Avoid exposure to extreme heat and/or cold.

The user must not adjust or tamper with the setup of the device.


Be aware of finger trap hazard at all times.

3 Construction

Principal Parts

• Hydraulic Body Assembly Including Pyramid	Aluminum, stainless steel, titanium
• Carrier Assembly	Aluminum, stainless steel
• Keel	Thermoplastic composite
• Keel Attachment Screws	Stainless steel
• Glide Sock	Ultra high molecular weight polyethylene (UHMW PE)
• Foot Shell	Polyurethane (PU)

4 Function

This device comprises a hydraulic body with a valve that can be adjusted to increase or reduce hydraulic resistance to plantar flexion and dorsiflexion simultaneously. The hydraulic body is connected to a carrier assembly via two pivot pins. A keel is attached to the carrier assembly using stainless steel screws. The foot is wrapped in a UHMW PE sock which is in turn surrounded by a PU foot shell.

5 Maintenance

Maintenance must be carried out by competent personnel.

We recommend the following annual maintenance:

- Visually inspect the foot shell and glide sock, check for damage or wear and replace if necessary.

There are no other serviceable parts in the foot assembly.

The user should be advised:

Any changes in performance of this device must be reported to the practitioner

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Foot products must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, it should be inverted and dried before further use.


Suitable for submersion

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

It is recommended that only Blatchford products be used in conjunction with the device.

7 Bench Alignment


Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

8 Biomimetic Alignment

The aim of alignment is to achieve a balance point while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved.

Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.


- * Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

9 Biomimetic Adjustment

Carry out static alignment while ensuring the user has some means of support such as parallel bars. Check the bench alignment and the heel height.

Note... This is standing alignment only.


The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

10 Dynamic Adjustment

The user should experience the ankle moving with the body through the gait cycle. There should be no effort exerted by the user to overcome the hydraulic resistance of the ankle.

Resistance Adjustment Procedure


Guidance

Following dynamic adjustment, trial the foot/ankle on ramps and stairs. Ensure the user is comfortable with the kind of terrain he/she may normally be expected to encounter. If the user reports any issues with comfort, usability or range of movement of the ankle, adjust accordingly. Caution: decreasing plantar flexion resistance may cause knee instability for transfemoral users.

11 Technical Data

Component Weight (size 26)	780 g (1 lb 11 oz)
Maximum User Weight	150 kg (330 lb)
Activity Level	2
Size Range	22 to 30 cm
Build Height (See diagram below)	115 mm
Heel Height	10 mm
Range of ankle movement	6 degrees plantar flexion to 6 degrees dorsiflexion
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)
Proximal connection	Male Pyramid (Blatchford)

Fitting Length


12 Assembly Instructions

 **Be aware of finger trap hazard at all times.**

 **Use appropriate health and safety equipment at all times including extraction facilities.**


1

Place foot in glide sock and slide assembly into foot shell.

2

Ensure the heel of the keel is fully engaged in the location feature in the foot shell.

Cutaway section through foot shell

Heel engagement

3

Abrade the top surface of the foot shell. Bond 12 mm Pelite interface to foot shell and shape as appropriate. Use Evostik 528 (926213) or equivalent.

12 mm Pelite

4

Bond foam cosmesis to top surface of the Pelite using Thixofix adhesive (926204), Evostik 528 (926213) or equivalent and shape to suit.

13 Ordering Information

Order Example

AV	25	L
	Size	Side

e.g. AV25L

Available from size 22 to size 30:

AV22L to AV30R

AV22LD to AV30RD

(add 'D' for a dark tone foot shell)

Foot Shell (for dark add 'D')

Small	Medium	Large	Extra Large
22L 539020	24L 539024	26L 539028	28L 539032
22R 539021	24R 539025	26R 539029	28R 539033
23L 539022	25L 539026	27L 539030	29L 539034
23R 539023	25R 539027	27R 539031	29R 539035
			30L 539036
			30R 539037

Glide Sock, Sizes 22–30

Allen 2.5 (Hex key, valve adjuster tool)

Part No. 405815

Part No. 940234

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months—glide sock 3 months—foot shell 12 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Avalon and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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© Blatchford Products Limited 2024. All rights reserved.
Endolite India Ltd.
Blatchford Europe GmbH