# Silcare Breathe Walk Cushion Liner

# Silcare Breathe Walk Locking Liner

## User Guide

SBWTTCP22 - SBWTTCP40  
SBWTTLP22 - SBWTTLP40

|    |                        |    |
|----|------------------------|----|
| EN | User Guide             | 2  |
| FR | Guide de l'utilisateur | 14 |
| DE | Bedienungsanleitung    | 26 |
| IT | Guida per il paziente  | 38 |
| ES | Guía del usuario       | 50 |
| NL | Gebruikershandleiding  | 62 |
| PL | Podręcznik użytkownika | 74 |
| PT | Guia do utilizador     | 86 |
| CS | Návod k obsluze        | 98 |

# Contents

- Contents .....2
- 1 Description and Intended Purpose .....3
- 2 Safety Information .....4
- 3 Construction .....5
- 4 Maintenance .....6
  - 4.1 Cleaning the Device .....7
  - 4.2 Cleaning the Valve (Locking Liner Only) .....7
  - 4.3 Cleaning the Residual Limb.....7
- 5 Limitations on Use .....8
- 6 Donning the Device .....8
- 7 Fitting the Valve .....9
- 8 Fitting the Spacer.....9
- 9 Fitting Advice.....10
  - 9.1 Cushion Liner.....10
  - 9.2 Locking Liner .....10
- 10 Technical Data .....11

---

# 1 Description and Intended Purpose

These instructions are for the user. Please keep these instructions.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Walk Cushion and Locking Liner unless otherwise stated.

Make sure that you understand all instructions for use, drawing particular attention to all maintenance and safety information sections.

## **Application**

This device is an interface component for use only as part of a lower limb prosthesis.

Intended for single user.

It is intended for a low to moderate activity user. If used on more active users device life may be compromised.

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

## **Cushion Liner**

The device is a biocompatible perforated, cushioned socket interface which allows moisture to pass through the perforations keeping the skin dry.

## **Locking Liner**

This device is a biocompatible, perforated locking liner that allows air and moisture to escape through its perforations. During the gait cycle, a silicone valve provides vacuum suspension by controlling airflow through the perforations in the distal cap.

## 2 Safety Information

This warning symbol highlights important safety information.

Inspect the device for damage/deterioration before and after each use.  
Report anything that might affect functionality.

Any deterioration in residuum condition or any change in sensation should be reported to the practitioner.

Make sure that any damaged skin or open wounds are properly and suitably dressed to prevent direct contact with the device.

Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.

For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.

Enlarged perforations can trap the skin and cause tissue damage. If the perforations enlarge, stop using the device

If distal oedematous swellings corresponding to distal perforations in the liner occur, use of the liner should be discontinued and the swellings reported to the practitioner.

Do not use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

Do not pull or stretch the device. Fingernails, sharp jewelry, and the locking pin can tear the device. If the device is torn, stop using it and contact a Blatchford sales representative.

Sockets with sharp proximal edges can tear the device.

Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.

When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.

To avoid danger of suffocation, keep the device away from babies and children.

Keep the device away from direct heat sources.

Only use the device in combination with corrosion resistant components.


**Do not** overtighten the locking pin.

# 3 Construction


## Principal Parts

- Fabric (polyamide and lycra)
- Main Body (silicone)
- Distal Cap (silicone)
- Valve (silicone)
- Washer (nylon)
- Umbrella (nylon)
- Spacer\* (silicone)

## Cushion Liner

## Locking Liner


\* For use with some locks.  
(See Section 8 *Fitting the Spacer.*)

## 4 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Distal oedematous swellings (see diagram below)
- Enlarged perforations
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device
  - Moisture in the distal end of the device
  - Loss of vacuum

Note... Perspiration discolors some socket materials.

## 4.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

**Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.**

**Dry thoroughly before use.**

**Do not tumble dry.**

### Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally.
4. Clean the valve and carefully clean the distal cap. (Locking liner only.)
5. Fill the device with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
6. Rinse the device with clean water to remove all residues.
7. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the device can stretch and become distorted.

### Machine Wash

Suitable for machine washing at 30 °C.

**Do not invert the device for machine washing.**

## 4.2 Cleaning the Valve (Locking Liner Only)

Make sure the distal end perforations are not clogged and that there is no contamination/debris trapped under the valve.

1. During hand wash squeeze water through the distal perforations.
2. Carefully, slightly lift valve and wipe clean underneath with a clean piece of cloth.
3. Check correct operation/function of valve.

## 4.3 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

**Any deterioration in residuum condition should be reported to the practitioner.**

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Dry residual limb.
4. Apply lotion to dry skin, as recommended by a practitioner.

**Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.**

## 5 Limitations on Use

### Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and stretch or distort.

Keep away from sharp objects (such as jewelry, fingernails).

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for shower use

## 6 Donning the Device





### Before you begin

Fit the valve to the locking liner. (See Section 7 *Fitting the Valve*.)

 **Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.**

 **Do NOT pull or stretch the device.**

1. Invert the device so that its silicone side points externally.
2. Align the locking pin with the long axis of the residual limb. (Locking liner only.)
3. Roll the device on the residual limb while releasing any trapped air.





If the user feels numbness, tingling or any unusual sensation while using the device remove it until normal sensation returns.

If this persists, contact your practitioner.

## 7 Fitting the Valve


**⚠ Do not overtighten the locking pin.**

We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.

## 8 Fitting the Spacer

Fit the spacer between the pin lock and the valve only if the pin lock prevents the valve from opening.

## 9 Fitting Advice

### 9.1 Cushion Liner

#### 9.1.1 Moisture Collects Inside the Distal end of the Device

| Cause                         | Solution                                                         |
|-------------------------------|------------------------------------------------------------------|
| The perforations are clogged. | Clean the device. (See Section 4.1 <i>Cleaning the Device</i> .) |

#### 9.1.2 Any Discoloration is Noticed on the Residuum

| Cause                       | Solution                   |
|-----------------------------|----------------------------|
| Incorrectly fitting device. | Contact your practitioner. |

#### 9.1.3 Oedematous Swellings Corresponding to Distal Perforations

| Cause                     | Solution                                                                 |
|---------------------------|--------------------------------------------------------------------------|
| Excessive distal contact. | Reduce end contact by addition of extra socks or lengthen/remake socket. |

### 9.2 Locking Liner

#### 9.2.1 Moisture in the Distal End of the Device

If moisture collects inside the distal end of the device, refer to the following table:

| Cause                                                                          | Solution                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The distal perforations are clogged.                                           | Clean the device. (See Section 4 <i>Maintenance</i> .)                                                                                                                                                                                    |
| The socket does not have enough space in its distal end for the valve to open. | Put the spacer between the pin lock and the valve. (See Section 8 <i>Fitting the Spacer</i> .)<br>If unsuccessful remake socket ensuring sufficient clearance for the valve to operate. (See Section 7 <i>Fitting the Valve</i> - box 4.) |
| The internal shape of the lock-body is stopping the valve from opening.        | Put the spacer between the pin lock and the valve. (See Section 8 <i>Fitting the Spacer</i> .)                                                                                                                                            |

#### 9.2.2 Loss of Vacuum

If the device loses vacuum, refer to the following table:

| Cause                                                                   | Solution                                                                            |
|-------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| The valve is damaged.                                                   | Stop using the device, and contact a Blatchford sales representative.               |
| Debris beneath the valve preventing correct seal.                       | Clean the valve. (See Section 4.2 <i>Cleaning the Valve (Locking liner only)</i> .) |
| The locking pin is not aligned with the long axis of the residual limb. | Re-align the locking pin with the long axis of the residual limb.                   |

### 9.2.3 Separation of the Valve from the Distal Cap

If the valve falls off the distal cap, refer to the following table:

| Cause                              | Solution                                                                                               |
|------------------------------------|--------------------------------------------------------------------------------------------------------|
| The valve is not fitted correctly. | Make sure that the washer is between the pin and the valve. (See Section 7 <i>Fitting the Valve</i> .) |

Note... Do not over-tighten the locking pin/lock rod.

### 9.2.4 Oedematous Swellings Corresponding to Distal Perforations


| Cause                     | Solution                                                                  |
|---------------------------|---------------------------------------------------------------------------|
| Excessive distal contact. | Reduce end pressure by addition of extra socks or lengthen/remake socket. |

## 10 Technical Data

|                                             |                                              |
|---------------------------------------------|----------------------------------------------|
| Principal Materials                         | polyamide, lycra, silicone                   |
| Shore Hardness                              | 30-35 Shore 00                               |
| Component Weight ( <i>size 28</i> )         | 695 g (1 lb 8 oz)                            |
| Activity Level                              | Low to Moderate                              |
| Size Range                                  | 22-40 cm                                     |
| Length<br>(See <i>Dimensions</i> )          | 435 mm                                       |
| Internal Length<br>(See <i>Dimensions</i> ) | 420 mm                                       |
| Operating and<br>Storage Temperature Range  | -15 °C to 50 °C (5 °F to 122 °F)             |
| Locking Liner Only                          |                                              |
| Matrix Length                               | Approx. 10 cm                                |
| Distal End Attachment                       | M10<br>Requires Locking Pin With a Shoulder* |
| Locking Pin Shoulder Diameter               | 13.5-19 mm                                   |

\*Locking Pin Not Supplied

Dimensions

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Warranty

The device is warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

## Trademark Acknowledgments

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

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© Blatchford Products Limited 2024. All rights reserved.
**Endolite India Ltd.**
Blatchford Europe GmbH