Silcare Breathe Walk Cushion Liner

Silcare Breathe Walk Locking Liner

Instructions for Use

SBWTTCP22 - SBWTTCP40

SBWTTLP22 - SBWTTLP40

EN	Instructions for Use	2
FR	Instructions d'utilisation	18
DE	Gebrauchsanweisung	34
IT	Istruzioni per l'uso	50
ES	Instrucciones de uso	66
NL	Gebruiksaanwijzing	82
PL	Instrukcje użytkowania	98
PT	Instruções de utilização	114
CS	Návod k použití	130

Contents

- Contents2
- 1 Description and Intended Purpose3
- 2 Safety Information5
- 3 Construction6
- 4 Function.....7
- 5 Maintenance7
 - 5.1 Cleaning the Device8
 - 5.2 Cleaning the Valve (Locking Liner Only)8
 - 5.3 Cleaning the Residual Limb.....8
- 6 Limitations on Use9
- 7 Choosing the Correct Size.....9
- 8 Casting/Scanning a Locking Liner..... 10
- 9 Casting/Scanning a Cushion Liner..... 11
- 10 Trimming the Device..... 12
- 11 Donning the Device 12
- 12 Fitting the Valve 13
- 13 Fitting the Spacer 13
- 14 Fitting Advice..... 14
 - 14.1 Cushion Liner 14
 - 14.2 Locking Liner 14
- 15 Technical Data 15
- 16 Ordering Information 16

1 Description and Intended Purpose

These instructions are for use by the practitioner.

Please read these instructions carefully before fitting the device.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Walk Cushion and Silcare Breathe Walk Locking Liner unless otherwise stated.

Ensure that the user understands all instructions for use, drawing particular attention to the maintenance and safety information.

Application

This device is an interface component for use only as part of a lower limb prosthesis.

Intended for single user.

The device is manufactured from biocompatible materials. The cushion liner provides a perforated cushioned socket interface which allows moisture to pass through the perforations keeping the skin dry. In addition to this, the locking liner's silicone valve provides vacuum suspension during the gait cycle by controlling the airflow through the perforations on the distal cap.

To optimize comfort ensure that the user is instructed in the correct way to handle and don/doff the device. See Section 11 *Donning the Device*.

Also ensure that they are made aware of how to maintain the device and to keep it clean for maintaining hygiene as set out in these instructions. See Section 5 *Maintenance*.

Activity Level

The device is recommended for Activity Level 1 and 2 users, although may be suitable for use at all four activity levels. If used on more active users device life may be compromised.

The silicone used in the Silcare Walk liners is softer and more flexible than that used in Silcare Active liners. The softer silicone offers more comfort and flexibility in all directions but may feel less connected or secure at higher activity levels. Silcare Walk liners are therefore generally recommended for the less active users (Activity Level 1 and 2).

The silicone used in Silcare Active liners is firmer and slightly less flexible, offering a more secure interface at higher activity levels. Silcare Active liners are therefore generally recommended for the more active users (Activity Level 3 and 4).

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications - Cushion Liner

- Conical residual limbs:
 - This device may not be optimally contoured to fit conical residual limbs
- Users whose suspension sleeve causes a bulk around their knee may be better suited with an alternative type liner
- Users with poor hand or cognitive function might find donning and cleaning difficult
- Poor hygiene

Contraindications - Locking Liner

- Conical residual limbs:
 - Can cause failure of the vacuum suspension
 - Reduce skin contact with the silicone inside the device and make a less effective seal
- Deep scars distally:
 - Can cause failure of the vacuum suspension
 - Reduce skin contact with the silicone inside the device and make a less effective seal
- Short residual limbs: These can cause failure of the silicone side of the device. After donning, if the unperforated fabric is on or above the patella tendon, the silicone is subject to more stress during knee flexion
- Users with poor hand or cognitive function might find donning and cleaning difficult.
- Poor hygiene.

Clinical Benefits

- Provides cushioning for the residual limb in the socket
- Distributes in-socket pressure more evenly, compared to other materials and alternative cushioning solutions
- Improvements in residual limb health and wound healing, compared to non-perforated liners
- Improved heat dissipation, compared to other temperature regulation solutions
- Removes sweat from skin-liner interface
- Patients reported a preference for their perforated liners, compared to non-perforated liners
- Reduces the need to remove prosthesis throughout the day to dry residual limb, compared to non-perforated liners
- The locking liner provides a means of suspension

2 Safety Information

This warning symbol highlights important safety information.

The user should be advised to contact their practitioner if their condition changes.

Any deterioration in residuum condition or any change in sensation should be reported to the practitioner.

Make sure that any damaged skin or open wounds are properly and suitably dressed to prevent direct contact with the device.

Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.

For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.

Enlarged perforations can cause tissue damage. If the perforations enlarge, stop using the device.

If distal oedematous swellings corresponding to distal perforations in the liner occur, use of the liner should be discontinued and the swellings reported to the practitioner.

Do not use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

Do not pull or stretch the device. Fingernails, sharp jewelry, and the locking pin can tear the device. If the device is torn, stop using it and contact a Blatchford sales representative

Sockets with sharp proximal edges can damage the device.

Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.

When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.

To avoid danger of suffocation, keep the device away from babies and children.

Keep the device away from direct heat sources.

Only use the device in combination with corrosion resistant components.


Do not overtighten the locking pin.

3 Construction


Principal Parts

- Fabric (polyamide and lycra)
- Main Body (silicone)
- Distal Cap (silicone)
- Valve (silicone)
- Washer (nylon)
- Umbrella (nylon)
- Spacer** (silicone)
- Retaining Screw* (nylon)
- Casting Dummy (silicone)

Cushion Liner

Locking Liner

938413

* For transportation only. **Do not** use for fitting.

** For use with some locks. (See Section 13 *Fitting the Spacer.*)

4 Function

Cushion

The device, used in connection with an air-tight suspension sleeve and corrosion resistant components, provides control, secure connection and cushioning between the limb and the socket. The perforations and the corrosion resistant expulsion valve of the device allow moisture to escape during stance.

Locking

The device provides control, secure connection and cushioning between the limb and the socket. The perforations of the device allow moisture to escape during stance.

During stance, the valve opens to allow moisture and air to escape through the perforations in the distal cap. During swing, the valve closes and stops air from re-entering the device, thus creating a vacuum, which both improves proprioception and reduces pistoning.

5 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Distal oedematous swellings
- Enlarged perforations
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device:
 - Moisture in the distal end of the device
 - Loss of vacuum

Note... Perspiration discolors some socket materials.

5.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.

Dry thoroughly before use.

Do not tumble dry.

Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally.
4. Clean the valve and carefully clean the distal cap. (Locking liner only.)
5. Fill the device with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
6. Rinse the device with clean water to remove all residues.
7. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the device can stretch and become distorted.

Machine Wash

Suitable for machine washing at 30 °C.

Do not invert the device for machine washing.

5.2 Cleaning the Valve (Locking Liner Only)

Make sure the distal end perforations are not clogged and that there is no contamination/debris trapped under the valve.

1. During hand wash squeeze water through the distal perforations.
2. Carefully, slightly lift valve and wipe clean underneath with a clean piece of cloth.
3. Check correct operation/function of valve.

5.3 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

Any deterioration in residuum condition should be reported to the practitioner.

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Dry residual limb.
4. Apply lotion to dry skin, as recommended by a practitioner.

Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and stretch or distort.

Keep away from sharp objects (such as jewelry, fingernails).


Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

7 Choosing the Correct Size

It is important that the correct size is specified to ensure proper fit and comfort in use.

1. Measure the circumference of the residual limb 4 cm up from its distal end.
2. If the circumference of the residual limb matches one of the sizes in the table, choose the next size down.
3. If the circumference of the residual limb is between one of the sizes in the table, choose the smaller size.

Circumference Measured (cm)	22.2-23.5	23.6-25	25.1-26.5	26.6-28	28.1-30	30.1-32	32.1-34	34.1-36	36.1-40	40.1-42
Available Sizes	22	23.5	25	26.5	28	30	32	34	36	40
Valve Size (Locking Liner)	Small			Medium		Large			Extra Large	

8 Casting/Scanning a Locking Liner

Before you begin

Allow the user to wear the device for 10 minutes.

(See Section 11 *Donning the Device*.)

If a scanning method is used


Fit/trim the casting dummy to the size of the distal end of the liner, see 938413.

To reduce the amount of distal pressure within the final socket, the model should be lengthened by 10-12 mm (depending upon redundant tissue etc.).

If a casting method is used

Fit/trim the casting dummy to the size of the distal end of the liner, see 938413.

We recommend using a casting method that accentuates weight bearing areas such as the medial condylar flare in conjunction with applying tension to the distal pin; otherwise excessive distal contact may result in the final socket, giving rise to oedematous swellings corresponding to the distal perforations as shown below.

Additionally, when fabricating a diagnostic socket, we recommend a 5 mm pelite spacer is used to further extend/lengthen the socket. When you are satisfied with the fit of the socket (by adding any necessary socks) the spacer can, if necessary, be removed if excess distal contact gives rise to distal oedematous swellings.

 **The softness of the silicone can mask the amount of distal contact that is occurring, this will only be evident when oedematous swellings are observed.**

 **The use of a diagnostic socket is strongly recommended to assess the amount of distal pressure/contact that is occurring.**

9 Casting/Scanning a Cushion Liner


Before you begin

Allow the user to wear the device for 10 minutes.

(See Section 11 *Donning the Device*.)

Casting/Scanning

We recommend using a casting method that accentuates weight bearing areas such as the medial condylar flare and displaces redundant tissue distally; otherwise excessive distal contact may result in the final socket, giving rise to oedematous swellings corresponding to the distal perforations as shown below.

Additionally, when fabricating a diagnostic socket, we recommend a 5 mm pelite spacer is used to further extend/lengthen the socket. When you are satisfied with the fit of the socket (by adding any necessary socks) the spacer can, if necessary, be removed if excess distal contact gives rise to distal oedematous swellings.

⚠ The softness of the silicone can mask the amount of distal contact that is occurring, this will only be evident when oedematous swellings are observed.

⚠ The use of a diagnostic socket is strongly recommended to assess the amount of distal pressure/contact that is occurring.

10 Trimming the Device

The device may be trimmed to suit as required but never below the socket trim line.

We recommend using a rotating disc cutter to achieve a bevelled edge.

Avoid trimming through the perforations.

Trim the proximal edge of the device to suit the shape of the user's socket.

If necessary, trim the proximal edge of the device in a wave pattern to reduce shear stress.

Where possible try to avoid abrupt changes of contour and sharp edges which could cause tears in the silicone and fabric.

11 Donning the Device

Before you begin

Fit the valve to the locking liner. (See Section 12

Fitting the Valve.)





Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.

Do NOT pull or stretch the device.

1. Invert the device so that its silicone side points externally.
2. Align the locking pin with the long axis of the residual limb. (Locking liner only.)
3. Roll the device on the residual limb while releasing any trapped air.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness, tingling or any unusual sensation within 10 minutes, do the following:


4. Doff the device, and wait until normal sensation returns.
5. Don the device again.
6. If the user feels numbness, tingling or any unusual sensation again, stop using the device.

12 Fitting the Valve

 **Do not overtighten the locking pin.**

We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.


1

Valve

Put the valve on the distal cap.

2


Washer

Locking Pin
Shoulder

Loctite 243


Put the washer between the locking pin and the valve. Apply Loctite 243 to the threads of the locking pin.

3

Torque the locking pin to 3 Nm, or tighten it with your fingers by 1/8 to 1/4 of a turn.


4

Confirm sufficient clearance in socket for valve to operate correctly.

13 Fitting the Spacer

Fit the spacer between the lock and the valve only if the lock prevents the valve from opening.

14 Fitting Advice

14.1 Cushion Liner

14.1.1 Moisture Collects Inside the Distal end of the Device

Cause	Solution
The perforations are clogged.	Clean the device. (See Section 5.1 <i>Cleaning the Device</i> .)

14.1.2 Any Discoloration is Noticed on the Residuum

Cause	Solution
Incorrectly fitting device.	Contact your practitioner.

14.1.3 Oedematous Swellings Corresponding to Distal Perforations

Cause	Solution
Excessive distal contact.	Reduce end contact by addition of extra socks or lengthen/remake socket.

14.2 Locking Liner

14.2.1 Moisture in the Distal End of the Device

If moisture collects inside the distal end of the device, refer to the following table:

Cause	Solution
The distal perforations are clogged.	Clean the device. (See Section 5 <i>Maintenance</i> .)
The socket does not have enough space in its distal end for the valve to open.	Put the spacer between the pin lock and the valve. (See Section 13 <i>Fitting the Spacer</i> .) If unsuccessful remake socket ensuring sufficient clearance for the valve to operate (See Section 12 <i>Fitting the Valve</i> - box 4.)
The internal shape of the lock-body is stopping the valve from opening.	Put the spacer between the pin lock and the valve. (See Section 13 <i>Fitting the Spacer</i> .)

14.2.2 Loss of Vacuum

If the device loses vacuum, refer to the following table:

Cause	Solution
The valve is damaged.	Stop using the device, and contact a Blatchford sales representative.
Debris beneath the valve preventing correct seal.	Clean the valve. (See Section 5.2 <i>Cleaning the Valve (Locking Liner Only)</i> .)
The locking pin is not aligned with the long axis of the residual limb.	Re-align the locking pin with the long axis of the residual limb.

14.2.3 Separation of the Valve from the Distal Cap

If the valve falls off the distal cap, refer to the following table:

Cause	Solution
The valve is not fitted correctly.	Make sure that the washer is between the pin and the valve. (See Section 12 <i>Fitting the Valve</i> .)

Note... Do not over-tighten the locking pin/lock rod.

14.2.4 Oedematous Swellings Corresponding to Distal Perforations


Cause	Solution
Excessive distal contact.	Reduce end pressure by addition of extra socks or lengthen/remake socket.

15 Technical Data

Principal Materials	polyamide, lycra, silicone
Shore Hardness	30-35 Shore 00
Component Weight (size 28)	695 g (1 lb 8 oz)
Activity Level	Low to Moderate
Size Range (see Section 7 <i>Choosing the Correct Size</i>)	22-40 cm
Length (See <i>Dimensions</i>)	435 mm
Internal Length (See <i>Dimensions</i>)	420 mm
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)
Locking Liner Only	
Matrix Length	Approx. 10 cm
Distal End Attachment	M10 Requires Locking Pin With a Shoulder*
Locking Pin Shoulder Diameter	13.5-19 mm

*Locking Pin Not Supplied

Dimensions

16 Ordering Information

Order Example

SB	W	TT	C	P	25
Silcare Breathe	Walk	Transtibial	C=Cushion L=Locking	Parallel	Size* (22-40)

Available from size 22 to size 40*:
SBWTTCP22 to SBWTTCP40
or SBWTTLP22 to SBWTTLP40

*SBWTTCP23 is for size 23.5. SBWTTCP26 is for size 26.5.

Part	Part Number
Small Valve Kit	559015
Medium Valve Kit	559016
Large Valve Kit	559017
Extra Large Valve Kit	559018
Casting Dummy	559019

Note... Spare parts in this table are to be used with locking liners only.

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: **www.blatchford.co.uk**

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

The device is warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgments

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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© Blatchford Products Limited 2024. All rights reserved.
Endolite India Ltd.
Blatchford Europe GmbH