Silcare Breathe Active Locking Liner

User Guide

SBATTLP22-SBATTLP40

EN	User Guide	2
DA	Brugervejledning	12
NO	Brukerhåndbok	22
FI	Käyttöopas	32
SV	Användarhandbok	42
RU	Руководство пользователя	52
ZH	用户指南	62
AR	دليل المستخدم	72

Contents2

1 Description and Intended Purpose3

2 Safety Information4

3 Construction5

4 Maintenance5

 4.1 Cleaning the device.....6

 4.2 Cleaning the Valve6

 4.3 Cleaning the Residual Limb.....6

5 Limitations on Use7

6 Donning the Device7

7 Fitting the Valve8

8 Fitting the Spacer8

9 Fitting Advice.....9

 9.1 Moisture in the Distal End of the Device9

 9.2 Loss of Vacuum9

 9.3 Separation of the Valve from the Distal Cap.....9

 9.4 Oedematous Swellings Corresponding to Distal Perforations9

10 Technical Data 10

1 Description and Intended Purpose

These instructions are for the user. Please keep these instructions.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Active Locking Liner.

Make sure that you understand all instructions for use, drawing particular attention to all maintenance and safety information sections.

Application

This device is an interface component for use only as a component of a lower limb prosthesis.

Intended for single user.

It is intended for moderate to high activity users with well covered residual limbs.

This device is a biocompatible, perforated locking liner that allows air and moisture to escape through its perforations. During the gait cycle, a silicone valve provides vacuum suspension by controlling airflow through the perforations on the distal cap.

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

2 Safety Information

This warning symbol highlights important safety information.

Inspect the device for damage/deterioration before and after each use.
Report anything that might affect functionality.

Any deterioration in residuum condition or any change in sensation should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.

Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.

For other medical conditions, please, follow the advice and recommendation of a physician or a medical practitioner regarding skin care.

Enlarged perforations can trap the skin and cause blisters. If the perforations enlarge, stop using the device

If distal oedematous swellings corresponding to distal perforations in the liner occur, use of the liner should be discontinued and the swellings reported to the practitioner.

Do not use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

Do not pull or stretch the fabric. Fingernails, sharp jewelry, and the locking pin can tear the device. If the fabric is torn, stop using it and contact a Blatchford sales representative.

Sockets with sharp proximal edges can tear the device.

Take care when handling the Device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.

When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.

To avoid danger of suffocation, keep the device away from babies and children.

Keep the device away from direct heat sources.

Only use the device in combination with corrosion resistant components.


Do not overtighten the locking pin.

3 Construction

Principal Parts

- Fabric (polyamide and lycra)
- Valve (silicone)
- Washer (nylon)
- Umbrella (nylon)
- Spacer* (silicone)


* For use with some locks.
(See *Fitting the Spacer* section.)

4 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Distal oedematous swellings (see diagram below)
- Enlarged perforations
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device, for example:
 - Moisture in the distal end of the device.
 - Loss of vacuum.

Note... Perspiration discolours some socket materials.

4.1 Cleaning the device

Wash the inside of the device daily to avoid any build up of bacteria.

Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.

Dry thoroughly before use.

Do not tumble dry.

Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally. Clean the valve and carefully clean the distal cap.
4. Fill the device with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
5. Rinse the device with clean water to remove all residues.
6. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the fabric can stretch and become distorted.

Machine Wash

Suitable for machine washing at 30 °C.

Do not invert the device for machine washing.

4.2 Cleaning the Valve

Make sure the distal end perforations are not clogged and that there is no contamination/debris is trapped under the valve.

1. During hand wash squeeze water through the distal perforations.
2. Carefully, slightly lift valve and wipe clean underneath with a clean piece of cloth.
3. Check correct operation/function of valve.

4.3 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin, as recommended by a practitioner.

Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

5 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and overstretch or distort.

Keep away from sharp objects (such as jewellery, fingernails).

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for shower use

6 Donning the Device

Before you begin

Fit the valve. (See *Fitting the Valve* section.)

Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.




Do NOT pull or stretch the device.

1. Invert the device so that its silicone side points externally.
2. Align the locking pin with the long axis of the residual limb.
3. Roll the device on the residual limb while releasing any trapped air.

Note: The unperforated fabric should be below the patella tendon.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness, tingling or any unusual sensation within 10 minutes, do the following:





4. Doff the device, and wait until normal sensation returns.
5. Don the device again.
6. If the user feels numbness, tingling or any unusual sensation again, stop using the device.

7 Fitting the Valve


⚠ Do not overtighten the locking pin.

We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.

8 Fitting the Spacer

Fit the spacer between the pin lock and the valve only if the pin lock prevents the valve from opening.

9 Fitting Advice

9.1 Moisture in the Distal End of the Device

If moisture collects inside the distal end of the device, refer to the following table:

Cause	Solution
The distal perforations are clogged.	Clean the device. (See <i>Maintenance</i> section.)
The socket does not have enough space in its distal end for the valve to open.	Put the spacer between the pin lock and the valve. (See <i>Fitting the Spacer</i> section.) If unsuccessful cast another socket.
The internal shape of the lock-body is stopping the valve from opening.	Put the spacer between the pin lock and the valve. (See <i>Fitting the Spacer</i> section.)

9.2 Loss of Vacuum

If the device loses vacuum, refer to the following table:

Cause	Solution
The valve is damaged.	Stop using the device, and contact a Blatchford sales representative.
Debris beneath the valve preventing correct seal.	Clean the valve. (See <i>Cleaning the Valve</i> section.)
The locking pin is not aligned with the long axis of the residual limb.	Re-align the locking pin with the long axis of the residual limb.

9.3 Separation of the Valve from the Distal Cap

If the valve falls off the distal cap, refer to the following table:

Cause	Solution
The valve is not fitted correctly.	Make sure that the washer is between the pin and the valve. (See <i>Fitting the Valve</i> section.)

Note... Do not over-tighten the locking pin/lock rod.

9.4 Oedematous Swellings Corresponding to Distal Perforations


Cause	Solution
Excessive distal contact	Reduce end pressure by addition of extra socks or lengthen/remake socket

10 Technical Data

Principal Materials	polyamide and lycra, silicone, nylon
Shore Hardness	40 shore 00
Component Weight (<i>size 28</i>)	695 g (1 lb 8 oz)
Activity Level	3-4
Size Range	22-40 cm
Length (See diagram below)	435 mm
Internal Length (See diagram below)	420 mm
Matrix Length	Approx. 10 cm
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)
Distal End Attachment	M10 Requires Locking Pin With a Shoulder*
Locking Pin Shoulder Diameter	13.5-19 mm

*Locking Pin Not Supplied

Dimensions

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Warranty

This device warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

Trademark Acknowledgments

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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© Blatchford Products Limited 2024. All rights reserved.
Endolite India Ltd.
Blatchford Europe GmbH