Tectus

Clinician's Manual

TEC-L, TEC-R

EN	Clinician's Manual	2
FR	Manuel de l'orthoprothésiste	67
DE	Gebrauchsanweisung für Techniker	132
IT	Manuale per il tecnico ortopedico	197
NL	Handleiding voor artsen	262
PT	Manual do profissional de saúde	327

Contents	2
Package Contents	3
1 Description and Intended Purpose	4
2 Safety Information	7
3 Component Identification	9
4 Function.....	11
5 Maintenance	12
5.1 Critical Service Indicator.....	13
5.2 Hard Reset.....	13
5.3 Cleaning	13
5.4 Replacing the Key Fob Battery	14
6 Limitations on Use	15
7 Battery Charging	16
7.1 Battery Charging Safety Information	16
7.2 Before Charging.....	17
7.3 On-Screen Battery Level Indicator	17
7.4 Connecting and Disconnecting the Charger	18
7.5 Battery Charger LED Indications.....	19
7.6 Charging and Operating Reference Times	19
8 Bench Alignment.....	20
9 How to Don and Doff the Tectus KAFO	21
9.1 Donning the Tectus KAFO.....	21
9.2 Doffing the Tectus KAFO	22
10 Turning the Device ON/OFF	22
11 Activating <i>Bluetooth</i> ® on the Device	24
12 Device Setup	25
12.1 Connecting to the Blatchford App.....	25
12.2 Tectus Programming Screens	27
12.3 Calibrating the Operating Modes for the Specific User	31
13 Key Fob Programming.....	37
13.1 Assigning Operating Modes to Key Fob Buttons.....	37
13.2 Clearing Key Fob Button Assignments	37
14 Tectus Knee Module Operation	38
14.1 User Modes.....	38
14.2 Safety While Switching Between User Modes	38
14.3 Device Over-Stress Warning	39
14.4 Activities and Associated User Modes.....	40
14.5 Operating Actions and Device Indications	44
14.6 Mode Switching on the Device Using Tectus Module Buttons A and/or B	45
14.7 Mode Switching on the Device Using the Key Fob.....	46
14.8 Additional Key Fob Functions.....	46
14.9 Mechanical Override Operation	47
15 Troubleshooting	48
15.1 Hydraulic System Optimization	48

16 Disassembly/Assembly Instructions..... 49

 16.1 Disassembly 49

 16.2 Assembly..... 50

17 Changing the Extension Assist Spring..... 54

18 Transport and Storage..... 55

19 Technical Data 56

Appendix 1 Installing the Blatchford Programming App 58

Appendix 2 Clinician Registration on Blatchford App 59

Appendix 3 Gait Visualizer 61

Manufacturer’s Declarations and Legal Information 63

 Please read all safety warnings in Section 2 and 7 before setting up or using the device.

Package Contents

- 1. Tectus knee module
- 2. Tectus follower joint
- 3. Knee joint protective covers
- 4. Battery pack, mount and charger
- 5. Set of three compression springs (one is pre-installed in the device)
- 6. Key fobs & lanyards (2 sets)
- 7. Clinician’s Manual
- 8. User Guide

Notes: i) Contents may vary from list depending on customer order.

 ii) Side-bars may be shipped separately.

 iii) The completed Tectus KAFO is supplied with a limb bag.

1 Description and Intended Purpose

These instructions are for use by the practitioner only.

The term *device* is used throughout these instructions for use to refer to Tectus, a knee system comprising Tectus knee module, follower joint, battery pack, side-bars and associated key fob controllers. The terms *Tectus KAFO* (Knee-Ankle-Foot Orthosis) or *brace* are used throughout these instructions for use to refer to a custom-made full leg orthosis incorporating the device and manufactured to a practitioner's prescription.

Note... The ankle joint is not manufactured by Blatchford and not covered by these instructions for use. Refer to manufacturer's documentation.

The term *programming device* is used throughout these instructions for use to refer to a tablet or mobile phone running an Android or iOS operating system.

The term *app* is used throughout these instructions for use to refer to **Blatchford** programming application which resides on the programming device.

Ensure that the user understands all instructions for use, drawing particular attention to the maintenance and safety information.

Application

This device is to be used exclusively as part of a lower limb orthosis.

Intended for a single user.

This device is a Class I CE marked medical device according to the Medical Device Regulation 2017/745 that is for use as a component of a lower limb orthosis. The brace is a custom-made medical device. The device, when assembled as part of a KAFO, is a microprocessor-controlled swing and stance device allowing automatic adjustment of flexion resistance and manual adjustment of extension assistance.

The device enables the user to undertake activities such as walking, standing, sitting and descending stairs. The user (wearer) can control the device by switching it between different operating modes for different activities, using control buttons on the device or control buttons on a key fob. The practitioner can also use these controls in conjunction with the **Blatchford** app to customize both swing and stance to the individual's requirements.

The practitioner can also fine tune the extension assistance by adjusting the pre-load/compression on the extension spring or a more coarse adjustment if needed by exchanging extension bias/assistance springs within the unit. Flexion resistance is applied hydraulically, and the level of flexion resistance applied depends on which operating mode is selected; the practitioner can adjust these resistance levels for individual operating modes via the **Blatchford** app.

Features

- Dynamic damping resistance improves gait characteristics, enabling safer walking
- Microprocessor control of flexion resistance in real time
- Mechanical adjustment of hydraulic extension assistance
- Device software adapts to individual user gait characteristics
- **Bluetooth®** communications between programming app and the device
- Power supplied via an integrated rechargeable battery pack
- Recharging through an easy-to-use magnetically self-aligning charger
- Device and battery pack splash proof to IP54

Indications

The device is intended solely for use as a lower limb orthosis and may only be prescribed by a practitioner trained and certified by Blatchford to prescribe Tectus. After considering all prescription parameters, the device may be suitable for users with the following conditions:

- Hip flexor muscle strength greater than power 2 on the Oxford scale and the user can initiate swing phase from their hip with muscle power with compensatory movements
- Quadriceps weakness
- Flaccid paralysis of lower limb
- Resultant paralysis from orthopaedic injury or procedure
- Insufficient knee control
- Spinal cord injuries
- Stroke patients (mild tone only)
- Poliomyelitis
- The user can stand independently

Note... The list above is not exhaustive. The device may be suitable for patients with other conditions after due consideration by a practitioner.

The Blatchford-certified practitioner will ultimately be responsible for deciding whether the Tectus KAFO is suitable for the intended user, and any prescription must take into account all indications listed above, and any standard outcome measures including but not limited to:

- The two-minute walk test (2MWT)
- The L test of functional mobility
- 6 minute walk test
- Pain score
- Berg Balance Test

In addition, the user's muscle strength, as referencing by the Oxford Scale shown below, must achieve grade 3 or above in their hip flexors (or 2 if they can initiate hip flexion with compensatory movements).

Oxford MRC (Medical Research Council) Scale

Grade	Muscle Activity
0	No contraction
1	Flicker/trace contraction
2	Active movement with gravity eliminated through full range
3	Active movement against gravity through full range
4	Active movement against gravity and resistance through full range
5	Normal power through full range

The user's activity level should also be considered using the descriptions below. The device is mainly intended for users who can achieve activity level 2 or 3 but can be used by activity level 1 or 4 users as well according to the recommendations of and limitations set out by a practitioner.

Orthotic User Activity Levels

Activity Level 1: Indoor User

Able, or may be able, to use an orthosis for transfers and walking slowly on level surfaces and for short distances indoor.

Activity Level 2: Restricted Outdoor User

Able, or may be able, to use an orthosis to walk indoors or outdoors with some variation in speed and able to traverse low-level environmental obstacles such as curbs, single steps, or uneven surfaces.

Activity Level 3: Unrestricted Outdoor User

Able, or may be able, to use an orthosis to walk at varied speeds indoors and outdoors and potentially traverse most environmental obstacles. Able, or may be able, to carry out other movements/exercises in addition to walking using an orthosis.

Activity Level 4: High Activity Unrestricted Outdoor User

Able, or may be able, to use an orthosis for activity levels greater than Activity level 3, creating higher impact loads, and increased stresses and strains on the orthosis, e.g. activity level of a child, active adult, or athlete.

Contraindications

- Not to be used bilaterally on both legs
- Unable to initiate hip flexion to initiate swing phase from muscle power or with compensatory movements
- Knee or hip flexion contracture greater than 10 degrees
- 10 plus degrees genu varus or valgus that can't be reduced
- Compromised trunk control
- A leg length discrepancy more than 15 cm
- Muscle conditions that prevent control of an orthosis including moderate or severe spasticity, or significant muscle tone
- Lack of sensory or cognitive ability necessary to operate and perceive device controls and indicators respectively
- Body weight is greater than specified in Section 19 *Technical Data*
- Insufficient manual dexterity to operate device controls

Clinical Benefits

- Supports and enhances the function of the affected limb
- More natural gait mechanics, compared to conventional KAFOs and simple stance KAFOs
- User will use less compensatory gait mechanisms (e.g. hip hitching), compared to fixed KAFOs and simple stance KAFOs
- More symmetrical loading between limbs during activities that use the flexion resistance of the device, for example transitioning from standing to sitting


2 Safety Information

-  This warning symbol highlights important safety information which must be followed carefully.
-  Advise the user to contact their practitioner if their condition changes.
-  Advise the user to regularly check the limb on which the Tectus KAFO is worn for signs of skin irritation, tissue damage, pressure sores or discomfort and report any such signs to their practitioner.
-  Assembly, programming, maintenance and repair of the device (including handling and replacement of the battery unit) must only be carried out by a suitably qualified practitioner or technician that has attended an approved training course.
-  Be aware of finger trap hazards at all times.
-  Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.
-  Avoid strong magnetic fields, radioactive environments, sources of electrical interference, atmospheres containing liquids and/or powders.
-  Do not wear the device during any scanning procedures (e.g. X-ray, CT or MRI scan).
-  The device is not intended for use when immersed in water or as a shower orthosis. If the orthosis comes into contact with water remove and dry immediately.
-  Any excessive changes in heel height after programming will adversely affect the function of the orthosis. Only footwear approved by the practitioner should be used with the orthosis.
-  Contact Blatchford or your service provider for advice on whether the device can be used by a user fitted with a pacemaker or any other electronic medical device.
-  Ensure no water enters the charger plug socket and the cover is only removed during charging.
See also additional safety information in Section 7.1 *Battery Charging Safety Information*.
-  Ensure that torque drivers used on the device are calibrated.
- Advise the user to observe the following warnings when using the device:
-  Any changes in the performance or function of the device e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your practitioner or service provider.
-  Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles. While wearing the orthosis, the leg on which the orthosis is worn must not be used to control the vehicle or its components (e.g. clutch pedal, brake pedal, accelerator).
-  The user must not operate heavy machinery while wearing the device.
-  Always use a handrail when ascending and descending stairs and at any other time if available.
-  Leg-over-leg stairs descent should only be undertaken by users who have received suitable training and have been approved to do so. Otherwise a step-to-step gait should be used descending stairs.

-  When starting to walk, the user should always lead with the unaffected leg. If walking is initiated with the braced limb, the device may not recognize the movement and as a safety measure will switch to *Lock* mode.
-  Confirm that *Stairs* mode is engaged before attempting to descend stairs.
-  Confirm that *Sit* mode is engaged before attempting to sit.
-  The device must not be used for unusual activities, including, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sport activities that involve jumping, sudden movements or fast sequences of steps (e.g., running, cycle racing, basketball, badminton, ice and snow sports) or extreme slopes and steps.
-  Avoid, where possible, using the device on slippery, uneven or unfirm surfaces (sand, snow, mud, cobbles etc.). If unavoidable, proceed with extreme caution.
-  If the critical service indicator is activated and cannot be reset (see Section 5.1 *Critical Service Indicator*), stop using the device as soon as is practicable and contact your service provider to arrange a service.
-  The user must not adjust or tamper with the setup of the device.
-  Do not modify or replace the KAFO straps or padding prescribed by the practitioner.
-  After continuous use, the knee module may become hot to the touch.
-  To minimize the risk of slipping and tripping, appropriate footwear that fits securely around the foot components must be used at all times.
-  Carrying loads might affect the user's posture and ability to meet the required parameters for the knee to release from lock.
-  Keep the key fob on a lanyard separate from any other keys to prevent accidental mode switching.
-  To prevent damage to clothing, take care not to trap clothing in the moving parts of the device.
-  Avoid any strangulation hazard that the key fob lanyard may present to the user when worn, such as catching the cord in rotating machinery.
-  The device is not designed for walking backwards. Do not attempt to walk backwards as the device may not support you.


In *Walk* mode, when backwards motion is detected, the device will lock in extension to provide standing support only. See Section 14.4 *Activities and Associated User Modes* for details.
-  Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
-  **WARNING:** Portable RF (radio frequency) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer; this may otherwise result in the degradation of performance of the device.

3 Component Identification


Note... For illustrative purposes, only the right leg Tectus KAFO is shown.


Tectus Knee Module (Right leg shown)


Follower Joint (Right leg shown)

Battery Module

Knee Side-bars

938467

4 Function

The device features a microprocessor-controlled system which uses measurements from integrated sensors as a basis to adjust a hydraulic unit influencing the flexion resistance of the stance phase. The device also regulates the extension of the swing phase using a manually adjustable extension assist spring.

During use, the device dynamically adapts in real time to the current motion situation (gait phase). An orthotist can make adjustments to the control system via the **Blatchford** App to adapt the device to the individual user's needs, see Section 12 *Device Setup*. The App runs on both Android and iOS programming devices (phone or tablet) and communicates with Tectus over a **Bluetooth®** wireless connection.

Extension assistance (extension bias) is controlled by a spring mechanism. The practitioner has a choice of 3 springs of different stiffness to allow coarse adjustment. (By default the device is fitted with the stiffest spring, part number 925042.) The practitioner can further adjust the extension assistance manually by tightening or loosening the selected spring with a hex key, see Section 12.3.5 *Tuning the Knee Extension Assist*.

In the event of any malfunction of the electronic or the hydraulic control system, or when the battery is fully discharged, the device automatically switches to a **Safety** mode. In this **Safety** mode the device reverts to pre-programmed resistance parameters and allows only limited functionality. Normal operation resumes when the device is first time powered on after the malfunction has been resolved, or the battery has been recharged, see Section 5.1 *Critical Service Indicator* and Section 7.3.2 *Low Battery Behavior*.

The device allows the user to switch operating modes manually using the 4-button key fob and/or A and B buttons to adapt device operation for specific circumstances. For example, to walk down stairs they would select **Stairs** mode just prior to starting to walk down the stairs. In **Stairs** mode the device will provide necessary additional support by increasing flexion resistance. Operation of the device and the operating modes available to the user are described in Section 14 *Tectus Knee Module Operation*.

A four-button key fob on a lanyard provides the user with a remote (Radio Frequency) control to switch the device between operating modes, see Section 14.7 *Mode Switching on the Device Using the Key Fob*. If required, the key fob functions can be reassigned via the app, see Section 13 *Key Fob Programming*.

A power ON/OFF button and two push buttons (A and B) on the device provide the local control through various switch pressing combinations. The device has a mechanical override button which allows knee flexion while pressed.

5 Maintenance

The Tectus KAFO requires a routine inspection at the following times after fitting:

6 months	12 months	24 months	36 months
----------	-----------	-----------	-----------

These inspections must be carried out by competent personnel according to the routine inspections schedule available from Blatchford ACS.

 **Completion of routine inspections according to the prescribed maintenance schedule and notification of Blatchford of the results is necessary to maintain the validity of the warranty.**

If the device is still under warranty and if the module alone is returned for maintenance or repair, than a replacement device may be made available to the user on loan while maintenance is being carried out. Contact Blatchford/your service provider for details.

The loaned device and the device returned after maintenance must be set up by the practitioner before use.

The user must be handed the User Guide supplied, and be advised of the following:

- Any changes in performance of this device must be reported to the practitioner. Changes in performance may include:
 - Increase or decrease in knee stiffness
 - Instability
 - Reduced knee support (free movement)
 - Any unrecognized beeps, warning lights, on-screen indications or device vibrations*
 - Any unusual noise
- Check the device every time before use for physical damage to the KAFO and:
 - Signs of wear or tear
 - Hydraulic leakage
 - Security of the joint covers
 - Security and condition of the straps

*The device emits vibrations when operating modes change. This is normal behavior.

 **Do not lubricate the device yourself, this will be carried out during servicing with the appropriate lubricants.**

 **Never use WD40 or any other water dispersant on any of the device moving parts in an attempt to fix unusual noises, squeaking or stiff joints.**

Advise the user to replace the key fob battery either every six months or:


- when switching between knee modes becomes inconsistent or
- if the key fob LED does not illuminate every time a button is pressed

See Section 5.4 **Replacing the Key Fob Battery**. If the user is unable to do this they must contact their service provider or clinician to arrange changing their battery.

 **Ensure that torque drivers used on the device are calibrated.**

5.1 Critical Service Indicator

In the event of the device detecting a malfunction, the critical service indicator graphic will be activated and flash on the device screen (see image) and the device will attempt go into **Safety** mode with higher than normal knee flexion and extension resistance. Simultaneously button B will start and continue flashing Red until the malfunction is resolved.

In this **Safety** mode the users can slowly extend or flex the knee to allow them to move to a safe position to remove the device.

If for any reason Safety mode does not engage, the user must stop using the device immediately.

If adjustment of the knee angle is required, the device can be allowed to flex or extend freely by operating the override button, please refer to Section 14.9.

WARNING! Only use the override button if it is safe to do so.

Inform the user that if the device malfunctions, and performing a hard reset (see Section 5.2) does not fix the fault, or if the Critical Service Indicator is displayed, they must contact their practitioner as soon as is possible, and must not use the device until the malfunction has been fixed and/or service has been carried out as necessary.

The critical service indicator will continuously flash and the knee will remain in this mode until the malfunction is resolved. If you power down with the critical service indicator on, it will extinguish however when you next power on the device the critical service indicator will begin to flash again and continue until the malfunction is resolved.

5.2 Hard Reset

In the case of a minor and non-recurring malfunction a hard reset may be attempted to allow device self-recovery. To conduct a hard reset, power the device OFF and ON via the power ON\ OFF button located on the knee module.

If the malfunction is persistent contact your service provider.

5.3 Cleaning

Blatchford recommends cleaning the device daily. Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers or bleach. Do not machine wash/tumble dry or use any hose, pressure washer/drier system.

Ensure water/liquids do not enter the charging socket. Thoroughly dry before use.

5.4 Replacing the Key Fob Battery

The key fob operates using a CR 2032 button battery. Replacement is necessary every six months or when switching between knee modes becomes inconsistent or if the key fob LED does not illuminate every time a button is pressed.

1

Place the key fob on a flat surface with the battery compartment lid facing upwards.

2

Using an appropriate cross head screwdriver, remove both screws and put aside for later.

3

Using an appropriate flat head screwdriver pry open the cover only from the key ring end as shown to prevent damage to the gasket.

4

Remove the discharged battery and dispose it according to the local recycling regulations.

Batteries should be kept out of reach of children and pets. In case of ingestion, seek medical attention.

5

With the side marked + facing outwards, place the new battery into the battery compartment by pressing it first against the side contact springs then into the compartment.

6

Ensure that the mating surfaces of the battery compartment and lid are clean. If required, use a dry cotton bud to clean these surfaces.

7

To ensure continued protection against ingress of moisture, make sure that the lid gasket is fitted without creases and free of damage.

8

Replace cover and screws as shown and tighten screws until hand tight. Do not overtighten!

After battery replacement, confirm all key fob functionality by testing that all buttons achieve the required knee mode change when pressed and that the LED on the fob illuminates once after each time a button has been pressed.

6 Limitations on Use

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment. The assessment should take into account that carrying loads might affect posture and ability to meet the required parameters for the knee to release from stance.

The combined user weight and lifting load must not exceed 100 kg in all cases.

Environment

Avoid exposing the Tectus KAFO to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Note... Only the knee module and battery pack are waterproof to IP54. The key fob is fully waterproof to IP68.

Exclusively for use between -10 °C and 50 °C
(14 °F to 122 °F).

Suitable for outdoor use

7 Battery Charging

 Please advise the user of the charging information in this section.

7.1 Battery Charging Safety Information

-  Never connect the charger to the device whilst still wearing it.
-  The battery should only be charged in temperatures between 0 °C to 35 °C (32 °F to 95 °F).
-  The battery charger can be used with electrical outlets rated between 100V and 240V AC, 50/60 Hz.
-  Only charge the battery of the device using the charger and charging cable supplied with the product. Do not use this charger to charge any other device.
-  Ensure the appropriate, country-specific mains adapter plug is fitted to the charger. Never force the plug into an incompatible mains socket.
-  Before charging ensure charging surfaces are dry, free from contamination and debris build-up, because this can prevent charging.
-  If any unusual noise or smell is noticed whilst charging, disconnect the charger and power off both charger and knee module. Advise user to report the problem to their practitioner.
-  The battery module and battery charger connectors contain strong magnets. Always keep a safe distance (at least 10 cm) between these magnets and all objects that can be affected by magnetism, such as watches, pacemakers, bank cards or any magnetically stored media.
-  Do not allow any ferrous material to attach itself to the charging surfaces.
-  Do not dismantle the adapter plug or battery charger.
-  Immediately replace damaged adapter plugs or battery chargers with parts supplied by Blatchford.
-  During a charging the device will power down. Once charged, the device will require powering on via the ON/OFF button located on the outer surface of the Tectus knee module.
-  Do not remove the device batteries, they are not field replaceable and will be replaced if necessary, during routine inspection.
-  If the device has been subjected to temperatures lower than -10 °C (14 °F) it must be returned to Blatchford for inspection for possible damage to the battery pack.

7.2 Before Charging

The device should not be charged while being worn; remove the Tectus KAFO before charging. Do not don the device while it is connected to the charger.

Always check the battery charge level before use, ensuring a sufficient amount of charge for daily activities. It is recommended to charge the device every day when using it on a daily basis.

It is recommended that you switch the device to **Lock** mode and, before charging, lay the Tectus KAFO down in a horizontal position (on a flat surface or floor), with the knee module facing upwards and with 5 degrees of knee flexion so that it is secure and does not roll over.

Note... As soon as the charger is connected, the knee module automatically switches off.


7.3 On-Screen Battery Level Indicator

7.3.1 Normal Indications

The battery level indicator will appear for 5 seconds on the side of the on-board graphic display:

- During mode switching.
- Long pressing button A and B together.

The following battery level indications are available:

7.3.2 Low Battery Behavior

Low Battery Warning

The device will warn the user when the battery charge level shown drops to 20% by emitting two short high-tone beeps and also displaying the warning symbol followed by one bar on the graphic display. This warning will be repeated every 10 minutes and/or when any button is pressed on the device until the charge in the battery reaches the critically low level.

Critically Low Battery Warning

A critically low battery level will be indicated by the warning symbol flashing and one long high-tone beep. The device will revert to a Safe Mode identical to the one described in Section 5.1

Critical Service Indicator.







When the Critically Low Battery warning is indicated, the user must stop using the device and not to reuse until the device has been recharged.


7.3.3 No Display after Connecting the Battery Charger

If after connecting and switching on the battery charger the device graphic display remains blank, the battery may be completely discharged. Leave the battery charger connected for at least 15 minutes and check the charge level by disconnecting and re-connecting the battery charger.

7.4 Connecting and Disconnecting the Charger

7.5 Battery Charger LED Indications


Note... Refer to the reverse side of the battery charger module and the instructions supplied with the charger for further information.

Start-Up Sequence

With the charger plugged into the mains socket, and connected to the device, switch on at the socket. The charger LED will illuminate in the following color sequence in approximately 10 seconds to indicate the device is charging.


7.6 Charging and Operating Reference Times

Process	Time
Time to fully charge battery. (From fully discharged state.)	Approximately 2 hours
Minimum charging time before use. (From fully discharged state.)	At least 1 hour
Maximum battery duration, from fully charged, under continuous device operation i.e. uninterrupted walking. (Subject to operating temperatures)	18 hours (see notes i and ii)

- Notes:
- i) For the maximum operating time with one battery charge, turning the device off when it is not being used is recommended.
 - ii) The operating time of the device may be reduced at low ambient temperatures or due to aging of the battery.


8 Bench Alignment

Bench alignment and fitting process should be completed before all other adjustment are made, and repeated by a practitioner after 6 weeks of use to allow for any adaptation/acclimatization that has taken place. Alignment confirmation may result in adjustment or rebuild of the Tectus KAFO.

The anatomical compromise knee centre is located as shown below, in line with conventional procedures as defined by M. Nietert^[1]. The mechanical centre of rotation of the device knee joint should fall slightly behind and below the anatomical compromise knee centre. The amount of additional posterior offset will depend on the sufficiency of knee and hip securing muscles.

The clinician provides Blatchford with user measurements such as height, weight, knee A-P width, etc. and then Blatchford calculates the exact position of the device's mechanical knee centre for them.

D Horizontal A-P depth of the knee.

● Device mechanical centre of rotation.

○ Anatomical compromise knee centre.

For the purpose of alignment Blatchford recommends using a 10-12 mm heel wedge to achieve 10 degrees anterior tilt of the lower leg. Additional heel wedge height can be added to accommodate any existing restrictions in ranges of motion at the ankle joint, such as a plantarflexion contracture. We recommend 7-12 degrees of toe/forefoot extension which will have been captured during the casting and manufacturing stages.

[1] Nietert, M. The Compromise Pivot Axis of the Knee Joint, Aachen: shaker, 2008.

9 How to Don and Doff the Tectus KAFO

Advise the user to don and doff the Tectus KAFO following the steps detailed below.

9.1 Donning the Tectus KAFO

Notes: i) The device should be charged before donning. (See Section 7 **Battery Charging**.)
ii) Advise the user to inspect daily for signs of deterioration or damage that might affect the safety or function of the device. If anything unusual is noticed, they should contact their practitioner for advice.

1. Completely undo and fully open all closure straps so they are clear of the shell openings.
2. Remove the shoe.
3. Sit on the front edge of a chair.
4. Power ON the device. (See Section 10.)
5. If powered ON with the device lying horizontal, the device should be in **Free** mode after power ON; if not, then set the device to **Free** mode via key fob or module buttons.
6. Flex the Tectus KAFO at the knee.
7. Slightly extend the leg and insert the thigh into the thigh shell.
8. Insert the foot into the foot component, positioning the heel first and insert the lower leg into the lower leg shell.
9. Fasten all straps. Do not overtighten.
10. Check the fit of the orthosis for comfort. (Alter as required e.g. appropriate trim lines.)
11. Put the shoe on over the foot component. (Some users may prefer donning the Tectus KAFO with the shoe already in place.)
12. Set the device to **Walk** mode and stand up (holding on to chair for support as necessary) and adjust all closures/securing straps as required.
13. If worn over clothing, ensure clothing is folded neatly to avoid discomfort.
If worn under clothing (i.e. dress, skirt, loosely fitting trousers), a sock or a thin under garment covering the full length of the leg should be worn.

To prevent damage to clothing take care not to trap clothing between the moving parts of the device.

14. Set the device to **Sit** mode and sit down again carefully.
15. Set the device to **Free** mode and gently flex and extend the knee to confirm that Tectus KAFO fits comfortably throughout the range of movement.
16. Set the device back to **Walk** mode before standing up.

9.2 Doffing the Tectus KAFO

1. Set the device to **Sit** mode.
2. Sit down on the front half of a chair (holding on to chair for support as necessary).
3. Set the device to **Free** mode.
4. Remove the shoe.
5. Open all fastening straps.
6. Remove Tectus KAFO.
7. Extend the Tectus KAFO to 5 degrees knee flexion and activate **Lock** mode on the device.
8. Power OFF the device.
9. Lay the Tectus KAFO down in a horizontal position (on a flat surface or floor), with the knee module facing upwards and with 5 degrees of knee flexion so that it is secure and does not roll over.

Note... It is recommended to put the device on charge every time after use.

Advise the user to regularly check the limb on which the Tectus KAFO is worn for signs of skin irritation, tissue damage, pressure sores or discomfort and report any such signs to their practitioner.

10 Turning the Device ON/OFF

The device is powered on or off via the Power ON/OFF switch on the device (see Section 3); there is no Power ON/OFF switch on the key fob.

Powering ON

To power on the device, press and hold the ON/OFF switch until LED illuminates green. Short beeps and vibrations will then occur during the start-up sequence, when they finish, then release.

Powering OFF

To initiate the shut-down sequence, press and hold the ON/OFF switch until two short beeps with two short vibrations are emitted, then release.


Start-up Sequence

Powering ON initiates the automatic start-up sequence during which the LED and graphic display indications change, audible beeps sound and the device vibrates as shown. The sequence lasts for approximately 4 seconds.

Note... If the programming of the device has not been previously completed, the BLUE LED will continue to flash on completion of the start-up sequence, indicating that *Bluetooth®* has been automatically enabled. If programming has been completed see operation in Section 11.

The device is powered ON and ready for use when the default *Free/Walk* mode (depending on device orientation during start-up) is displayed on the on-board graphic display, (see sequence shown on right).

Note... The display goes into sleep mode (display goes blank) after 5 seconds of inactivity but is reactivated after every user mode change initiated by either buttons A/B on the device, or buttons on the key fob.


11 Activating *Bluetooth*® on the Device

After the Programming Complete button has been pressed on the Programming screen of the *Blatchford* app, the device will not switch on its *Bluetooth*® at startup. Before connecting to the app on the programming device the *Bluetooth*® function must be manually activated.

It is recommended for the user to sit down while *Bluetooth*® connection is being established.

1


On-screen and Audible Timing Aids

Button B


To initiate *Bluetooth*® connection press and release button B after the above sequence of screen indications and accompanying beeps have occurred (approx. 3 s). Successful *Bluetooth*® activation is indicated by two low tone beeps and vibrations.*

2


The Power button LEDs illuminate flashing blue, indicating that *Bluetooth*® is active on the device and the device is discoverable by the programming app.

3

Ensure the Tectus module is in range of the programming device (within 5 m). Follow the procedures in Section 12.1 to connect via *Bluetooth*® to the programming app. The LED flashes blue while establishing connection, or until time-out has been reached.


4

After successful connection of the app to the device and when the Programming button on the app is tapped, the blue LED will stop flashing and remain illuminated steady blue until the connection is terminated or lost.

If connection is lost the LED indicator will return to the discoverable state (flashing blue).

5

Live *Bluetooth*® Connection Indicator

Successful *Bluetooth*® connection is also indicated by the on-board graphic display showing the above symbol accompanied by one medium tone beep and vibration.

6

Termination of *Bluetooth*® connection can be initiated via the programming app by selecting DISCONNECT in the main menu.

Following disconnection, the device's *Bluetooth*® functionality can be turned off by:

- pressing button B until all three on-screen and audible timing indications (see box 1) have occurred.
- or waiting for it to automatically time-out after a period of inactivity.

In all cases the Power button LED extinguishes.

Powering the device off will result in both loss of *Bluetooth*® connection and turning off the device's *Bluetooth*® functionality.

*See description of audio and visual cues in Section 14.5 *Operating Actions and Device Indications*.


12 Device Setup

The device is set up and adjusted primarily using the **Blatchford** app running on an Android or an iOS mobile device (phone or tablet). In addition, the correct extension assist spring must be selected and spring compression must be manually adjusted to set the correct amount of extension assist (see Section 12.3.5).


12.1 Connecting to the Blatchford App

Prerequisites:


- The **Blatchford** app is installed on the programming device (see Appendix 1) and
 - **Bluetooth**® is enabled on the programming device.
 - The practitioner or other user of the application has been granted appropriate authorization to use the **Blatchford** app by Blatchford. Authorization is provided as a result of attending and passing an approved training course.
1. Power on the device (see Section 10).
 2. Activate **Bluetooth**® on the knee module (see Section 11).
 3. Tap the **Blatchford** icon on the menu screen of your programming device.

4. On the **WELCOME** page click **LOG IN**. (For registration click **CREATE AN ACCOUNT**. See Appendix 2 for details.)


5. Enter login details then Press **SIGN IN**.


6. Tap **SELECT A DEVICE** to show devices available for you.


7. Make sure the device is switched on and **Bluetooth®** is enabled then choose Tectus.


8. Choose the Tectus device with the connection code matching the device you intend to connect to. (There may be more than one available nearby.)

9. Enter **Connection Code** (a 4-digit **Bluetooth®** number printed on the Tectus KAFO label) then tap **CONNECT**.

10. The connection status updates as the **Blatchford** app connects to the device.


11. When the **Bluetooth®** connection has been successfully established, tap **PROGRAM** to navigate to the main programming screens. (See Sections 12.2.1 and 12.2.2.) Options displayed depend on authorisation level.

12.2 Tectus Programming Screens


On selecting **START PROGRAMMING** the **BASIC SETUP** screen appears (scroll down on the app to reveal all). Its functions and controls are described in Section 12.2.1.

An **ADVANCED SETUP** screen is also available by selecting **ADVANCED SETUP** from the **BASIC SETUP** screen. The **ADVANCED SETUP** screen and its additional functions and controls are described in Section 12.2.2.


12.2.1 BASIC SETUP Screen

- a) **Battery.** Current battery level indicator (percentage remaining) displayed in 20% increments.
- b) **State.** Real-time indication of device operating mode sub-states during **Walk** or **Sit** operating modes.
- c) **Mode.** Indication of the current user mode activated.
- d) **Leg Side.** Leg selected for user is displayed.
- e) **Thigh Angle.** Real-time indication of stride profile. The black needle shows the current thigh angle reading in degrees while the red shows the average of the maximum thigh angle values reached. (Readout only – no adjustment.)
- f) **Hydraulic Resistance.** Real-time hydraulic valve position indicative of hydraulic resistance in knee flexion. (Readout only – no adjustment.)
- g) **User Standing Calibration.** To be used as part of the standing calibration process described in Section 12.3.1 *Standing Calibration in Walk Mode*. To confirm completion, tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen.
- h) **Sitting Resistance.** Sets the amount of flexion resistance that will support the user as they sit down in **Sit** mode. Adjust value up or down with +/- buttons. **Important:** Press **CONFIRM** to confirm value then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Value displayed changes to green while being altered. After confirmation it turns black again.
- i) **Stance Release Angle.** Sets thigh angle (extension) that must be achieved for the knee to release into swing. **Important:** Press **CONFIRM** to confirm value then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Value displayed changes to green while being altered. After confirmation it turns black again.
- j) **Stairs Resistance.** Sets the amount of flexion resistance applied in **Stairs** mode. **Important:** Press **CONFIRM** to confirm value then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Value displayed changes to green while being altered. After confirmation it turns black again.
- k) **Swing Adjustment.** Adjusts the hydraulic swing resistance. This also complements the spring-operated knee extension assist feature. An adjustment to the swing control value will instruct the microprocessor to make an adjustment to the hydraulic valve resulting in a reduction in knee velocity during the swing phase. The higher the value the greater the hydraulic resistance. **Important:** Press **CONFIRM** to confirm value then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Value displayed changes to green while being altered. After confirmation it turns black again.

- l) **Pressure Release Threshold.** Sets the pressure at which the knee releases during walking. Pressure can be set between 20 and 30 bar. **Important:** Press **CONFIRM** to confirm value then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Value displayed changes to green while being altered. After confirmation it turns black again.
- m) **Set Device Mode.** Enables selection of user modes (**WALK, SIT, LOCK, STAIRS, FREE**) on the device via the app. Tap the button corresponding to the required device mode change then tap **CONFIRM AND CONTINUE** button on the **Safety alert** screen. Button illuminates teal (blue-green) when selected operating mode has been confirmed by the device.
- n) **COMPLETE.** When you have finished calibration, this button should be pressed and confirmed on the **COMPLETE** screen by tapping **CONFIRM** (selection results in the button illuminating green). It serves as a reminder that you have completed all calibration tasks. If after pressing this button you start adjusting any of the parameters again the button will switch off (display grey).
- Note... As a result of completion of programming, the automatic **Bluetooth®** activation during start-up on the device will be deactivated.
- o) **TECTUS DOCUMENTATION.** Provides link to Tectus documentation.
- p) **CONTACT SUPPORT.** Provides link to contact support.
- q) **ADVANCED SETUP.** Navigates to the **ADVANCED SETUP** screen.
- Note... **Safety alert** screen. Appears after making a selection and requests confirmation of selections made, such as changing operating mode.

12.2.2 ADVANCED SETUP Screen

- a) **Training Mode.** Provides the user with audible and vibration cues in this mode while the user learns to use the device. Tap **ENABLE/DISABLE** button then **CONFIRM AND CONTINUE** button on the **Safety alert** screen to enable/disable **Training Mode**. Button illuminates teal (blue-green) when selected operating mode has been confirmed by the device.
- b) **Silent Mode.** Silences the audible cues that normally sound during mode switching, however the vibrational cues will still be active.

Notes: i) Enabling **Silent Mode** does not silence the audible cues emitted when **Training Mode** is active.
ii) Enabling **Silent Mode** does not silence any critical warning sounds or vibrations when a critical warning is active.


The selection is retained by the device even after the device is powered OFF and ON. Tap **ENABLE/DISABLE** button then **CONFIRM AND CONTINUE** button on the **Safety alert** screen to enable/disable **Silent mode**. Button illuminates teal (blue-green) when selected operating mode has been confirmed by the device.

- c) **Max Hydraulic Pressure.** The maximum hydraulic pressure readout is a monitoring tool to help the practitioner determine if the user is or has been over-stressing the device. Indicates the maximum pressure reached during any one operating session up to maximum 200 bar, resets at power-off.
- d) **Device step count since last reset.** Real-time indication of total number of steps taken by the device since last reset.
- e) **BASIC SETUP.** Press to return to the **BASIC SETUP** screen.

Tectus programming is completed by returning to the **BASIC SETUP** screen by tapping **BASIC SETUP** then tapping **COMPLETE**.

12.2.3 SERVICING Screen

This screen is only available to practitioners who have the relevant authorization granted when registering to use the *Blatchford* app. This is usually granted as a result of completing the prerequisite certification training provided by Blatchford. Authorized practitioners can access additional app functions by choosing the *SERVICING* menu option on the main menu screen.

- a) *Key Fob Programming*. By pressing the **CLEAR KEY FOBS** button all previous key fob button assignments are cleared.
 **WARNING. All key fobs paired to the device will be affected.**
- b) *Pair Key Fob Buttons*. This option allows for the assignment or reassignment of operating modes to key fob buttons. On the app screen press the mode you wish to assign to a key fob button (**WALK, SIT/STAIRS, FREE, LOCK**) then press the corresponding button on the key fob. Successful assignment is confirmed by the app button color turning teal (blue-green). See Section 13.1 *Assigning Operating Modes to Key Fob Buttons* for details.
- c) *Selected Leg*. The leg on which the Tectus KAFO is worn is assigned here. If the selected leg is changed, power the Tectus device off and on again to register the change.
- d) **COMPLETE**. After completing all required service actions click **COMPLETE** to return to main menu.

12.3 Calibrating the Operating Modes for the Specific User

Carry out calibration of the operating modes sequentially following order given in this sections 12.3.1 to 12.3.11.

12.3.1 Standing Calibration in Walk Mode

Always use parallel bars for standing calibration.

- Switch the device to **Walk** mode on the app.
- The user should stand up and assume a natural comfortable position; legs aligned and arms to the side, not holding onto the bars.
- Ensure the Tectus KAFO alignment is accurate by confirming:
 - The ankle range and shoe heel pitch is accommodated (add heel wedges if required).
 - The lower leg has approximately a 10 degree anterior tilt.
 - The Tectus mechanical joint location is accurate (this will be posterior and distal to the anatomical compromise knee centre and as prescribed).
- Perform standing calibration by pressing **CALIBRATE** on the **BASIC** setup screen on the app.
- The app confirms completion of standing calibration by **changes saved** displaying at the bottom of the screen.
- After the calibration is completed the user should be able to:
 - Stand in a normal position.
 - Feel comfortable and safe when standing.

12.3.2 Sit Resistance Calibration

Walking rails or a chair with arms should be used for support during this procedure.

- Switch the device to **Sit** mode on the app.
- Ask user to stand up and slowly sit back down.
- Confirm if they can feel the resistance of the unit.
- The device should provide sufficient support but should not prevent sitting down.
- Repeat the procedure and adjust the **Sitting Resistance** on the **BASIC** setup screen in small increments, sitting down, standing back up and re-selecting **Sit** mode each time until the patient is satisfied with the result.

12.3.3 Walk Mode Stance Release Angle Calibration (Training Mode Option)

Always use parallel bars for walking calibration.

- Switch the device to **Walk** mode on the app.
- With the user standing ensure even weight bearing on the lower limbs and good posture.
- The user should take a few steps in the default **Stance Release Angle** setting.
- Ensure the user is able to progress through each stage of the stance phase with initial contact of the heel, progression through mid stance and loading the forefoot throughout third rocker in gait. If not;
- Adjust the **Stance Release Angle** slightly up or down as necessary and ask the user to take a few steps again. Note any improvement in gait i.e. when the knee releases in stance phase.
- Continue adjustment until the user is able to progress through each stage of the stance phase and feels safe and comfortable to walk unassisted.
- If satisfactory knee release cannot be achieved by adjusting **Stance Release Angle** alone, refer to the following **Refining Walk Mode Calibration** table.

Refining Walk Mode Calibration

Advise the user if pressure points are observed to stop using Tectus KAFO and to inform their practitioner as soon as possible.

Symptom	Solution
Tectus module reaching and staying in a high yield position and knee not releasing in late stance or swing phase.	1. Check the patient has a step-over-step gait with a full step cycle and loading into late stance to ensure the thigh inclines. It may take some repetition for the patient to adapt to this pattern. 2. Make the Stance Release Angle less negative and repeat the standing calibration procedure. 3. If this does not work, try reducing the Swing Adjustment and repeat the standing calibration procedure. 4. If this still does not work, try reducing the Pressure Release Threshold.
Tectus module is slow in swing phase and patient is catching their toe or stumbling.	1. Check the set-up alignment with the ankle joint set with a 10 to 12 mm heel wedge/pitch so there is approximately a 10 degree anterior tilt of the lower leg. This will ensure you are not pitched with the foot/ankle complex pointing downwards. 2. Reduce the hydraulic swing resistance by reducing the Swing Adjustment setting on the app (please refer to Section 12.3.4). 3. Increase spring stiffness by adjusting the knee extension assist spring. Please refer to tables in Section 12.3.5 and 12.3.6.
Tectus module is quick in swing phase with excessive heel rise and/or terminal swing.	1. Adjust the knee extension assist spring. Please refer to tables in Section 12.3.5 and 12.3.6.
Tectus module is slow in swing phase with insufficient heel rise.	1. Adjust the knee extension assist spring. Please refer to tables in Section 12.3.5 and 12.3.6.
Tectus module does not reach full extension. Note... This should not happen if there is sufficient control to initiate swing phase of the effected limb.	1. If the patient loses balance and shortens the swing phase to then bear weight on the Tectus module early, the device will provide a high yield resistance for stumble recovery.
Sitting and stairs stability.	1. Please ensure sitting and stairs calibration guidance is followed, see Sections 12.3.2 and 12.3.8 respectively. If the patient does not feel stable during sitting or stairs descent, then increase the associated Sitting or Stairs resistance on the app accordingly.

Training Mode Option

- **Training** mode can be enabled from the **ADVANCED SETUP** screen on the app.
- With **Training** mode enabled,
 - When in **Walk** mode the direction of walking is determined by the first step from standing. Two short beeps are emitted every time the device enters the **Standing State** of **Walk** mode.
 - One long beep is emitted when the device enters the **Backward Walk State** of **Walk** mode.
 - No beeps are emitted when the device enters **Forwards Walk State** of **Walk** mode.
- **Training** mode does not alter the operation of the knee module in any other way.
- Once the required walking confidence is reached the practitioner can disable **Training** mode.

12.3.4 Swing Adjustment—Tuning the Hydraulic Swing Resistance


The hydraulic system in **Walk** mode exerts resistance both during swing phase and stance phase. The app allows fine tuning of the swing resistance in motion. The value can be adjusted using **Swing Adjustment** setting on the **BASIC Programming** screen. (See Section 12.2.1.)

Gait Deviation Experienced During Walk	Swing Adjustment (in App)
Insufficient heel rise at initial swing phase.	Reduce hydraulic resistance.
Insufficient knee flexion at initial swing phase.	Reduce hydraulic resistance.
Excessive heel rise at initial swing phase.	Increase hydraulic resistance.
Excessive knee flexion movement at initial swing phase.	Increase hydraulic resistance.
Mid-swing and/or terminal swing moving too slow.	Reduce hydraulic resistance.
Mid-swing and/or terminal swing moving too fast.	Increase hydraulic resistance.
Catching toes on ground in swing phase.	Reduce hydraulic resistance. (After confirming that ankle joint position and heel pitch wedge height is correct.)

12.3.5 Tuning the Knee Extension Assist

The device is equipped with an internal spring to assist with knee extension during the swing phase of the gait. This feature provides extension assist and should not prevent knee flexion. Depending on spring rate and pre-load, flexion resistance in swing phase may be affected. The device is delivered with three springs rated according to their stiffness (see Section 12.3.6 *Spring Selection and Adjustment Guide*) and it comes with the stiffest spring factory fitted as standard (Spring 3, part number 925042). If required, it can be changed to one of a lower stiffness rating as required (see Section 17 *Changing the Extension Assist Spring*).

1


Adjuster Cover

Adjuster


Pry open the extension assist adjuster cover from the side opposite the hinge and rotate the cover to expose the adjuster feature.

2

End stops on either end of the adjustment range ensure that the adjuster is not removed and remains intact.

3


Clockwise
Increases
Assistance

Anticlockwise
Decreases
Assistance

Adjust the amount of knee extension assist using the 5 mm key provided. Never force the screw beyond its highest and lowest settings.

4

After completing the adjustments make sure the adjuster cover is put back into the original closed position.

12.3.6 Spring Selection and Adjustment Guide

The following guide is provided as an approximate indication which spring out of the three stiffness options should be used as a starting point of the tuning process.

The unit is delivered with all three springs. Spring 3 is fitted at the factory as default.

See Section 12.3.5 *Tuning the Knee Extension Assist* for description of tuning process and Section 17 *Changing the Extension Assist Spring* on how to change springs.

	Spring 1	Spring 2	Spring 3
Part Number	925040	925041	925042
Stiffness	Low	Medium	High
Activity Level	2, 3 and 4	2, 3 and 4	2, 3 and 4
Quadriceps Power	Grade 2 and above	Grade 1 and above	Grade 0 and above

Fine tuning of the amount of assistance the knee unit provides during swing phase can be done by:

- Adjusting spring tension/pre-load (see Section 12.3.5, box 3).
- Replacing it with a spring of different stiffness when further adjustment is not possible (see Section 17 *Changing the Extension Assist Spring*).

Gait Deviation Experienced During Walk	Spring Extension Assistance Adjustment
Insufficient heel rise at initial swing phase.	Reduce spring stiffness.
Insufficient knee flexion at initial swing phase.	Reduce spring stiffness.
Mid-swing and/or terminal swing moving too slow.	Increase spring stiffness.
Mid-swing and/or terminal swing moving too fast.	Reduce spring stiffness.
Catching toes on ground in swing phase.	Increase spring stiffness. (After confirming that ankle joint position and heel pitch wedge height is correct.)

12.3.7 Adjusting the Ankle Joint

The Tectus KAFO is fitted with an ankle supplied by a manufacturer other than Blatchford. For setting up and adjusting instructions, please consult the corresponding fitting instructions delivered with the ankle.

Blatchford recommends a 16 mm or 20 mm wide lateral ankle joint system that provides a dynamic and adjustable dorsiflexion and plantar flexion stop. The stops are set by the spring units.

The medial supporting ankle joint is to be free in motion and we recommend a 16 mm or 20 mm wide ankle joint is used with parallel placement to the lateral joint.

12.3.8 Stairs Mode Calibration (Stairs Descent)

Always use handrails for *Stairs* mode calibration.

- Set the ***Stairs Resistance*** to the same value as set for ***Sitting Resistance***.
- Ask the user to stand at the top of a flight of stairs.
- Ask the user to set the device to ***Stairs*** mode on the key fob or you set it to ***Stairs*** mode on the ***BASIC Setup*** screen.
- Holding onto the handrails the user should start stairs descent.
- When stepping down with the unaffected leg, the affected leg should pivot over the midpoint of the foot during every step.
- If resistance is too high the knee will not bend fast enough resulting in the patient leaning forward excessively. Reduce the ***Stairs Resistance*** and ask the user to retry stairs descent until resistance is satisfactory.
- If the resistance is too low the knee will bend too quickly resulting in the patient leaning backwards. Increase the ***Stairs Resistance*** and ask the user to retry stairs descent until resistance is satisfactory.

To ensure patient safety, always switch to *Lock* mode while adjustments are made to *Stairs Resistance*.

Note... Alternatively a short ramp or a single step of appropriate height can be used instead of a flight of stairs.

12.3.9 Confirming Key Fob Pairing and Functionality

- Press each key fob button in turn (***Walk***, ***Sit/Stairs***, ***Free*** and ***Lock***) and confirm the device selects the associated mode (or selects the step counter), and indicates the mode status change accordingly on the device on-board graphic display.

Note... Be aware that there are some mode changing sequences that are, for safety reasons, not allowed. Use only permitted mode changing commands for testing the outcome of key fob pairing. See Section 14.2 ***Safety While Switching Between User Modes***.

12.3.10 Confirming Silent Mode Availability

- Select ***Silent*** mode on the ***ADVANCED*** screen of the programming app by pressing the ***ENABLE/DISABLE*** button.
- With ***Silent*** mode on, test that no sounds are emitted when button operations and mode changes are made.
- After confirmation ***Silent*** mode can be left enabled or disabled according to user preference.

Note... Audible warnings are not disabled by activation of Silent mode and will still sound.

12.3.11 Registering Completion of Calibration and Setup

When you have completed all calibration operations press ***PROGRAMMING COMPLETE***. Button illuminates green when pressed and acts as a reminder that you have completed all calibration and setup stages.

Note... All settings on the BASIC SETUP screen need to have been changed before programming can be successfully completed.

13 Key Fob Programming

Tectus is delivered with two pre-programmed and paired key fobs ready for use.
Additional key fobs can be ordered and paired to the device using the *Blatchford* app when:

- Use of a third key fob is required by the user.
- One or more of the key fobs need replacement.

Note... a maximum of three key fobs can be paired to one Tectus module at a time.

Prerequisites:

- Service Level login authorisation for the app.
- The presence of both the Tectus Module and the key fob in close proximity of the programming device with *Bluetooth®* connection established between app and Tectus module.

13.1 Assigning Operating Modes to Key Fob Buttons



1. Navigate to the *SERVICING* screen.
2. Press the *WALK* button on the app.
3. Press the *WALK* button on the key fob. When assigned, *WALK* on the app changes to green.
4. In the same way apply the *SIT/STAIRS*, *FREE*, *LOCK* functions to buttons *Sit/Stairs*, *Free*, *Lock* on the key fob.
5. Test all key fob buttons' functionality by pressing all of them in turn and confirming they achieve the required mode change on the device.

These changes can be seen:

- On the on-board graphic display of the device where a mode indicator symbol appears.

- At the bottom of the *BASIC* setup screen of the app where the corresponding mode button turns blue.

13.2 Clearing Key Fob Button Assignments

The key fob button assignments are cleared by pressing the *CLEAR* button on the app screen.

 **WARNING.** All key fobs paired to the device will be affected.

14 Tectus Knee Module Operation

14.1 User Modes

The device has the following user selectable operating modes:

- *Walk*
- *Sit*
- *Stairs*
- *Lock*
- *Free*
- *Silent*

The device may enter automatically the following modes:

- *Critical Service*
- *Safe*

The device has the following practitioner-only selectable operating mode:

- *Training*

14.2 Safety While Switching Between User Modes

14.2.1 Default User Mode

The device after start-up, by default, switches to **Free** mode if the user is seated or **Lock** state of **Walk** mode if the user is standing. During start-up **Free** mode is only intended to be used for donning the device. Another mode (typically **Walk** mode) should be selected before standing after donning the device.

14.2.2 Safety Restrictions to User Mode Switching

The device will not allow switching to **Free** mode when the user is walking, descending stairs or standing with the knee locked. This is for safety reasons as the consequent sudden loss of support poses a considerable risk of falling. Therefore the following mode switching combinations will not be executed by the device:


Walk mode to *Free* mode


Stairs mode to *Sit* mode

Stairs mode to *Free* mode


Sit mode to *Stairs* mode

Lock mode to *Free* mode

Free mode to *Stairs* mode

14.3 Device Over-Stress Warning

Very rarely, and under special circumstances, the device may emit a rapid beeping audio and vibration alarm, for up to 15 seconds. This indicates the device is over-stressed, through not being used as intended, or due to an internal fault.

When Stress Warning Occurs

If the user (or the device itself) takes successful remedial action to remove the over-stress condition before the end of the 15 seconds alarm period, then the alarms stop and the device resumes normal operation.

If the user (or the device itself) does not take successful remedial action to remove the over-stress condition before the end of the 15 seconds alarm period, then the device interprets this as a fault situation and switches automatically to Safety mode.



Occurrence and Avoidance of Over-stressing the Device



The following table summarizes the scenarios in which the device might be over-stressed and the possible remedial actions.

Alarm Occurrence	Reason	Remedial Action
During Mode Switching	The user may have attempted to change between user modes when the previous activity had not been fully completed and the hydraulic system was not ready to adjust to the new settings.	The user should try to reduce the load on the device until it can execute the mode switching command. The rapid beeping should stop immediately to indicate that normal operation has been resumed. Advise the user to always allow sufficient time for the device to respond to mode switching commands.
During Sitting Down	During the process of sitting down the device automatically switches to a lower resistance after about a 45° of movement. If the user attempts to sit down too quickly and forcefully, it over-stresses the hydraulic system during mode switching.	The user should try to reduce the load on the device until it is able to switch to Free mode. Advise the user on the correct way and speed of sitting down while wearing the device.
In Case of a Malfunction	If the remedial action does not fix the problem within the 15 seconds alarm period at the end of this period, the device will assume there is a malfunction and revert to Safety mode and the critical service indicator will appear on the device screen.	Please refer to Section 5.1 Critical Service Indicator on how to deal with a malfunction and for the description of Safety mode.

14.4 Activities and Associated User Modes

Activity	Mode Activated and Displayed	Function and Procedure Guidance
Walking 	Walk 	Automatically provides support during stance and release during swing. Although not microprocessor controlled, knee extension assist is always active.
Standing 	Walk  Lock 	If the user stands still while in Walk mode, the user will be supported. From standing, the device will be alert to any movement, such as a forward or backward step and provide support accordingly. To provide appropriate support, using Lock mode is recommended when standing for an extended period of time (e.g., queueing) or in an environment where the user's balance can be compromised (e.g., standing in a crowd). Note... Lock mode is a fixed mode, and the device does not monitor for movement. Lock mode needs to be set and released manually.
Backstep 	Walk 	If the user steps backwards or leans backwards while in Walk mode, the device goes into a locked state called Backstep. This gives the user support while they rebalance. Note... Backstep is a state within Walk mode and is not the same as the user-selectable Lock mode. If Training mode is enabled, the device indicates it has gone into this Backstep state by emitting one long high-tone beep and one vibration. The device reverts to standing Lock state only after the user has returned to and remained in the standing position for at least two seconds. If training mode is enabled, this is indicated by the device emitting two beeps and vibrations.
Lock 	Lock 	Knee locked against flexion but free to extend until maximum extension is achieved. Lock mode will remain active until another mode is selected or until the mechanical override button is pressed.
Free 	Free 	The device provides minimal resistance in flexion and minimum amount extension assist resistance in extension. Free mode will remain active until another mode is selected.  This mode does not provide sufficient support for standing.

Activity	Mode Activated and Displayed	Function and Procedure Guidance
Stand-to-Sit 		Fixed resistance at knee to provide support while moving to seated position, the geometry of the module provides deep yield functionality. To sit, the user must manually activate Sit mode via mode switching buttons on the device or on the key fob. The user should be aware of the surroundings before sitting down. Both legs should be loaded equally when sitting down. The user should bend forward at the waist, push buttocks rearward, place hands on armrests (if available) and bend the legs, sitting down in a smooth motion.
Sitting 	 	In Sit mode after 2 seconds of achieving a resting position the device automatically switches to Free mode to allow, if necessary, adjustment of limb position while sitting. Note... The device will start up in this Free mode when start-up is initiated while the user is in a seated position. While in Sit mode the device will detect when the user attempts to stand and locks against flexion as they rise, so they cannot sit back down again until either the leg is fully extended or the override button is pressed.
Sit-to-Walk 		When in Sit mode the device detects when the user attempts to stand up and start walking and enters into Walk mode automatically. To stand up, the user's feet should be placed evenly, they should move their trunk forward and stand up with the support of their arms. The user should try to load weight through the orthosis and begin walking.  If the user does not complete the motion and attempts to sit back down again the device will lock against flexion.
Stairs Ascent 		 For safety, the user must switch manually to Lock mode before starting to ascend stairs.  Where possible, handrails must always be used. When Lock mode is used the knee will not yield (although it can extend), so a "step-to" gait should be used leading with the sound limb going up stairs.

Activity	Mode Activated and Displayed	Function and Procedure Guidance
Stairs Decent 	<i>Stairs</i>  <i>Lock</i> 	⚠ In <i>Stairs</i> mode, the user must not attempt to descend stairs using leg-over-leg gait until they have had adequate physiotherapy and completed supervised training signed off by their practitioner! ⚠ Ensure <i>Stairs</i> mode is actually selected before descending stairs. ⚠ The user should always approach stairs with caution and be aware of their surroundings. ⚠ Where possible, handrails must always be used. ⚠ For stairs decent, either <i>Stairs</i> mode or <i>Lock</i> mode must be enabled before starting descent and once selected, the user must not change it to another mode until the stairs decent is complete. <i>Stairs</i> mode provides appropriate support for leg-over-leg descent with programmed flexion resistance combined with appropriate extension bias. Note... The user should place only the posterior half of their foot on the step then load orthosis and push it in to knee flexion, step down with the sound limb and allow knee flexion on the orthotic limb. When <i>Lock</i> mode is used the knee will not yield, it is locked in full extension, so a "step-to" gait should be used, leading with the locked fitted limb going down the stairs.
Entering Vehicle 	<i>Sit</i> 	⚠ The user must always approach vehicles with caution and be aware of their surroundings. To sit into the vehicle, the user should activate <i>Sit</i> mode. The device will enter <i>Free State</i> of <i>Sit</i> mode when the thigh is horizontal. (Not identical to <i>Free</i> mode.) This allows the user to manoeuvre the limb into a comfortable and safe position. Note... If the thigh angle is not horizontal e.g. due to the angle of the vehicle seating surface, then the device will enter <i>Sit-to-Stand State</i> of <i>Sit</i> mode and can only be positioned by using the mechanical override button.
Exiting Vehicle 	<i>Sit</i> 	Before exiting the vehicle the user should make sure the device is in <i>Sit</i> mode so that when they exit the vehicle they are supported as they move from sitting to standing and is ready to switch automatically to <i>Walk</i> mode.

Activity	Mode Activated and Displayed	Function and Procedure Guidance
Driving 	Free 	⚠ The user must not operate the device controls or key fob while driving. ⚠ Warnings in this document and any national guidance must be followed when considering driving while wearing a Tectus KAFO. ⚠ Only the unaffected leg may be used to operate controls of a car. For driving, the device should be manually switched to Free mode by the user via the key fob or the device buttons.
Safe Mode		The device will revert into Safe mode every time a malfunction or critically low battery level is detected, or when the service mode warning period has expired. In Safe mode the device will switch to high yielding knee flexion resistance. See Section 5.1 <i>Critical Service Indicator</i> and Section 7.3.2 <i>Low Battery Behavior</i> for more details.
Critical Service		
Training	Training	Note... Training mode can only be enabled through practitioner (clinician) level access to the app. In Training mode, aural and vibrational cues are provided to indicate when an operating state change within Walk mode has been made to help the user learn how to use the device. For example, when the user makes the appropriate first step movement from stationary to activate Backwards Walk state, the device emits one long beep and vibrates to confirm Backwards Walk state has been selected. Similarly, if the user stops walking to stand still, the device emits two short beeps and vibrates to confirm Lock state has been selected.

14.5 Operating Actions and Device Indications

Device and Key Fob Button Operation

The device is manually operated by pressing either buttons A and/or B on the device, or by pressing buttons on one of the key fobs. The button or button combination pressed, together with the hold duration (Short, Long or Very Long) of the press determines which operation is selected.

Note... Key fob buttons have only two options; Short or Long press.

To help the user distinguish between the three different press durations, for each press duration the device emits a corresponding unique pattern of beeps and vibrations, and displays a corresponding icon on the screen, as shown in the following table:

Name	Duration	Timing Beeps with Vibrations	On-screen Icon	Symbol Used in Document
Short Press	Max. 2 seconds			
Long Press	2-4 seconds			
Very Long Press	Longer than 4 seconds			


Device Indications

The device will acknowledge that the requested operation, e.g. change of mode, has been registered by the device by emitting audible and/or vibrational cues as shown in the following table:

Type of Indication		Symbol Used in Document	
Audible (Beep)		Short	Long
	High Tone		
	Medium Tone		
	Low Tone		
Vibration			
On Screen		See detailed description in Section 14.4, 14.6 and 14.7.	
ON/OFF Switch LED Light		See detailed description in Section 10.	

14.6 Mode Switching on the Device Using Tectus Module Buttons A and/or B


Once start-up routine has successfully completed, the device can be controlled with button combinations on the Tectus Knee Module or Key Fob.

Modes	Initiation (On Device)	Timing Aid		Device Response		
		On Screen	Beep & Vibrate	Beep Signal	Vibrate Signal	Displayed Symbol
<i>Walk</i>	Short Press button A 		1 	2xLow tone short 	2xShort 	
<i>Sit</i>	Long Press button B 		1+2 	1xMedium tone short 	1xShort 	
<i>Stairs</i>	Short Press button B 		1 	2xMedium tone short 	2xShort 	
<i>Free</i>	Short Press button A&B 		1 	1xMedium tone long 	1xLong 	
<i>Lock</i>	Long Press button A 		1+2 	2xMedium tone long 	2xLong 	
<i>Silent</i>	Very Long Press button A 		1+2+3 	2xMedium tone short 	2xShort 	

Note... Do not press button B for more than 4 seconds or you may accidentally switch on **Bluetooth®** as indicated by the power button illuminating flashing blue and emitting three low tones and vibrations. To switch off **Bluetooth®**, power off the device and power on again using the ON/OFF switch (see Section 10), or press button B for longer than 4 seconds.

14.7 Mode Switching on the Device Using the Key Fob

Modes	Initiation (On Fob)	Timing Aid		Device Response		
		On Screen	Beep & Vibrate	Beep Signal	Vibrate Signal	Displayed Symbol
Walk	Short Press Walk button			2xLow tone short 	2xShort 	
Sit	Short Press Sit \ Stairs Button			1xMed. tone short 	1xShort 	
Stairs	Long Press Sit \ Stairs Button		1+2 	2xMed. tone short 	2xShort 	
Free	Short Press Free button			1xMed. tone long 	1xLong 	
Lock	Short Press Lock button			2xMed. tone short 	2xShort 	

14.8 Additional Key Fob Functions

Modes	Initiation (On Fob)	Timing Aid		Device Response		
		On Screen	Beep & Vibrate	Beep Signal	Vibrate Signal	Displayed Symbol
Step Counter	Long Press Walk button		1+2 	-	-	
Connection Code	Long Press Lock Button		1+2 	-	-	

14.9 Mechanical Override Operation

The knee can be flexed in any mode by pressing the manual override button located on the top surface of the control module. The button must be pressed and held down fully until the desired knee flexion angle has been achieved.

Do not use any object (e.g. screwdriver, pen) to press the override button.


This operation should only be used in exceptional circumstances (e.g. device malfunction, battery depletion).

Once the manual override button is released the knee flexion reverts to the same resistance it had before the manual override button was pressed.

Note... The hydraulic module has been designed to always permit free leg extension.

Be aware of finger trap hazard at all times when flexing the device.

15 Troubleshooting

Advise the user if pressure points are observed to stop using Tectus KAFO and to inform their practitioner as soon as possible.

Symptom	Solution
Tectus module beeping when walking.	1. Please ensure the Training mode has been turned off and check the graphic display for any warning symbols. 2. Unrecognized beeps or any unusual noise must be reported to the support team or practitioner.
Warning symbols on the Tectus module on-board graphic display.	1. Please refer to Section 5 for details and instructions.
Red LED flashing.	1. If accompanied by one high tone beep, it indicates a critically low battery level or a critical error. 2. If accompanied by two high tone beeps, it indicates a low battery level. In both cases, please, refer to Section 7.3.2 Low Battery Behavior and 5.1 Critical Service Indicator. If the error persists, please, contact your practitioner or service provider.
Slight play felt by user when device is in Lock mode.	If no obvious reason is detected by visual inspection, attempt Hydraulic System Optimization procedure. If fault detected or Hydraulic System Optimization does not remedy symptom, contact Blatchford or your service provider.

15.1 Hydraulic System Optimization

The slight play experienced by the user does not constitute a fault but if it causes inconvenience this Hydraulic System Optimization procedure can be performed.

Ask the user to remove the KAFO for the duration of the procedure. Then activate Lock mode on the device and repeat the following actions 3-5 times:

1. Press and hold down the manual override button.
2. Slowly flex the KAFO until full flexion is achieved.
3. Release the manual override button.
4. Let the KAFO extend slowly until full extension is achieved.

After completing the Hydraulic System Optimization, instruct the user to don the KAFO and confirm that the previously experienced symptoms have been resolved.

If Hydraulic System Optimization needs to be applied frequently, the user should stop using the device and Blatchford or your service provider should be contacted.

16 Disassembly/Assembly Instructions

These instructions are for the trained practitioner and/or technician only.

This section describes how to remove and/or replace the knee module, battery pack or follower joint.


16.1 Disassembly

The disassembly process should generally follow the assembly steps in reverse order.

Special care should be given during removal of the battery pack connecting cable. There are two versions of cable routing the manufacturer might have chosen to apply:

- a) Routed on the outside of the thigh shell and secured by P-clips.
- b) Routed on the inside held in place by the inner lining of the thigh shell.

The connecting cable has to be detached using one of the following two methods (1a or 1b), before the battery pack can be disconnected:

16.2 Assembly

1

Knee Module

Proximal Plate Cover

1

Tectus Retainer Screws (x2)
M4 x 6 mm Hex

2

Distal Plate Cover

Tectus Cover Retainers

Tectus Cover Placement Screws (x8)
M3 x 8 mm Hex

Proximal Follower Cover

3

Follower Joint Cover Retainers

Follower Joint Retainer Screws (x2)
M4 x 8 mm Hex

4

5

Distal Follower Cover

6

Follower Joint Cover Placement Screws (x6)
M3 x 6 mm Hex

A new Tectus module and follower joint is delivered loosely assembled; all covers, retainers, retainer screws and cover placement screws are fitted only hand tight.


Remove all retainers, retainer screws, placement screws and covers and keep them for later.

Notice the slightly different sizes required for the Tectus joint and the follower joint.


- 2** Fit the knee module onto the proximal and distal knee side-bars on the lateral side of the Tectus KAFO.

Fit the follower joint onto the proximal and distal side-bars on the medial side.


The rotational angle of the control module of the device can be adjusted by loosening the clamp screw on the proximal plate. After any adjustment, tighten the screw to 3 Nm.

- 3** Fit all side-bar connecting screws (M5×10 mm) and use a T25 Torx torque driver, set to a maximum torque 5 Nm.


- 4** Fit the distal plate cover and the follower cover. Then reassemble and tighten cover retainers and retainer screws to 3 Nm torque using a 3 mm Allen key.

The side-bar connecting and retainer screws can be removed and refitted a maximum of 5 times. After that, a set of new screws should be used.

5


Fit all Tectus cover placement screws (×8) and all follower joint cover placement screws (×6). Tighten all screws only finger tight using a 2 mm Allen key. Do not overtighten so the cover would distort.

6


Depending on the design chosen during KAFO manufacturing, route the connecting cable from the actuator to the battery pack through:

- the refitted P clips.
(Only tighten the retaining screws loosely at this stage.)
- or the thigh shell internal cable channel.
(Do not refit the inner lining yet.)

7

Before connecting the cable to the battery pack, note the orientation of the connector required for correct alignment and engagement.


To prevent damaging the connector pins, always connect the battery cable to a battery pack before connecting the battery pack to the to the upper thigh enclosure.


- 8** Connect the knee power cable to the battery pack via the connector and retention nut.

Ensure the connectors are aligned correctly and then engage them. Then, holding only the battery pack and connector retention nut, secure them by tightening the retention nut.

Do not rotate the cable or connector.


- 9** Attach battery pack to Tectus KAFO as shown and secure it using the three M2.5 screws supplied. Tighten the screws to 0.5 Nm torque.


- 10a** Make sure that there are no bends in the cable with a radius of less than 15 mm.

Then tighten the P clip retaining screws hand tight. Apply Loctite 243, or equivalent threadlocker, to secure them.

- 10b** Refit inner lining of thigh shell using adhesive specified by manufacturer.


17 Changing the Extension Assist Spring


Before disassembling the extension assist spring assembly, make sure the Tectus module is fully extended to reduce the hydraulic pressure on the spring. Failure to do so could result in forceful ejection of the spring and the adjuster cap on release.

- 1** Pry open the extension assist adjuster cover from the side opposite the hinge. Using the 4-pin installation tool fully unwind the spring adjuster cap. The 4-pin installation tool should be used in conjunction with a 13 mm socket.


Note... Align slit pin with recess in installation tool, if necessary, by turning the adjuster screw with a 5 mm hex key.

4-pin
Installation Tool
(829061)


2


Remove the adjuster cap, the compression spring and the spring coupling.

3

Replace the compression spring with one that has the required stiffness and reassemble the spring coupling, the compression spring and the adjuster cap.

4

Using the 4-pin installation tool and a calibrated torque driver, tighten the adjuster cap to 4 Nm then close adjuster cover.

18 Transport and Storage


Transport

1. For transport purposes the unit must be switched off, and packed securely so that
 - a) it cannot be inadvertently switched on and become operative while in transit,
 - b) the Tectus KAFO/device is thoroughly protected from damage in transit.
2. If the device is shipped on its own it should be placed into protective foam S-0000481 and packaged securely in a suitable box.
3. If the device is shipped built into a Tectus KAFO it should be put into Lock mode in an extended position. The device should be covered with protective foam S-0000481. The Tectus KAFO should be placed into its own limb bag (when available) and securely packaged in a suitable box.
4. If a device is to be returned, contact Blatchford prior to sending if any obvious physical damage has occurred to the battery pack.

Storage

1. It is recommended to store the device upright with the device put into **Lock** mode before switching off.
2. If the Tectus KAFO is to be stored upright then secure the device to prevent it falling over.
3. For extra protection during storage, it is recommended that the device is covered with protective foam S-0000481.

Note... Protective foam S-0000481 is part of the original packaging.

19 Technical Data

Operating and Storage Temperature Range: -10 °C to 50 °C
(14 °F to 122 °F)

Charging Temperature Range: 0 °C to 35 °C
(32 °F to 95 °F)

Component Weights:

- Tectus Module (with joint covers and screws) 732 g (1 lb 9.8 oz)
- Battery Module 107 g (3.8 oz)
- Four side bars and screws 330 g (11.6 oz)
- Follower Joint (with joint covers and screws) 185 g (6.5 oz)

Activity Level: (1), 2, 3, (4)

Maximum Combined User and Load Weight: 100 kg (220 lb)

Minimum User Weight: 50 kg (110 lb)

Range of Motion: 110°

Build Height: See *Fitting Length* diagram on page 57.

Device Battery

Battery Type: Rechargeable Ni-MH 950 mAh 6V

Battery Life: Variable, but max. 18 hrs from full charge. Charge daily.

Battery Charger

Charging time to full charge: 2 hours

Input Voltage: 100V to 240V AC

Frequency: 50 Hz to 60 Hz




Max. Output Voltage: 21 V DC

Max. Output Current: 0.35 A

Key Fob Battery

Battery Type: CR2032, Disposable Lithium 3V

Fitting Length

Maximum cable length from actuator housing to battery pack assembly: approx. 250 mm

Appendix 1 Installing the Blatchford Programming App

A1-1 Android Devices

1. Open the **Google Play Store** app on your Android programming device.
2. Search for the **Blatchford App** in the search bar at the top of the screen.
3. Tap on the app name and verify the information of the distributor is 'Blatchford Products Ltd'.
4. Tap the **Install** button on the app.
This will begin the process of downloading and installing the **Blatchford App** on your programming device.
5. When successfully installed, return to your **Home Screen** -> **App menu** and you should have the **Blatchford App** as an option to open (see Fig. A1-1).

System Requirements

Android 11, Red Velvet Cake (API level 30) and above.


A1-2 Apple iOS Devices

1. Open the **App Store** app on your Apple programming device.
2. Tap the search icon at the bottom right of the screen.
3. Search for the **Blatchford App** in the search bar at the top of the screen.
4. Tap on the app name and verify the information of the distributor is 'Blatchford Products Ltd'.
5. Tap the **GET** button on the **Blatchford App**.
This will begin the process of downloading and installing the **Blatchford App** on your programming device.
6. When successfully installed, return to your **Home Screen** -> **App menu** and you should have the **Blatchford App** as an option to open (see Fig. A1-1).


System Requirements

iOS 14.0 and above.


Appendix 2 Clinician Registration on Blatchford App


1. On the *Blatchford Welcome* screen, select **CREATE AN ACCOUNT**.


2. To register as a clinician, select *Register standard Blatchford account* under **For professionals**.


3. Enter the e-mail address you regularly use to communicate with Blatchford. Click on **SEND VERIFICATION CODE**.


4. Check the inbox of your e-mail account. Enter the verification code you received and click **VERIFY CODE**.


5. After successful verification of your e-mail address, click **NEXT**.


6. Complete the registration form.

7. Click **CREATE** to finalize your registration.


8. Select the preferred language for the application.

9. Select the region where your main practice is located.

10. After reading all app-related warnings and permission requests, click **ACCEPT AND CONTINUE**.

11. The app opens at your **Blatchford** home page and is ready for use.

Appendix 3 Gait Visualizer

The *Blatchford* Gait Visualizer app function provides real-time visualization of the gait patterns of the user and the parameters monitored by the device.

The Gait Visualizer can record, simultaneously with a video image, the following parameters:

1. *Thigh Angle*


Similarly to the Thigh Angle dial on the main programming screen, this trace shows in degrees how the thigh angle changes during an activity. Its value can vary between -15° and 15°.

2. *Pressure*


This trace shows how the pressure is changing during an activity. Its value varies between 0 and 200 bars.

3. *Motor Position*


This trace shows values corresponding to the actuator motor position during an activity. Its value varies in between 1000 (full lock) and 2000 (free swing).


1. After expanding the main menu, choose the *Gait Visualizer* option. **a**


2. The main *Gait Visualizer* page gives the option to record (tap *START A NEW RECORDING* **b**) or review a session (tap *REVIEW* **c**).


3. Several charts can be opened simultaneously (**1**, **2**, ...) by tapping the *ADD CHART* button **d**.


4. Tap the **Add Trace** **e** button and choose from the **Available Traces** **f** list the parameter or parameters to be displayed on a chart (max. 3).


5. The video can be started or stopped by tapping the middle of the image **g**. The chosen parameter traces will be displayed in real time.

6. Adjusting the start/stop sliders **h** on the video editing strip **h** enables review of a selected segment of the recorded session.

7. The X axis always shows the elapsed time, while the Y axis is always scaled to fit the largest value displayed. (Indicated by arrow.)

8. Press **←** to return to the main **Gait Visualizer** page. The main menu reappears **i** and allows access to other app features.

Manufacturer’s Declarations and Legal Information

Electromagnetic Environment

The device is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Electromagnetic Emissions

Interference Measurements	Compliance	Electromagnetic Environment Directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighboring electronic devices is unlikely.

Electromagnetic Interference Immunity

Phenomenon	EMC Basic Standard or Test Procedure	Interference Immunity Test Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
High-frequency electromagnetic fields	IEC 61000-4-3	10V/m 80 MHz to 6 GHz 80% AM at 400 Hz 100% AM at 1 kHz
	FCC CFR 47 part 15 (b)	30 MHz to 12.5 GHz
	ANSI C63.4:2014	
	AIM Standard 7351731 (Rev 2.00) 02:2017 Radiated RF Immunity	433.92 MHz, 3 V/m 860-960 MHz, 54 V/m 2.45 GHz, 54 V/m
	AIM Standard 7351731 (Rev 2.00) 02:2017 Magnetic Field Immunity	134.2 kHz, 65 A/m 13.56 MHz, 12 A/m
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

Bluetooth® Module

Part No.	BGM220PC22HNA2
FCC ID	QOQ-GM220P
IC	5123A-GM220P
Frequency	2.4 GHz
Effective Radiated Power	8 dBm

FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth®** module: FCC ID: QQQ-GM220P IC: 5123A-GM220P.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342

800-548-3534 (toll free) | 937-291-3636

customerservice@blatchfordus.com

RSS Compliance Statement

This device complies with RSS 210 of ISED Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices.

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 36 months.

A routine inspection of the device is required after 2 weeks, 3 months, 6 months, 12 months and 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Tectus and Blatchford are registered trademarks of Blatchford Products Limited.

The **Bluetooth®** word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under License. Other trademarks and trade names are those of their respective owners.

Android™ trademark is a registered trademark of Google LLC. All rights reserved.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

The MRC Muscle Scale is used with the permission of the Medical Research Council, © Crown Copyright.

TORX is a registered trademark of Acument Intellectual Properties LLC in the United States and other countries.

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09:54 Loading...
 Confirm New Password
 What is this?
 Display Name
 What is this?
 Given Name
 What is this?
 Surname
 What is this?
 CREATE →

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09:15
 Blatchford
 Select your preferred language
 English
 CONTINUE →

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09:16
 Blatchford
 Select your region
 EUROPE →
 UNITED →
 UNITED →

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Blatchford
 Thanks, you're setup and ready to use the Blatchford app.
 The features in this app must only be used by clinicians who have completed a Blatchford approved training course.
 Using this app without the Blatchford certification will void the product warranty.
 We also require Bluetooth and location permission to configure any of your Blatchford devices.
 Accept and →

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Blatchford
 Hi
 CONNECT TO A NEW DEVICE
 PRODUCT DOCUMENTATION
 REQUEST PRODUCT ACCESS
 Home Search Settings

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L'échelle musculaire MRC est utilisée avec l'autorisation du Medical Research Council, © Crown Copyright.
Die MRC-Muskel-Skala wird mit Genehmigung des Medical Research Council verwendet, © Crown Copyright.
La scala muscolare MRC viene utilizzata con il permesso del Medical Research Council, © Crown Copyright.
De MRC-spierschaal wordt gebruikt met toestemming van de Medical Research Council, © Crown Copyright. TORX is een gedeponiseerd handelsmerk van Acument Intellectual Properties LLC in de Verenigde Staten en andere landen.
A Escala de avaliação da força muscular do MRC é utilizada com a autorização do Medical Research Council, © Crown Copyright.
© Blatchford Products Limited 2026. All rights reserved.