



Paediatric Footwear

Issue 8

KSA

 **Mendivil**

DISTRIBUTED BY
Blatchford

Mendivil Paediatric Footwear



Mendivil Paediatric Footwear is exclusively distributed by Blatchford.

The range is available in sizes 16 to 45, with a wide variety of colour options and a choice of lace or reverse hook-and-loop fastenings. Two back heights are offered: UK and Continental.

You can also select from standard, anti-varus, or orthoses-accommodating width fittings, along with three different stiffener options.

Please note: not all features are available with every style, refer to the specific details in the catalogue for full specifications.

All Mendivil Paediatric footwear can be fully customised with floats, wedges, and raises to support your clinical needs.

Made to Order footwear allows for up to three upper last adaptations per last.

All footwear contained in this catalogue is distributed by Blatchford in accordance with supplier agreements. Mendivil is a trade name of Biomedical footwear who have a designated (UKRP) United Kingdom Responsible Person for product compliance.

Contents

	Page
Mendivil Paediatric Footwear	1
KSA Lazarus Boots and Bars	3
Made to Order Footwear	5
Specification Guide/Adaptations and Insoles	5-6
Colour Options	7
Ordering Information	8
Shoes	9-10
Boots	11-24
Measurement Charts	25-26
Made to Order Orthoses Boots	27-28
Orthoses Measurement Chart	29

KSA Lazarus Boots and Bars



Plantar Flexion Stop	
Size	Product Code
0000	LP.PFS.0000
000	LP.PFS.000
00	LP.PFS.00
0	LP.PFS.0
1	LP.PFS.1
2	LP.PFS.2
3	LP.PFS.3
4	LP.PFS.4
5	LP.PFS.5
6	LP.PFS.6
7	LP.PFS.7
8	LP.PFS.8

Lazarus Boot 2.0 - Grey

Size	Product Code Left	Product Code Right	Product Code Pair
6X0	LP.GR.S.000000.L	LP.GR.S.000000.R	LP.GR.S.000000
5X0	LP.GR.S.00000.L	LP.GR.S.00000.R	LP.GR.S.00000
0000	LP.GR.S.0000.L	LP.GR.S.0000.R	LP.GR.S.0000
000	LP.GR.S.000.L	LP.GR.S.000.R	LP.GR.S.000
00	LP.GR.S.00.L	LP.GR.S.00.R	LP.GR.S.00
0	LP.GR.S.0.L	LP.GR.S.0.R	LP.GR.S.0
1	LP.GR.S.1.L	LP.GR.S.1.R	LP.GR.S.1
2	LP.GR.S.2.L	LP.GR.S.2.R	LP.GR.S.2
3	LP.GR.S.3.L	LP.GR.S.3.R	LP.GR.S.3
4	LP.GR.S.4.L	LP.GR.S.4.R	LP.GR.S.4
5	LP.GR.S.5.L	LP.GR.S.5.R	LP.GR.S.5
6	LP.GR.S.6.L	LP.GR.S.6.R	LP.GR.S.6
7	LP.GR.S.7.L	LP.GR.S.7.R	LP.GR.S.7
8	LP.GR.S.8.L	LP.GR.S.8.R	LP.GR.S.8
9	LP.GR.S.9.L	LP.GR.S.9.R	LP.GR.S.9
10	LP.GR.S.10.L	LP.GR.S.10.R	LP.GR.S.10
11	LP.GR.S.11.L	LP.GR.S.11.R	LP.GR.S.11
12	LP.GR.S.12.L	LP.GR.S.12.R	LP.GR.S.12

Lazarus Boot 2.0 - Royal Blue

Size	Product Code Left	Product Code Right	Product Code Pair
6X0	LP.RB.S.000000.L	LP.RB.S.000000.R	LP.RB.S.000000
5X0	LP.RB.S.00000.L	LP.RB.S.00000.R	LP.RB.S.00000
0000	LP.RB.S.0000.L	LP.RB.S.0000.R	LP.RB.S.0000
000	LP.RB.S.000.L	LP.RB.S.000.R	LP.RB.S.000
00	LP.RB.S.00.L	LP.RB.S.00.R	LP.RB.S.00
0	LP.RB.S.0.L	LP.RB.S.0.R	LP.RB.S.0
1	LP.RB.S.1.L	LP.RB.S.1.R	LP.RB.S.1
2	LP.RB.S.2.L	LP.RB.S.2.R	LP.RB.S.2
3	LP.RB.S.3.L	LP.RB.S.3.R	LP.RB.S.3
4	LP.RB.S.4.L	LP.RB.S.4.R	LP.RB.S.4
5	LP.RB.S.5.L	LP.RB.S.5.R	LP.RB.S.5
6	LP.RB.S.6.L	LP.RB.S.6.R	LP.RB.S.6
7	LP.RB.S.7.L	LP.RB.S.7.R	LP.RB.S.7
8	LP.RB.S.8.L	LP.RB.S.8.R	LP.RB.S.8
9	LP.RB.S.9.L	LP.RB.S.9.R	LP.RB.S.9
10	LP.RB.S.10.L	LP.RB.S.10.R	LP.RB.S.10
11	LP.RB.S.11.L	LP.RB.S.11.R	LP.RB.S.11
12	LP.RB.S.12.L	LP.RB.S.12.R	LP.RB.S.12

Lazarus Boot 2.0 - Pink

Size	Product Code Left	Product Code Right	Product Code Pair
6X0	LP.PI.S.000000.L	LP.PI.S.000000.R	LP.PI.S.000000
5X0	LP.PI.S.00000.L	LP.PI.S.00000.R	LP.PI.S.00000
0000	LP.PI.S.0000.L	LP.PI.S.0000.R	LP.PI.S.0000
000	LP.PI.S.000.L	LP.PI.S.000.R	LP.PI.S.000
00	LP.PI.S.00.L	LP.PI.S.00.R	LP.PI.S.00
0	LP.PI.S.0.L	LP.PI.S.0.R	LP.PI.S.0
1	LP.PI.S.1.L	LP.PI.S.1.R	LP.PI.S.1
2	LP.PI.S.2.L	LP.PI.S.2.R	LP.PI.S.2
3	LP.PI.S.3.L	LP.PI.S.3.R	LP.PI.S.3
4	LP.PI.S.4.L	LP.PI.S.4.R	LP.PI.S.4
5	LP.PI.S.5.L	LP.PI.S.5.R	LP.PI.S.5
6	LP.PI.S.6.L	LP.PI.S.6.R	LP.PI.S.6
7	LP.PI.S.7.L	LP.PI.S.7.R	LP.PI.S.7
8	LP.PI.S.8.L	LP.PI.S.8.R	LP.PI.S.8
9	LP.PI.S.9.L	LP.PI.S.9.R	LP.PI.S.9
10	LP.PI.S.10.L	LP.PI.S.10.R	LP.PI.S.10
11	LP.PI.S.11.L	LP.PI.S.11.R	LP.PI.S.11
12	LP.PI.S.12.L	LP.PI.S.12.R	LP.PI.S.12



Sherling Pad - Universal One Size

Product Code	Description
LP.ZPT.GR	Grey Shearling Padded Tongue Pack 10
LP.ZPT.PI	Pink Shearling Padded Tongue Pack 10
LP.ZPT.BU	Burgundy Shearling Padded Tongue Pack 10
LP.ZPT.RB	Royal Blue Shearling Padded Tongue Pack 10
LP.ZPT.NB	Navy Blue Shearling Padded Tongue Pack 10

- All bars are fully assembled
- Simple quick release mechanism
- Set the width of the bar and the angle to your patient preference
- Fasten the boots to the bar with the click system.



Product Code: LP.SB.0.S
Static Bar Small - 0 Degree



Product Code: LP.SB.0.L
Static Bar Large - 0 Degree



Product Code: LP.SB.15.S
Static Bar Small - 15 Degree



Product Code: LP.SB.15.L
Static Bar Large - 15 Degree



Product Code: LP.AB.0.S
Dynamic Bar Small



Product Code: LP.AB.0.L
Dynamic Bar Large




Size	Length	Width
000000	5.5-6.5	2.5-4.0
00000	5.8-6.9	3.0-4.5
0000	6.0-7.3	3.3-4.8
000	6.3-7.8	3.0-5.5
00	6.8-8.3	3.0-5.5
0	7.5-9.0	3.0-6.0
1	8.0-9.5	3.5-6.0
2	9.0-10.5	3.8-6.5
3	10.0-11.5	4.5-7.0
4	11.0-12.5	4.7-7.0

Size	Length	Width
5	12.0-13.5	4.8-7.0
6	13.0-14.5	5.0-7.5
7	14.0-15.5	5.5-8.0
8	15.0-16.5	6.7-8.2
9	16.0-17.5	7.0-8.5
10	17.0-18.5	7.4-8.9
11	18.0-19.5	7.7-9.2
12	19.0-20.5	8.0-9.5
Small Bar	145mm-212mm	
Large Bar	217mm-355mm	

Made to Order Specification Guide

Please clearly state the following options when ordering your chosen style. Not all options are applicable to every style, please check the specification provided alongside each shoe/ boot in the catalogue pages.

- 1 Choose from either a **lace** or **reverse hook and loop** fastening (hook and loop incurs an additional charge).
- 2 Specify the correct height:




- 3 **5 standard, 3 anti-varus** and **3 orthoses accommodating last** options are available. See page 25 for measurement charts. (Split sizes are available but will incur an additional charge.)
- 4 **3 stiffener** options are available:



- 5 **Last additions: a maximum of 3 additions per last from the 7 illustrated** are available with modifications in multiples of 3mm, allowing you to choose any combination as illustrated below. Total modifications must not exceed 18mm per pair, each modification is charged as one unit with the code LADP-PAEDS.

If the limits are exceeded then a fully drafted Made To Measure order is required. If the footwear is to be size 32 or smaller, a cast is also required.
All footwear is supplied complete fittings are not available.
Last adaptions are not available in Accommodating Orthoses footwear.



- 6 **5 sole options** are available in the size ranges shown below:



Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.


- 7 **Various colour and material options** can be specified for each style, as illustrated on page 7.
- 8 **All 5 standard last width MTO footwear** is supplied with the following internal inlay.



Adaptations and Insoles

A17-PAEDS	High boot uppers additional 25mm to footwear
B1	Outside compensating raise up to 25mm with a rand
B2	Each subsequent 25mm or part thereof to B1
B3	Outside compensating raise up to 25mm without a rand
B4	Each subsequent 25mm or part thereof to B3
B5	Raising the heel in a single item of footwear up to 25mm
B7	Addition of a heel seat or valgus/varus wedge fitted to shoe (each)
B8-PAEDS	Stirrup, heel or caliper retaining strap (each)
B9	Any additional padding in a single item of footwear
B15-PAEDS	Addition of orthosis stiffener
XSTRG-PAEDS	Extra Strength Stability Stiffener
B16	Addition of a metatarsal bar fitted to the sole (each)
B17	Inclusion of a rocker sole (each)
B18	Inclusion of wedge, fitted into/onto the sole (each)
B19	Inclusion of a floated sole in a single item of footwear
B20	Inclusion of a sole strip or plate in composite material (each)
B21	Inclusion of a sole strip or plate in non-composite material (each)
B22	Addition of a rubber tip to the sole (each)
B23	Addition of a tip, with flange, to protect the welt (each)
B24	Addition of an external protective toe cap (each)
B25	Full/half-elongated, oblique, floated heel or heel wedge or torque heel in a single item of footwear
B30	Socket to one/both sides of the heel, without front/backstop to single item of footwear
B31	Socket to one/both sides of the heel, with front/backstop to single item of footwear
B32	Inclusion of a toe raising socket to single item of footwear
B34-PAEDS	Addition of a lined and padded double-decker T or Y strap (each)
B35	Rebuilding the heel to allow socket fixation in a single item of footwear
B40a	3/4 or full length non-moulded inlay, non-shock absorb material e.g. Leather/Neoprene/Plastazote® (each)
B40b	3/4 or full length non-moulded inlay, shock absorb material e.g. Poron®/PPT® (each)
B40c	3/4 or full length non-moulded inlay in silicone type material (each)
B41	Metatarsal/valgus or heel pad to inlay (each)
B42	Metatarsal/valgus or heel pad fitted direct to shoe (each)
B43	Inclusion of a heel elevator in cork, HDP or EVA covered (each)
B50	Filler to compensate for toe amputation(s) (each)
B51	Filler to compensate for amputation of forepart of foot (each)
LADP-PAEDS	Last supplement per adjustment (each)
NSB030	Flat carbon fibre insole (each)
PEADHL	Paediatric Footwear H&L Fastening Option
SURMENDZR	Paediatric Footwear Split Size Charge up to Size 38
SURMENDS	Paediatric Footwear Split Size Charge Over Size 39
WLN-PAEDS	Addition of lambswool lining in a single item of footwear (each)

(The dark blue lines on the diagram to the right illustrate the area lined.)




Colour Options

Not all colour options are applicable to every style, please check the specification provided alongside each shoe/boot in the catalogue pages.




Leather - Natural
Colour - Black



Leather - Natural
Colour - Navy



Leather - Nubuck
Colour - Tartufo



Leather - Box
Colour - White



Leather - Nubuck
Colour - Black



Leather - Nubuck
Colour - Chine



Leather - Val
Colour - Linen



Leather - Nubuck
Colour - Pink



Leather - Patent
Colour - Black



Leather - Natural
Colour - Royal




Leather - Nubuck
Colour - Light Brown




Leather - Bolonia
Colour - Pink



Leather - Calgary
Colour - Black



Leather - Nubuck
Colour - Navy




Leather - Val
Colour - Soan



Leather - Nubuck
Colour - Fuchsia



Leather - Nubuck
Colour - Grey



Leather - Antic
Colour - Jean



Leather - Val
Colour - Saddale



Leather - Bolonia
Colour - Fuchsia



Leather - Bolonia
Colour - Lilac

Ordering Information

1. Ordering and Enquiries:

Ordering Paediatric footwear from Blatchford is simple and easy. You can contact our Customer Services directly:

Telephone: +44 (0) 114 2637900

Email: orthoticscs@blatchford.co.uk

Facsimile: +44 (0) 114 2637901

Post: Blatchford Limited, 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

2. Delivery Promise:

Our goal is to provide fast and efficient delivery. We aim to despatch all orders for Paediatric Made to Order footwear from receipt in 14 working days. This does not apply where specific delivery schedules have been previously agreed.

3. Catalogues/Product Information:

If you require additional catalogues or product information please contact our customer service advisors who will be happy to help with your enquiry.

4. Single Patient Use:

All footwear featured in this catalogue is prescribed/ distributed by Blatchford for single patient use only and should be fitted by a qualified healthcare professional. The reuse of any footwear can potentially cause harm to the user and will invalidate the UKCA/CE or any warranty agreement. Please refer to the IFU for more information.

5. Delivery Shortages:

If there are shortages in the delivery this must be notified in writing to customer services within 5 working days of receipt of goods.

6. Returns:

In the unlikely event of Made to order footwear being supplied faulty or having been incorrectly supplied against the order specification, please contact customer services within 20 working days of delivery to report the issue so a credit can be processed.

Blatchford reserves the right to reject items returned for credit if the above conditions are not met. Blatchford is unable to accept return of items that are 'made to order' for credit, unless the goods are faulty or it can be demonstrated they have been damaged in transit. Products and prices are subject to variation at any time and are subject to availability. E & OE.



Need footwear for school?

No problem!

The majority of styles shown in this catalogue are available in a range of black leathers, making them perfect for school days.

Kiwi



- Last Options**
- 5 standard widths
- Sizes**
- 16-45 (subject to sole option)
- Anti-Varus Last**
- 3 widths - sizes 16-33 (subject to sole option)
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

Rubber Kick

(Sizes 20-45)

Black Light Beige

PU Sole

(Sizes 16-45)

Black White Beige Grey Pink



Leather - Natural
Colour - Navy




Leather - Natural
Colour - Black



Leather - Patent
Colour - Black



Leather - Nubuck
Colour - Fuchsia



Leather - Bolonia
Colour - Pink



Leather - Val
Colour - Linen

Skylark



- Last Options**
- 5 standard widths
- Sizes**
- 16-45
- Anti-Varus Last**
- 3 widths - sizes 16-33
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

Sports

(Sizes 16-45)

Black Grey White

Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Colour Options as illustrated:
White with Pink Trim or White with Blue Trim

Roadrunner



- Last Options**
- 5 standard widths
- Sizes**
- 16-45
- Anti-Varus Last**
- 3 widths - sizes 16-33
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

Sports

(Sizes 16-45)

Black Grey White

Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Colour Options as illustrated:
Blue or Pink

Albatross



- Last Options**
- 5 standard widths
- Sizes**
- 16-45
- Anti-Varus Last**
- 3 widths - sizes 16-33
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

Sports

(Sizes 16-45)

Black Grey White

Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Colour Options as illustrated:
White with Pink Trim or White with Blue Trim

Shikra



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Orthoses Accommodation

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Rubber Kick (Sizes 20-45) Black Light Beige

Not suitable with extra strength stability stiffener.






Leather - Bolonia Colour - Pink



Leather - Natural Colour - Royal



Leather - Bolonia Colour - Lilac



Leather - Nubuck Colour - Grey



Leather - Nubuck Colour - Chine



Leather - Nubuck Colour - Tartufo

Owl



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Orthoses Accommodation

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Not suitable with extra strength stability stiffener.





Leather - Nubuck Colour - Fuchsia



Leather - Nubuck Colour - Black



Leather - Nubuck Colour - Chine



Leather - Nubuck Colour - Tartufo



Leather - Natural Colour - Black



Leather - Natural Colour - Royal

Raptor



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45) Black White Beige Grey Pink


Rubber Montana (Sizes 18-45) Black Beige



Leather - Bolonia Colour - Lilac



Leather - Calgary Colour - Black



Leather - Nubuck Colour - Tartufo



Leather - Natural Colour - Black



Leather - Natural Colour - Royal

Eagle



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Rubber Montana (Sizes 18-45) Black Beige

Rubber Kick (Sizes 20-45) Black Light Beige



Leather - Nubuck Colour - Pink



Leather - Calgary Colour - Black



Leather - Nubuck Colour - Grey



Leather - Bolonia Colour - Fuchsia



Leather - Nubuck Colour - Fuchsia

Kite



- Last Options** 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Rubber Montana (Sizes 18-45) Black Beige

Rubber Omega (Sizes 16-42) Black

Rubber Kick (Sizes 20-45) Black Light Beige



Leather - Val
Colour - Linen Upper/
Cream Stitching



Leather - Nubuck
Colour - Fuchsia Upper/
White Stitching



Leather - Natural
Colour - Black Upper/
Cream Stitching



Leather - Nubuck
Colour - Black Upper/
Cream Stitching



Leather - Natural
Colour - Royal Upper/
White Stitching



Leather - Val
Colour - Saddale/
Brown Stitching

Osprey



- Last Options** 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Rubber Kick (Sizes 20-45) Black Light Beige



Leather - Nubuck
Colour - Fuchsia



Leather - Patent
Colour - Black



Leather - Calgary
Colour - Black



Leather - Natural
Colour - Royal



Leather - Bolonia
Colour - Lilac

Condor



- Last Options** 5 standard widths and 3 orthoses widths
Sizes 18-45
Anti-Varus Last 3 widths - sizes 18-33
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height Continental
Stiffeners Standard, Orthoses Accommodation

Rubber Montana (Sizes 18-45) Black (with yellow stitch effect)

Not suitable with extra strength stability stiffener.



Leather - Calgary
Colour - Black



Leather - Bolonia
Colour - Lilac



Leather - Patent
Colour - Black



Leather - Natural
Colour - Royal



Leather - Natural
Colour - Black

Harrier



- Last Options** 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45) Black White Beige Grey Pink

Rubber Montana (Sizes 18-45) Black Beige

Rubber Omega (Sizes 16-42) Black



Leather - Nubuck
Colour - Chine



Leather - Patent
Colour - Black



Leather - Calgary
Colour - Black



Leather - Natural
Colour - Royal



Leather - Natural
Colour - Black

Bateleur



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole
(Sizes 16-45) Black White Beige Grey Pink

Rubber Montana
(Sizes 18-45) Black Beige

Rubber Kick
(Sizes 20-45) Black Light Beige

Leather - Patent
Colour - Black

Leather - Calgary
Colour - Black

Leather - Natural
Colour - Black

Leather - Natural
Colour - Royal

Leather - Val
Colour - Saddale

Leather - Bolonia
Colour - Pink

Gryfalcon



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole
(Sizes 16-45) Black White Beige Grey Pink

Rubber Omega
(Sizes 16-42) Black

Rubber Kick
(Sizes 20-45) Black Light Beige

Only available in colour options below.

Leather - Nubuck
Colour - Black Upper/Box
White Toe Cap and Detail

Leather - Nubuck
Colour - Black Upper/Toe
Cap/Detail

Leather - Nubuck
Colour - Fuchsia Upper/Toe
Cap/Detail

Leather - Nubuck
Colour - Pink Upper/Toe
Cap/Detail

Leather - Val
Colour - Saddale Upper/
Val Soan Toe Cap and Detail

Leather - Antic
Colour - Jean Upper/Box
White Toe Cap and Detail

Kestrel



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation

PU Sole
(Sizes 16-45) Black White Beige Grey Pink

Rubber Omega
(Sizes 16-42) Black

Rubber Kick
(Sizes 20-45) Black Light Beige

Rubber Montana
(Sizes 18-45) Black Beige

Not suitable with extra strength stability stiffener.

Leather - Nubuck
Colour - Fuchsia

Leather - Patent
Colour - Black

Leather - Bolonia
Colour - Pink

Leather - Nubuck
Colour - Black

Leather - Natural
Colour - Black

Leather - Val
Colour - Saddale

Kookaburra



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole
(Sizes 16-45) Black White Beige Grey Pink

Rubber Kick
(Sizes 20-45) Black Light Beige

Only available in colour options below.

Leather - Nubuck
Colour - Chine Upper/Natural
Royal Detail/Box White Toe Cap

Leather - Nubuck
Colour - Pink Upper/Bolonia
Pink Detail and Toe Cap

Leather - Nubuck
Colour - Black Upper/Patent
Black Detail/Box White Toe Cap

Leather - Calgary
Colour - Black Upper/Detail/
Toe Cap

Leather - Natural
Colour - Black Upper/Detail/
Toe Cap

Leather - Val
Colour - Saddale Upper and
Toe Cap/Soan Detail

Merlin



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Only available in colour options below.



Leather - Antic
Colour - Jean Upper/
Nubuck Chine Cuff



Leather - Natural
Colour - Black Upper/
Black Nubuck Cuff




Leather - Nubuck
Colour - Fuchsia



Leather - Natural
Colour - Royal



Leather - Val
Colour - Saddale



Leather - Box
Colour - White

Peregrine



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Only available in colour options below.



Leather - Natural
Colour - Black Upper/
Black Nubuck Cuff



Leather - Val
Colour Saddale Upper/
Val Soan Cuff



Leather - Nubuck
Colour - Chine Upper/
Orchid Nubuck Cuff



Leather - Nubuck
Colour - Tartufo



Leather - Val
Colour - Linen



Leather - Nubuck
Colour - Pink

Besra



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability




Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.



Leather - Natural
Colour - Black Upper/
Box White Stripe



Leather - Antic
Colour - Jean Upper/
Box White Stripe



Leather - Box
Colour - White Upper/
Natural Black Stripe



Leather - Val
Colour - Saddale Upper/
Soan Stripe

Bazza




Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability




Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.



Leather - Natural
Colour - Black Upper/
Box White Stripe



Leather - Nubuck
Colour - Grey Upper/
Box White Stripe



Leather - Box
Colour - White Upper/
Natural Black Stripe



Leather - Nubuck
Colour - Fuschia Upper/
Box White Stripe



Leather - Nubuck
Colour - Chine Upper/
Box White Stripe

Hobby



- Last Options

5 standard widths and 3 orthoses widths
- Sizes

16-45
- Anti-Varus Last

3 widths - sizes 16-33
- Fastening

Lace or Reverse Pull Hook and Loop
- Height

UK or Continental
- Stiffeners

Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Sports

(Sizes 16-45)

Black

Beige

Grey

White

Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.

Leather - Val

Colour - Linen Upper/Inyet White Stripe

Leather - Antic

Colour - Jean Upper/Box White Stripe

Leather - Natural

Colour - Black Upper/Box White Stripe

Leather - Nubuck

Colour - Fuschia Upper/Inyet White Stripe

Leather - Nubuck

Colour - Grey Upper/Inyet White Stripe

Leather - Nubuck

Colour - Black Upper/Inyet White Stripe

Goshawk

Black rubber heel and toe protector is standard on design.



- Last Options

5 standard widths and 3 orthoses widths
- Sizes

16-45 (subject to sole option)
- Anti-Varus Last

3 widths - sizes 16-33 (subject to sole option)
- Fastening

Lace or Reverse Pull Hook and Loop
- Height

UK or Continental
- Stiffeners

Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Rubber Montana

(Sizes 18-45)

Black

Beige

Rubber Omega

(Sizes 16-42)

Black

Only available in colour options below.

Leather - Natural

Colour - Black

Leather - Natural

Colour - Navy

Leather - Antic

Colour - Jean

Leather - Patent

Colour - Black

Leather - Val

Colour - Linen

Buzzard



- Last Options

5 standard widths and 3 orthoses widths
- Sizes

16-45 (subject to sole option)
- Anti-Varus Last

3 widths - sizes 16-33 (subject to sole option)
- Fastening

Lace or Reverse Pull Hook and Loop
- Height

UK
- Stiffeners

Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Rubber Kick

(Sizes 20-45)

Black

Light Beige

Leather - Natural

Colour - Navy

Leather - Natural

Colour - Black

Leather - Patent

Colour - Black

Leather - Nubuck

Colour - Fuchsia

Leather - Bolonia

Colour - Pink

Leather - Val

Colour - Linen

Jay



- Last Options

5 standard widths and 3 orthoses widths
- Sizes

16-45 (subject to sole option)
- Anti-Varus Last

3 widths - sizes 16-33
- Fastening

Lace only with medial zip
- Height

UK or Continental
- Stiffeners

Standard, Orthoses Accommodation

PU Sole

(Sizes 16-45)

Black

Rubber Omega

(Sizes 16-42)

Black

Not suitable with extra strength stability stiffener.

Only available in colour options below.

Leather - Natural

Colour - Black

Leather - Patent

Colour - Black

Eider



- Last Options**
- 5 standard widths and 3 orthoses widths
- Sizes**
- 16-45 (subject to sole option)
- Anti-Varus Last**
- 3 widths - sizes 16-33 (subject to sole option)
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Height**
- UK
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Rubber Omega

(Sizes 16-42)

Black

Rubber Kick

(Sizes 20-45)

Black

Light Beige

Rubber Montana

(Sizes 18-45)

Black

Beige

Leather - Natural Colour - Black

Leather - Patent Colour - Black

Leather - Val Colour - Soan

Leather - Val Colour - Saddale

Leather - Nubuck Colour - Pink

Raven



- Last Options**
- 5 standard widths and 3 orthoses widths
- Sizes**
- 16-45 (subject to sole option)
- Anti-Varus Last**
- 3 widths - sizes 16-33 (subject to sole option)
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Height**
- UK
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Rubber Omega

(Sizes 16-42)

Black

Rubber Kick

(Sizes 20-45)

Black

Light Beige

Leather - Natural Colour - Black

Leather - Patent Colour - Black

Leather - Nubuck Colour - Black

Leather - Nubuck Colour - Grey

Leather - Nubuck Colour - Fuchsia

Leather - Val Colour - Linen

Cub Starter Boot



- Last Options**
- 5 standard widths
- Sizes**
- 16-26
- Anti-Varus Last**
- 3 widths - sizes 16-26
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Height**
- UK
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole

(Sizes 16-45)

Black

White

Beige

Grey

Pink

Leather - Nubuck Colour - Black

Leather - Nubuck Colour - Navy

Leather - Nubuck Colour - Fuchsia

Leather - Nubuck Colour - Pink

Leather - Val Colour - Soan

Leather - Nubuck Colour - Chine

Slickers Boot II



- Last Options**
- 5 standard widths
- Sizes**
- 18-45
- Anti-Varus Last**
- 3 widths - sizes 18-33
- Fastening**
- Lace or Reverse Pull Hook and Loop
- Height**
- UK
- Stiffeners**
- Standard, Orthoses Accommodation and Extra Strength Stability

Rubber Montana

(Sizes 18-45)

Black

Beige

Leather - Natural Colour - Black

Leather - Val Colour - Saddale

Leather - Nubuck Colour - Black

Leather - Val Colour - Soan

Leather - Natural Colour - Royal

Leather - Nubuck Colour - Pink


Sprinter Boot



Last Options 5 standard widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Trainer sport soles are not suitable for use on footwear where modifications to the last would alter the original tread-line. This includes any length or width adjustments at the joints.



Catapult Boot II



Last Options 5 standard widths
Sizes 18-45
Anti-Varus Last 3 widths - sizes 18-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Saker Sandal



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Fastening Reverse Pull Hook and Loop
Height Continental
Stiffeners Standard, Orthoses Accommodation

When ordering size down by 2 sizes.




Not suitable with extra strength stability stiffener. Only available in colour options below.



Measurement Charts

Fitting Principles

The measurement charts lay out a step by step procedure which can help the clinical assessment and measurement process achieve a right first time fit. It is important to note that when the foot weight bears, the foot posture changes and can have an effect on the shoe size required. Please consider this to ensure the desired toe room is added to the overall length.



Standard Last

Size		Narrow				Medium				Wide				X Wide				XX Wide			
Size	Length	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth
16	106	46	122	104	14	49	130	106	14	52	138	119	15	55	146	120	15	58	154	122	16
17	113	47	126	109	15	50	134	112	15	53	142	124	16	56	150	125	16	59	158	128	17
18	119	48	130	113	15	51	138	117	16	54	146	129	16	57	154	129	17	60	162	133	17
19	126	50	134	118	16	53	142	122	16	56	150	134	17	59	158	134	17	62	166	138	18
20	133	51	138	122	16	54	146	128	17	57	154	138	17	60	162	138	18	63	170	144	18
21	139	53	142	128	17	56	150	134	17	59	158	143	18	62	166	144	18	65	174	149	19
22	146	54	146	134	18	57	154	139	18	60	162	148	19	63	170	149	19	66	178	154	20
23	153	55	150	140	18	58	158	145	19	61	166	153	19	64	174	155	20	67	182	159	20
24	159	56	154	146	19	59	162	150	19	62	170	158	20	65	178	160	20	68	186	164	21
25	166	58	158	150	19	61	166	156	20	64	174	162	20	67	182	166	21	70	190	168	21
26	173	60	162	154	20	63	170	161	20	66	178	167	21	69	186	171	21	72	194	173	22
27	179	62	166	158	21	65	174	166	21	68	182	172	22	71	190	176	22	74	198	178	23
28	186	64	170	162	21	67	178	171	22	70	186	177	22	73	194	181	23	76	202	183	23
29	193	65	174	166	22	68	182	176	22	71	190	182	23	74	198	186	23	77	206	188	24
30	199	66	178	170	22	69	186	180	23	72	194	188	23	75	202	190	24	78	210	193	24
31	206	67	182	175	23	70	190	185	23	73	198	193	24	76	206	196	24	79	214	198	25
32	213	69	186	180	24	72	194	190	24	75	202	198	25	78	210	202	25	81	218	205	26
33	219	71	190	185	24	74	198	195	25	77	206	203	25	80	214	208	26	83	222	212	26
34	226	72	194	190	25	75	202	200	25	78	210	208	26	81	218	214	26	84	226	219	27
35	232	74	198	196	25	77	206	205	26	80	214	218	26	83	222	220	27	86	230	226	27
36	239	76	202	202	26	79	210	210	26	82	218	224	27	85	226	227	27	88	234	232	28
37	246	77	206	207	27	80	214	215	27	83	222	229	28	86	230	232	28	89	238	237	29
38	252	79	210	213	27	82	218	220	28	85	226	235	28	88	234	238	29	91	242	243	29
39	259	80	214	219	28	83	222	225	28	86	230	240	29	89	238	244	29	92	246	248	30
40	266	82	218	224	28	85	226	230	29	88	234	246	29	91	242	252	30	94	250	254	30
41	272	84	222	233	29	87	230	238	29	90	238	253	30	93	246	258	30	96	254	260	31
42	279	85	226	242	30	88	234	246	30	91	242	259	31	94	250	263	31	97	258	269	32
43	286	87	230	252	30	90	238	254	31	93	246	266	31	96	254	269	32	99	262	277	32
44	292	89	234	260	31	92	242	252	31	95	250	272	32	98	258	274	32	101	266	286	33
45	299	90	238	268	31	93	246	270	32	96	254	278	32	99	262	280	33	102	270	294	33

Measurements for all sections are shown in mm.



Anti-Varus Last

Size		Medium			Wide			X Wide		
Size	Length	Joint Width	Instep Circ.	Ankle Circ.	Joint Width	Instep Circ.	Ankle Circ.	Joint Width	Instep Circ.	Ankle Circ.
16	106	49	130	106	52	138	119	55	146	120
17	113	50	134	112	53	142	124	56	150	125
18	119	51	138	117	54	146	129	57	154	129
19	126	53	142	122	56	150	134	59	158	134
20	133	54	146	128	57	154	138	60	162	138
21	139	56	150	134	59	158	143	62	166	144
22	146	57	154	139	60	162	148	63	170	149
23	153	58	158	145	61	166	153	64	174	155
24	159	59	162	150	62	170	158	65	178	160
25	166	61	166	156	64	174	162	67	182	166
26	173	63	170	161	66	178	167	69	186	171
27	179	65	174	166	68	182	172	71	190	176
28	186	67	178	171	70	186	177	73	194	181
29	193	68	182	176	71	190	182	74	198	186
30	199	69	186	180	72	194	188	75	202	190
31	206	70	190	185	73	198	193	76	206	196
32	213	72	194	190	75	202	198	78	210	202
33	219	74	198	195	77	206	203	80	214	208

Measurements for all sections are shown in mm.

Orthoses Accommodating Footwear

The styles below feature an additional zip fastening for ease of access and donning.

Falcon



Last Options 3 orthoses widths
Sizes 20-44
Fastening Lace with Zip
Height UK
Stiffeners Orthoses Accommodation
No last adaptation available on Accommodating Footwear.

Rubber Montana (Sizes 18-45) Black (with yellow stitch effect)
Rubber Kick (Sizes 20-45) Black (with yellow stitch effect)

Not suitable with extra strength stability stiffener.



Leather - Calgary
Colour - Black



Leather - Nubuck
Colour - Black



Leather - Natural
Colour - Black



Leather - Natural
Colour - Royal



Leather - Nubuck
Colour - Brown


Hawk



Last Options 3 orthoses widths
Sizes 20-45
Fastening Lace or Reverse Pull Hook and Loop with Zip
Height UK
Stiffeners Orthoses Accommodation
No last adaptation available on Accommodating Footwear.

Sports (Sizes 16-45) Black Beige Grey White

Not suitable with extra strength stability stiffener.




Leather - Box
Colour - White



Leather - Natural
Colour - Black



Leather - Natural
Colour - Royal



Leather - Bolonia
Colour - Pink



Leather - Bolonia
Colour - Fuchsia

Swan Rear Opening Boot



Last Options 3 orthoses widths
Sizes 20-45
Fastening Lace Only
Height UK
Stiffeners Orthoses Accommodation Only
No last adaptation available on Accommodating Footwear.

Sports (Sizes 16-45) Black Grey White

Not suitable with extra strength stability stiffener.

Colour Options as illustrated:
White with Blue Trim or White with Red Trim



Orthosis Accommodating Last

Size		Orthoses Narrow Last				Orthoses Wide Last				Orthoses Extra Volume			
Size	Length	Joint Width	Joint Circ.	Ankle Circ.	Long Heel	Joint Width	Joint Circ.	Ankle Circ.	Long Heel	Joint Width	Joint Circ.	Ankle Circ.	Long Heel
20	133	57	170	161	225	60	178	168	235	60	184	168	241
21	139	59	174	167	233	62	182	174	243	62	188	174	249
22	146	60	178	173	241	63	186	180	251	63	192	180	257
23	153	61	182	179	249	64	190	186	259	64	196	186	265
24	159	62	186	185	257	65	194	192	267	65	200	192	273
25	166	64	190	191	265	67	198	198	275	67	204	198	281
26	173	66	194	197	273	69	202	205	283	69	208	205	289
27	179	68	199	204	281	71	207	212	291	71	213	212	297
28	186	70	203	211	289	73	211	219	299	73	217	219	305
29	193	71	207	218	297	74	215	226	307	74	221	226	313
30	199	72	212	225	305	75	220	233	315	75	226	233	321
31	206	73	217	234	315	76	225	242	325	76	231	242	331
32	213	75	221	243	323	78	229	251	333	78	235	251	339
33	219	77	225	252	331	80	233	260	341	80	239	260	347
34	226	78	229	261	339	81	237	269	349	81	243	269	355
35	232	80	233	268	347	83	241	278	357	83	247	278	363
36	239	82	237	279	355	85	245	287	365	85	251	287	371
37	246	83	240	288	365	86	248	296	375	86	254	296	381
38	252	85	244	297	375	88	252	305	385	88	258	305	391
39	259	86	248	306	385	89	256	314	395	89	262	314	401
40	266	88	252	315	395	91	260	323	405	91	266	323	411
41	272	90	255	324	405	93	263	332	415	93	269	332	421
42	279	91	258	333	415	94	266	341	425	94	272	341	431
43	286	93	261	342	425	96	269	350	435	96	275	350	441
44	292	95	264	351	435	98	272	359	445	98	278	359	451
45	299	96	268	360	444	99	275	368	454	99	281	368	460

Measurements for all sections are shown in mm.

To meet the clinical requirement of a patient Blatchford have the right to modify the UKCA marked footwear and by default this will be labeled as Custom-made device. Blatchford maintains a QMS system which is ISO 9001 and ISO 13485 accredited. Blatchford products and systems meet the essential requirements of EU Medical Devices Regulation 2017 (MDR).

Distributed by



blatchfordmobility.com | orthoticscs@blatchford.co.uk | [@blatchfordgrp](https://twitter.com/blatchfordgrp)