

Tel: +44 (0) 114 2637900 | Fax: +44 (0) 114 2637901 | orthoticscs@blatchford.co.uk

Patient Information

Patient Initials _____

Clinic _____

Clinician _____

☐ Left leg ☐ Right leg ☐ Bilateral

Trade/Orthotist/Practitioner Information

(Person who took the measurements)

Clinic/Customer Name _____

Customer PO/Order No. _____

Telephone No. _____

Email Address _____

Was a Trial Kit Used?

☐ Yes ☐ No ☐ Small
☐ Medium
☐ Large

☐ Assembled
XTERN Summit Required
☐ Small ☐ Medium
☐ Large

☐ Non-assembled
XTERN Summit Required
☐ Small ☐ Medium
☐ Large

Assembled and Non-assembled braces are the same price.

Measurements

Shoe Size _____

☐ Male ☐ Female

Approx Height _____

Approx weight in Kg _____

A

Forefoot Total Width:

☐ Change width of
forefoot section

☐ _____ mm narrower

☐ _____ mm wider

B

Heel Total Width:

☐ Change width of
heel section

☐ _____ mm narrower

☐ _____ mm wider

C

Front Length (tick one):

Caution - the XTERN
is set up to the same
measures on both sides.

☐ 0 ☐ 1 ☐ 2 ☐ 3

☐ 4 ☐ 5 ☐ 6 ☐ 7


If the form is not completed correctly it could lead to a possible delay in processing the device for shipment. Once a brace has been manufactured, this is a custom-made device and is non-returnable.

D

Calf Height: Slide the calf band at the desired position on patient calf. On the back of the orthosis, measure the distance between the 'Y' intersection and the top of the plastic calf band. _____ mm

Accessories and Options (Extra charge for all items, see price list for details.)

Item	Size	Part No.	Quantity
☐ Lace Clip Attachment	Small	200000-L00-S	
☐ Lace Clip Attachment	Med/Large	200000-L00-ML	
☐ Extra Cable Ties	Pack of 1000	700000-ZZZ-002	
☐ Extra Cable Ties	Pack of 20	700000-ZZZ-001	
☐ Summit XTERN Moulded Calf Pad/Strap	Med/Large	501302-K00-ML	
☐ Summit XTERN Moulded Calf Pad/Strap	Small	501302-K00-S	
☐ XTERN Magnet Buckle Strap	S/M/L	501433-MBS-SML	
☐ Ankle Stabilisation Strap	Universal	500100-AST-SML	
☐ Heel Bumper Kits (8 per kit)	S/M/L	500250-HBP-SML	
☐ Extension Stopper 40mm	S	500202-EXT-40	
☐ Extension Stopper 45mm	M/L	500202-EXT-45	



Additional Information

For technical support please email david.bownes@blatchford.co.uk
@blatchfordgrp | turbomedfootdrop.co.uk

Tel: +44 (0) 114 2637900 | Fax: +44 (0) 114 2637901 | orthoticscs@blatchford.co.uk
Blatchford Limited, 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

3320451012SF Iss8 09/26

Distributed by

Blatchford

TurboMed
ORTHOTICS